

Reynard IV access and clinical skin preparation. Our extensive range means you're prepped for anything.

Superior solution skin preparation products.

Reynard manufactures a broad range of intravenous access and skin preparation products giving you the perfect solution for every clinical situation.

Innovative packaging with 'easy-peel' and 'tear-open' options combined with unsurpassed value for money ensure our prep pads and swabsticks are used every day throughout Australia. Providing every health system provider in the country with the highest quality, TGA registered product solutions.



Reynard is reassurance.



Reynard
Health Supplies

Reynard IV access & clinical skin preparation.

The organisms may be micro
but the impact can be big.

Surgical site infections (SSIs) are a major challenge for patients and health sector resources. Effective preparation, cleansing, and disinfection using appropriate Reynard Prep Pads and Swabsticks can help minimise the risk of surgical site infections and improve overall patient outcomes.

- SSIs negatively impact patient physical and mental health.
- Increased patient morbidity, mortality and loss of earnings during recovery are some of the indirect costs associated with infection.
- Intangible costs may also be incurred by the patient, such as pain and anxiety. In addition, patients may experience delayed wound healing and be more susceptible to secondary complications, such as bacteraemia.
- Distress may also be caused to the patient and family members if the patient is absent from home and work for a prolonged period.
- Accordingly, prolonged hospitalisation and increased morbidity as a result of developing an SSI, have been shown to negatively impact on patient health-related quality of life.
- Increase in healthcare costs per patient and additional stress on healthcare systems and staff.

The perfect solution for
every problem.

Reynard provides the broadest range of prep pads and swabsticks in a variety of formulations and configurations to suit your specific clinical need and appropriate application.

- 0.5% chlorhexidine digluconate (CHG) and 70% isopropyl alcohol (IPA)
- 2% chlorhexidine digluconate (CHG) and 70% isopropyl alcohol (IPA)
- 70% isopropyl alcohol (IPA) - sterile
- 0.5% Chlorhexidine digluconate (CHG)
- 2% Chlorhexidine digluconate (CHG)



Reynard Health Supplies was established in New Zealand in 1999 and today is the number one supplier of surface cleaning & disinfection wipes, specialist patient care wipes and skin preparation products to hospitals in New Zealand and Australia.

Used extensively in the private sector and every public healthcare system in Australia.

At Reynard, we strive to provide products and support services that provide you with the reassurance you need to be able to focus on optimum patient outcomes.

- **Quality & range** - a wide range of quality products to suit every need and application
- **Cost savings** - cost effective and pricing stability provide real value
- **Reliability** - a reliable, proven supply chain from point of manufacture to delivery at your facility
- **Support** – technical expertise, training, education, and compliance support services
- **Research & development** - focus drives continuous product development and a pipeline of innovative solutions
- **Flexibility** - the capability and desire to design new products and supply chain options



Reynard
Health Supplies

Our extensive range means you're prepped for anything.

SIZE AND TYPE OF PRODUCT													
Formula	PREP PADS					SWABSTICKS							
	Substrate only					Cotton head				Foam head			
	6 x 3cm	6 x 6cm	4.5 x 8.5cm	16 x 20cm	Easy peel	Regular	Triple pack	Easy peel	Tinted	Regular	Maxi	Reg tinted	Maxi tinted
70% IPA only	RHS451G	RHS417	RHS438										
2% chg + 70% IPA		RHS401	RHS437	RHS422	RHS422EP	RHS408	RHS408T	RHS408EP	RHS408TR	RHS473		RHS473TR	RHS472
0.5% chg + 70% IPA										RHS470		RHS470TR	RHS471
2% chg only			RHS421							RHS475			
0.5% chg only		RHS479								RHS477	RHS478		

Skin preparation considerations.

PRODUCT FORMULATION	RECOMMENDED PRODUCT APPLICATION
0.5% chg + 70% IPA	General use, also suitable for paediatrics and areas where higher concentrations of CHG are not preferred due to potential allergies, skin reactivity and risk.
2% chg + 70% IPA	General use skin preparation where sensitivity etc. are not a concern. This solution provides broad capability for most patients.
70% IPA	General use, also suitable for paediatrics and areas where higher concentrations of CHG are not preferred due to potential allergies, skin reactivity and risk.
0.5% chg only	General use, also suitable for paediatrics and areas where higher concentrations of CHG are not preferred due to potential allergies, skin reactivity and risk. Alcohol (IPA) free, appropriate where use and application of IPA is not suitable.
2% chg only	General use skin preparation where CHG Concentration sensitivity is not a concern. Alcohol (IPA) free, appropriate where use and application of IPA is not suitable.

All Reynard clinical skin preparation products are suitable for a wide range of skin preparation requirements including but not limited to:

- Blood culture collections
- Central venous catheters insertion
- General skin preperation, pre-procedure, injection or incision
- Insertion of percutaneous devices
- Percutaneous device site maintenance
- Peripheral intravascular catheters insertion
- Site maintenance of central lines

Note: Only use on intact skin before any invasive procedures



Directions	Product should be applied with gentle friction using a repeated back-and-forth motion for manufacturer-recommended application time.
Dry Time	Allow product to air dry completely; dry time can be affected by volume of product, body site, presence or absence of hair, humidity etc.

Reynard Prep Pad and Swabstick Product Application

Instructions for use:

1. Check allergy adverse reaction status with your patient
2. Check that the chosen formula is compatible with the catheter material
3. Thoroughly clean and treat the area surrounding the insertion site with chlorhexidine and/or alcohol using a back-and-forth scrubbing motion for 20 to 30 seconds
4. Apply a “window frame” of skin preparation solution around the catheter site, if needed, to help the dressing stay in place (If you plan on placing a chlorhexidine-impregnated sponge dressing at the insertion site, do not apply the skin protectant to the area that the sponge will cover)
5. Use one or two pads or swab sticks dependent on the clinical procedure
6. Allow the skin area to fully air dry for approx 30 seconds for a dry site or up to 2 minutes for a moist site, to be most effective against bacteria growth
7. Discard of prep pads or swab sticks after single use, do not flush
8. Dispose of prep pads and swab sticks in adherence with your facility Infection Control Guidelines

Is the product flammable?

- Refer to manufacturer's instructions for use, and any associated hazard warnings
- Follow AORN guidelines for perioperative practice



Prep pads

CHLORHEXIDINE & ALCOHOL



Product code	RHS401	RHS422	RHS422EP	RHS437
Swabs per box	200	100	50	100
Boxes per carton	25	10	20	10
Wipe size	6 x 6cm	16 x 20cm	16 x 20cm	4.5 x 8.5cm
Solution	70% alcohol & 2% chg	70% alcohol & 2% chg	70% alcohol & 2% chg	70% alcohol & 2% chg
Volume	0.7ml	4.5ml	5ml	1.6ml
GSM	45gsm	40gsm	40gsm	100gsm

ALCOHOL




Product code	RHS417	RHS451G	RHS438
Swabs per box	200	200	100
Boxes per carton	25	25	10
Wipe size	6 x 6cm	6 x 3cm	4.5 x 8.5cm
Solution	70% alcohol	70% alcohol	70% alcohol
Volume	0.7ml	0.5ml	1.6ml
GSM	45gsm	45gsm	100gsm

CHLORHEXIDINE



Product code	RHS421	RHS479
Swabs per box	100	200
Boxes per carton	10	25
Wipe size	4.5 x 8.5cm	6 x 6cm
Solution	2% chg	0.5% chg
Volume	0.8ml	0.7ml
GSM	45gsm	45gsm

Skin surface area coverage for Reynard Prep pads*




22cm

Surface area Coverage for Reynard Prep Pads		
	Total Area	Products
Small Swabs	42cm ²	451G
Medium Swabs	100cm ²	401, 417, 479
Large Swabs	225cm ²	421, 437, 438
X-Large Swabs	484cm ²	422, 422EP

*note- Coverage area may vary in accordance with variations in atmospheric conditions and condition of the skin. Tests were conducted on hairless male skin at 21 degree Celsius. The skin area was wetted sufficiently to ensure it was wet for a minimum 30 seconds before allowing to dry. In each test enough liquid was left in the product to allow up to a further 20% area to be covered.

Swabsticks

COTTON

				
Product code	RHS408	RHS408T	RHS408TR	RHS408EP
Sticks per sachet	1	3	1	1
Sticks per box	50	75	50	50
Boxes per carton	10	10	10	10
Stick type	Cotton	Cotton	Cotton	Cotton
Solution	70% alcohol & 2% chg	70% alcohol & 2% chg	70% alcohol & 2% chg	70% alcohol & 2% chg
Tint	No	No	Yes	No
Volume	1.6ml	4.5ml	1.6ml	1.6ml


FOAM

						
Product code	RHS470	RHS470TR	RHS473	RHS473TR	RHS475	RHS477
Sticks per box	25	25	25	25	25	25
Boxes per carton	10	10	10	10	10	10
Stick type	Foam	Foam	Foam	Foam	Foam	Foam
Solution	70% alc & 0.5% chg	70% alc & 0.5% chg	70% alc & 2% chg	70% alc & 2% chg	2% chg	0.5% chg
Tint	No	Yes	No	Yes	No	No
Volume	2.5ml	2.5ml	2.5ml	2.5ml	2.5ml	2.5ml

MAXI FOAM

			
Product code	RHS471	RHS472	RHS478
Sticks per box	25	25	25
Boxes per carton	10	10	10
Stick type	Maxi Foam	Maxi Foam	Maxi Foam
Solution	70% alcohol & 0.5% chg	70% alcohol & 2% chg	0.5% chg
Tint	Yes	Yes	No
Volume	5ml	5ml	5ml

Skin surface area coverage for Reynard Swabsticks*



Surface area Coverage for Reynard Swab Stocks		
	Total Area	Products
Regular Cotton Swabsticks	72cm ²	408, 408T, 408TR, 408EP
Regular Foam Swabsticks	144cm ²	470, 470TR, 473, 473TR, 475, 477
Large Foam Swabsticks	484cm ²	471, 472, 478

*note- Coverage area may vary in accordance with variations in atmospheric conditions and condition of the skin. Tests were conducted on hairless male skin at 21 degree Celsius. The skin area was wetted sufficiently to ensure it was wet for a minimum 30 seconds before allowing to dry. In each test enough liquid was left in the product to allow up to a further 20% area to be covered.

Reynard focus on killing germs so you can focus on saving lives.

Germs.

They can proliferate in the cleanest of places. A constant bugbear, literally, to be eradicated, first time every time.

Reynard constantly adopt new technology, new research to develop new products in all our ranges.

Working to protect patients, and everyone and everything they come into contact with.

An ever growing Reynard product range means you can choose a product that meets every operational need, whether cleaning, disinfecting, IV Therapy or patient care.

Many applications developed in consultation with health professionals, in response to urgent clinical needs.

So you can concentrate on doing what you do best. Caring for patients. Visit www.reynardhealth.com.au

Reynard is reassurance.



Reynard Product ranges

Personal Care



Our diverse range of personal care products provide gentle and effective cleansing while minimising cross-infection risks.

Cleaning and disinfection



Developed with convenience in mind, this range provides fast and effective solutions for every cleaning needs.

Personal protective equipment



Reynard P.P.E. range aids in infection control and works towards minimising the spread of bacteria, dust and pathogens.

Dry wipes



Disposable, practical and absorbent, Reynard dry wipes are one of the softest and strongest wipes on the market.

Hand Hygiene



Kills 99.9% of germs while leaving hands soft and refreshed. Available with and without chlorhexidine.

Baby Care



Chemical and fragrance-free, our baby care range is designed to be ultra-soft and gentle on little ones.

References

1. Anderson, D., et al. Strategies to Prevent Surgical Site Infections in Acute Care Hospitals: 2014 Update. Infection Control and Hospital Epidemiology. 2014; 35(6), 605-627. doi:1. Retrieved from <http://www.jstor.org/stable/10.1086/676022> doi:1
2. Chaiyakunapruk N, Veenstra DL, Lipsky BA, Saint S. Chlorhexidine compared with povidone-iodine solution for vascular catheter-site care: a meta-analysis. Ann Intern Med 2002;136:792-801.
3. Hefford M, Black M, Wyatt S, Cost benefit analysis of the proposed national surgical site infection surveillance and response programme, 2011
4. Maki DG, Ringer M, Alvarado CJ. Prospective randomised trial of povidone-iodine, alcohol, and chlorhexidine for prevention of infection associated with central venous and arterial catheters. Lancet 1991;338:339-43.
5. O'Grady NP, Alexander M. Guidelines for the Prevention of Intravascular Catheter-related Infections. Recommendations of the Healthcare Infection Control Practices Advisory Committee (HICPAC). American Journal of Infection Control 2011; 39:4 (1 Suppl): S1- 34.
6. Prevention and treatment of surgical site infection. NICE Clinical Guidelines <https://www.nice.org.uk/guidance/cg74/evidence/full-guideline-242005933> Published October 22, 2008. Accessed December 13, 2016
7. Whitehouse et al. The impact of surgical-site infections following orthopedic surgery at a community hospital and a university hospital: adverse quality of life, excess length of stay, and extra cost. Infect Control Hosp Epidemiology. 2002; 23(4):183-189.
8. World Health Organization (WHO). Global Guidelines for the Prevention of Surgical Site Infection. 2016. <http://www.who.int/gpsc/ssi-prevention-guidelines/en/> (Accessed February 3, 2017).
9. Zimlichman E, Henderson D, Tamir O, et al. Health Care-Associated Infections: A Meta-analysis of Costs and Financial Impact on the U.S. Health Care System. JAMA Internal Medicine, Dec. 9/23, 2013. 173(22), 2042-2044.



Australia

44 Hampden Road
Artarmon, NSW, 2064
enquiries@reynardhealth.com.au
1300 667 199

New Zealand

PO Box 8470
Havelock North
Sales@reynardhealth.com
0508 494 737

UK

13 White Horse Business Park,
Stanford in the vale SN7 8NY
uk@reynardhealth.com
+44 (0)845 263 7335

International

PO Box 8470
Havelock North
Sales@reynardhealth.com
+64 (0)6 650 0709