

Tristel Jet™ FOAM

SPORICIDAL FOAM

Jet is a powerful sporicidal disinfectant that can be used on surfaces and medical equipment located near the patient.

As a foam, Jet's unique properties deliver chemistry activated at point of use without creating aerosols. This ensures fresh chemistry is delivered directly to the surface, without the risk of aerosols associated with sprays.

Jet contains powerful cleaning agents that destroy DNA and RNA in seconds. Its cleaning efficacy meets the requirements of a medical device cleaner, and non-ionic surfactants lift and remove organic soiling and prevent them from clinging back on surfaces, optimising disinfection action.


Jet utilises Tristel chlorine dioxide chemistry and is effective against microorganisms of concern. Jet passes bactericidal, virucidal, fungicidal and mycobactericidal tests in 30 seconds, and sporicidal tests in 2 minutes. Disinfection is completed well within the time taken for the chemistry to dry.

Jet is the best safeguard for patients and staff wherever the risks of infection are highest, and is ideal for areas such as isolation wards, operating theatres, intensive care units, oncology wards and neonatal wards.

EFFECTIVE
AGAINST
CORONAVIRUSES



Ideal for use in patient areas such as:



Jet achieves bactericidal, virucidal, fungicidal and mycobactericidal efficacy in 30 seconds, and sporicidal efficacy in 2 minutes, including against:

| BACTERIA | CONTACT TIME | VIRUSES | CONTACT TIME |
|---|--------------|---|--------------|
| <i>Staphylococcus aureus</i> | 30 seconds | Poliovirus Type 1 | 30 seconds |
| <i>Enterococcus hirae</i> | 30 seconds | Adenovirus Type 5 | 30 seconds |
| <i>Pseudomonas aeruginosa</i> | 30 seconds | Murine Norovirus | 30 seconds |
| Carbapenem Resistant Enterobacteriaceae (CRE) | 30 seconds | Influenza Virus (H1N1) | 30 seconds |
| <i>Klebsiella pneumoniae</i> | 30 seconds | Coronavirus (SARS-CoV-2)* | 30 seconds |
| Vancomycin Resistant Enterococci (VRE) | 30 seconds | Vaccinia Virus | 30 seconds |
| <i>Enterococcus faecium</i> | 30 seconds | Human Papillomavirus (using Polyoma virus SV40 Surrogate) | 30 seconds |
| | | Human Papillomavirus** | 30 seconds |
| | | Herpes Simplex Virus Type 1 | 30 seconds |
| MYCOBACTERIA | CONTACT TIME | | |
| <i>Mycobacterium terrae</i> | 30 seconds | | |
| <i>Mycobacterium avium</i> | 30 seconds | | |
| FUNGI | CONTACT TIME | SPORES | CONTACT TIME |
| <i>Candida albicans</i> | 30 seconds | <i>Clostridium sporogenes</i> | 2 minutes |
| <i>Aspergillus brasiliensis</i> | 30 seconds | | |

Easy to use with any available wipe, to reduce your plastic use!




=



Saving space, reducing waste

Jet can be used with any sustainably sourced wipe, avoiding the use and waste of plastic wipes.



Powered by chlorine dioxide foam

Jet is built upon the same core chlorine dioxide technology as the Tristel medical device disinfectant products, but specially adapted for surface disinfection where cleaning and biofilm destruction are so important.

Chlorine dioxide is widely regarded as one of the most effective disinfectants for biofilm removal and prevention.

Jet delivers chlorine dioxide as a foam and – unlike sprays – generates no aerosols.

Complete compliance with EN 14885:2018 and the latest efficacy regulatory requirements

| | | | |
|--------------|----------|---------------|----------|
| SPORES | EN 17126 | YEASTS/ FUNGI | EN 13624 |
| | EN 13704 | | EN 14562 |
| MYCOBACTERIA | EN 14348 | | EN 16615 |
| | EN 14563 | BACTERIA | EN 13727 |
| VIRUSES | EN 14476 | | EN 14561 |
| | EN 16615 | | EN 16615 |

**AVAILABLE
TO ORDER
NOW!**



Each pack contains: 2x Jet Heads and 2x 800ml Jet Packs

One 800ml Jet pack delivers 570 doses of foam

Quote: TSL0518 14

Jet is classified as biocide under the EU Biocides Regulation. Use biocides safely. Always read the label and product information before use.

Tristel™

Australia: Tristel Pty Ltd,
40/328 Reserve Road, Cheltenham, VIC 3192
T 1300 680 898 - F +61 (0)3 9533 6193 - E mail-au@tristel.co.uk

New Zealand: Tristel New Zealand Limited,
P.O Box 15451, Tauranga, New Zealand
T +64 (0)7 5771560 - F +64 (0)7 5771567 - E info@tristel.co.nz

Copyright © Tristel Solutions
Mkt-Bro-1906-1
23/FEB/2022

For Tristel patent information please visit: <http://www.our-patents.info/tristel>
*A representative sample of Tristel chlorine dioxide chemistry has been tested in accordance with EN 14476:2013+A2:2019, at a concentration of 20 parts per million (ppm). JET has a chlorine dioxide concentration greater than 20ppm at point of use. **Based on native infectious HPV Type 16 and Type 18 as tested by Meyers et al., (2020).