

THE WORLD'S MOST WIDELY APPROVED, VALIDATED AND USED HIGH-LEVEL DISINFECTION WIPES SYSTEM

For heat-sensitive, non-lumened medical devices, such as:

- Nasendoscopes
- Transvaginal and transrectal ultrasound probes
- Transoesophageal echocardiogram (TOE) probes
- Laryngoscopes
- Manometry catheters
- Ophthalmic devices



DESTROYS HPV TYPES 16 & 18 IN 30 SECONDS*

*Meyers, C., Milici, J., Robison, R. (2020) 'The Ability of Two Chlorine Dioxide Chemistries to Inactivate Human Papillomavirus-contaminated Endocavitary Ultrasound Probes and Nasendoscopes' DOI: 10.1002/jmv.25666.

THE TRISTEL TRIO WIPES SYSTEM IS EASY TO USE!



CLEAN

Tristel Pre-Clean Wipe

A non-woven wipe impregnated with a triple enzymatic detergent and surfactant used for a thorough cleaning of the medical device's surface to remove soil and organic matter.



HIGH-LEVEL DISINFECT

Tristel Sporidical Wipe

This wipe combines with the Tristel Activator Foam to generate chlorine dioxide, a highly effective biocide that eliminates bacteria, viruses, fungi, mycobacteria and spores from the surface of the medical device in 30 seconds.



RINSE

Tristel Rinse Wipe

A sterile packed, non-woven wipe impregnated with deionised water and a low antioxidant content; used to remove any remaining chemical residues from the surface of the medical device.



TRACE

Tristel 3T or the Quality Audit Trail Record Book

Intended for the simple registration of all steps of the process, including validation by the responsible person.

Tristel Trio™
Wipes System



Tristel Pre-Clean Wipe
Class I Medical Device



Tristel Sporidical Wipe
Class IIb Medical Device



Tristel Rinse Wipe
Class I Sterile Device

STERILE R



UNRIVALLED EFFICACY

The Tristel Sporicidal Wipe is proven effective as a high-level disinfectant according to BS EN 14885 and the latest efficacy standards published, in 30 seconds. The product's fast-acting and uniform contact time leaves no room for doubt or error.



FLEXIBLE

The Tristel Trio Wipes System is a mobile solution which can either be used in the consultation or decontamination room. Its flexibility eliminates the need for the instrument to be sent to central sterilisation units for reprocessing, saving valuable time.



COMPATIBLE

The Tristel Trio Wipes System has been tested and confirmed to be compatible with medical devices from over 30 manufacturers, including:

- Advanced Tactile Imaging, Inc.
- Airtraq*
- Arthrex
- ATMOS MedizinTechnik
- Becton Dickinson (Bard Access, CareFusion)
- BK Medical
- Butterfly Network
- Canon (Toshiba)
- Carestream
- Clarius Mobile Health
- Clarus Medical
- Crystal Photonics
- Echosens
- Esaote*
- Exact Imaging
- Fisher & Paykel
- GE Healthcare
- Haag-Streit
- Hitachi
- Karl Storz
- Laborie (Cogentix Medical, MMS)
- Machida Endoscope Co.
- MCube
- Medcom
- Medtronic
- Mindray
- Neoplas Tools
- Ocular Instruments Inc
- Optim LLC
- Orlvision
- Pentax
- Philips
- Polydiagnost
- Samsung Healthcare
- SapiMed
- Siemens Healthineers
- SonoScape*
- FUJIFILM SonoSite
- Unisensor
- Verathon

* Compatible with the Tristel Sporicidal Wipe as part of the Tristel Trio Wipes System



FAST

Full decontamination cycle, including traceability, in a matter of minutes.



SMART TRACEABILITY

The Tristel Trio Wipes System includes either paper-based or digital traceability. The widely-used Quality Audit Trail Record Book documents each disinfection event with pen and paper. Tristel 3T digitalises this process for a more efficient, accurate and smart approach to traceability.



ECONOMICAL

No investment required in machines and maintenance contracts.




SAFE

The Tristel Sporicidal Wipe is non-toxic to embryos and can therefore be used in IVF environments.

The Tristel Trio Wipes System offers complete decontamination of heat-sensitive, non-lumened medical devices. Thanks to the concept of cleaning, high-level disinfection and rinsing, a validated result is achieved with a simple and fast three-step process.



Microbiological test reports

| | European norms (EN) | Uniform Contact Time |
|---------------|---------------------|---|
| Bacteria | EN 13727 |  |
| | EN 14561 | |
| | EN 16615 | |
| Yeasts/ Fungi | EN 13624 | |
| | EN 14562 | |
| | EN 16615 | |
| Viruses | EN 14476 | |
| | EN 16615 | |
| Mycobacteria | EN 14348 | |
| | EN 14563 | |
| Spores | EN 16615 | |
| | EN 17126 | |

TRISTEL TRIO WIPES SYSTEM 50

- 1 box contains 50 decontamination procedures, available with or without the Quality Audit Trail Record Book.
- Tristel 3T is sold separately.



More information about the Tristel Trio Wipes System, such as safety data sheets, reports, publications on materials and studies, is available on request or online at www.tristel.com

MEDICAL DEVICE TRANSPORTATION

Tristel protect bags are designed specifically for the transportation and short-term storage of medical devices, with traceability in mind.

TRISTEL PROTECT

- 50 bags per box.
- Available in small (21cm x 34cm) for laryngoscopes and large (31cm x 71cm) for nasendoscopes.



Tristel
WE HAVE CHEMISTRY.

Manufactured in the United Kingdom by:
Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY
T +44 (0) 1638 721500 · E mail@tristel.com · W www.tristel.com

Hong Kong & Taiwan:
Tristel Asia Limited, 21st Floor, 168 Electric Road, Fortress Hill, Hong Kong
T +852 2895 6968 · F +852 2869 4388 · E customerservicehk@tristel.com

New Zealand: Tristel New Zealand Limited, 23 Birch Avenue, Judea, Tauranga
T +64 (0)7 5771560 · F +64 (0)7 5771567 · E mail-nz@tristel.com

Australia: Tristel Pty Ltd, 40/328 Reserve Road, Cheltenham, VIC 3192
T 1300 680 898 · F +61 (0)3 9533 6193 · E mail-au@tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>