

Quark RMR
Indirect Calorimetry

The gold standard for metabolic measurements in
applied human physiology research

“ Assess, Measure,
Improve my
Performance ”



COSMED
The Metabolic Company

“Quark RMR was demonstrated to be unbiased, precise, reproducible, and accurate device for measuring oxygen consumption⁽¹⁾”

- Measurement of Oxygen Consumption over the entire human physiological range
- Top of the range O₂/CO₂ gas analyzers (Paramagnetic, NDIR)
- Three flowmeters available for Resting and Exercise applications
- Gas Exchange measurements via Breath by Breath or optional Mixing Chamber
- Powered by OMNIA, the most intuitive software in the industry
- Independently validated at rest and exercise conditions



The Quark RMR is a state-of-the-art metabolic cart for gas exchange analysis (VO₂, VCO₂) either during resting or exercise.

The number of available configurations make the Quark RMR the most versatile metabolic cart for applied research in human physiology.

Quark RMR's accuracy and reliability have been validated against Gold Standard methods either with spontaneously breathing subjects (at rest and during exercise) and mechanically assisted patients.

Design

Unsurpassed reliability. Fast-response stable and durable paramagnetic technology for O₂ sensor, and rapid infrared for the CO₂. Both analyzers can ensure reliable data for a long time without requiring their replacement.

Breath by Breath & Mixing Chamber. Quark RMR is provided with Breath by Breath analyzers however the system is also available with an optional Mixing Chamber.

Modular architecture allows to configure Quark RMR according to the different metabolic testing requirements. This cost-effective solution gives the opportunity to scale at any time to a more complex configuration.

Low running costs and easy maintenance. Quark RMR design has been conceived to reduce ordinary maintenance and to easily and rapidly solve any possible technical problem through parts replacement.

Powered by OMNIA software innovative user interface, touch screen ready, easy-to-use and self-explanatory.

Calibrations and Verifications procedures are available to ensure that main measurement components perform according to their specifications.

Independently validated technology. Quark RMR is the only metabolic cart in the market that has been validated both on different gas exchange methods (Breath by Breath and Mixing Chamber), and on the whole human physiological range (from resting to a wide range of exercise intensities).

Main Features and Tests

Resting Energy Expenditure (REE) with Canopy Hood	Standard
---	----------

Resting Energy Expenditure (REE) with Mask	Standard
--	----------

Resting Energy Expenditure (REE) with Ventilator (for mechanically assisted patients)	Option
---	--------

“Breath by Breath” Cardio Pulmonary Exercise Testing (CPET) with Face mask	Option
--	--------

“Mixing Chamber” Cardio Pulmonary Exercise Testing (CPET) with Face mask/Mouthpiece	Option
---	--------

Spirometry (FVC, SVC, MVV etc.)	Option
---------------------------------	--------

⁽¹⁾ Ashcraft C.M. et al. “A Test of Validity of a New Open-Circuit Indirect Calorimeter.” *J Parenter Enteral Nutr.* 2014 Mar 10

Resting Energy Expenditure (REE)

The Quark RMR in its standard configuration provides the following features:

Breath by Breath Gas exchange Measurement of oxygen consumption (VO_2), carbon dioxide production (VCO_2) and related ventilatory and metabolic parameters.

Assessment of either spontaneously breathing or mechanically ventilated subjects.

Intended for testing patients above 15kg of weight or 6 years of age.

Available with "Low Flow" Turbine Flowmeter for canopy and mask tests and a single-use pneumotach for tests with mechanical ventilated patients.

REE by Canopy Dilution

Provided with an Adult Canopy hood, a paediatric version of canopy hood is available as an option.

Canopy blower is integrated in the device and it is easily controlled through software.

The software prompts an intuitive widget to help the operator in maintaining a stable CO_2 expired fraction (F_{eCO_2}) during dilution.

The Canopy veil is easy to mount and made in medical grade LDPE. It's a single-use item in order to avoid any possible cross contamination between subjects.

Cleaning the hood is easy and can be done with easily accessible solutions.

REE by Mask and Mouthpieces

REE tests can also be done by wearing multi-use silicone oro-nasal face masks (available in 5 sizes: 3 adult, 2 pediatric).

In addition to Canopy and Mask, users can also use mouthpiece with Antibacterial filters, together with a nose clip.

REE on mechanically ventilated patients

The ICU Kit is an optional module available for measuring REE in patients undergoing mechanically assisted ventilation in intensive care units.

Flow and Volume is measured with a single-use pneumotach flowmeter (Flow-REE), to be positioned in line between the endotracheal tube and the "Y" connector of the ventilator circuit.

All parts required during testing (Flowmeter, sampling line and HME filter) are single patient, with no need for cleaning and disinfecting after a test.

Quark RMR allows to assess ventilated patients up to $FiO_2 \leq 70\%$.

The patient setup makes Quark RMR completely independent from any type of ventilator in use.



REE by dilution with canopy hood



REE breath by breath by face mask and with mouthpiece and AB filter



REE on mechanically ventilated patients

Cardio Pulmonary Exercise Testing

On top of the standard features, the optional module for Cardio Pulmonary Exercise Testing (CPET) extends the possibility to perform full exercise protocols during exercise efforts.

Fast response analyzers provide accurate, reliable, breath-by-breath gas exchange data at any exercise intensity.

CPET made easy thanks to OMNIA, the new generation of COSMED software. The intuitive, beautiful, and innovative user interface brings complex CPET procedures to a new simpler stage.

CPET by breath by breath

BxB is the standard configuration of the CPET Module. It includes a "High-Flow Range" flow reader with 2 extra turbines.

Tests are conducted using ergonomic multi-use silicone oro-nasal face masks (available in 5 sizes: 3 adult, 2 pediatric) for comfortable testing in any condition.

Masks are also available with 2 inspiratory valves, to reduce inspiratory resistance and to prevent moisture accumulation especially at high intensity exercise.

CPET by Mixing Chamber

This optional module includes a physical mixing chamber (7 liters) with 2/way valve and adapters.


Ideal for gas exchange analysis when testing athletes ventilating at a frequency over 60 breaths per minute.

Simplified patient set-up, with turbine flowmeter placed at the exhalation port of the mixing chamber, avoids the use of the cumbersome conventional helmet.

The software provides flowmeter calibration specific for Mixing Chamber test to linearize response at its best.



Mixing chamber



Spirometry

Software module for performing FVC, SVC, MMV and Pre/Post Bronchial Provocation.

Real time acquisition and capture of Exercise Flow/Volume loops (EFVL) with comparison of resting FVC for evaluating ventilatory limitation.


Trial Selection and Quality Control in compliance with ERS/ATS guidelines.

Paediatric incentivations with user defined effort grade on both volume and flow.

Full compliance with "2005 ATS/ERS consensus" (Interpretation, QC, etc.).

GOLD COPD Interpretation on FVC PostBD.

Latest Global Lung Initiative (GLI) predicted (including Z-score).



FVC Test

Options and Accessories

Carts. Full range of carts either medical-graded with isolation transformer (available either with 230 or 120 VAC) or not electrified cart. Both 1 and 3-cylinder holder carts can be equipped with 1 or 2 monitors.

High FiO₂ kit. Gas exchange measurements using hypoxic and hyperoxic gas mixtures.

Ethanol burning Kit. The kit consists of a lamp, parts and connectors to be wired to the Quark RMR. Burning ethanol generates a predictable ratio of VO₂ and VCO₂ and it can be used to verify the Quark RMR accuracy of Respiratory Quotient measurement.

Wide selection of ergometers, available from COSMED, including treadmills, cycle-ergometers, arm-ergometers and recumbent bikes, suitable for any kind of tested subject.

Interfaced Devices

Diagnostic quality 12-lead Stress ECG (available either in wireless or patient cable configuration) with full disclosure and scroll back during test. High resolution ECG processing produces an exceptionally clear on-screen display and allows detailed, reliable analysis of ST segments and minimal arrhythmia changes. Available with Resting and Exercise ECG interpretation software.

Pulse Oximeter. High quality monitors (Nonin® technology) with a broad range of sensors (finger, earlobe or forehead/reflectance).



Ethanol burning kit



COSMED stress test ECGs (wireless or patient cable)



Pulse oximeter

Non-Invasive Blood Pressure (Suntech Tango). Cardiac stress blood pressure monitor specifically designed to overcome noise, motion and physical difficulties associated with cardiac stress and exercise testing

Cardiac Output (Physioflow Enduro) Portable, battery powered, non-invasive hemodynamic monitor for reliable and repeated cardiac output measurements during exercise.

Philips IntelliBridge compatibility. It allows to transfer REE parameters during testing (VO₂, VCO₂, RQ, REE, VE, Rf) through Philips monitoring solution directly to the Hospital Information System or to Philips IntelliSpace Critical Care (ICCA) systems.



Cardiac Output monitor (Physioflow)



Blood Pressure Monitor (Tango)

Data Management & Software

Quark RMR comes with **OMNIA Metabolic Module**, the new software designed by COSMED, compatible with the entire COSMED product range, OMNIA allows the user to operate different equipment in a single software environment.

Easy-to-use touch-screen graphic user interface with intuitive workflow and hierarchy.

Manage and display data and charts through standard (9 panel plot, etc.) or user defined Dashboards.

Select and define charts, data and widgets to define your preferred working environment.

Powerful chart creation (up to 4 Y axis and one X axis) with full control on settings.

Easy, quick and fully assisted calibration for high accuracy measurements, either for flowmeters (calibration and linearity check) or for gas sensors (zero, gain and delay).

Powerful post-test editing phase for data filtering, calculation of thresholds (AT, RCP), VO_2 max, EFVL, VE/VCO_2 slope, intercept and other parameters requested for interpretation.

Comprehensive interpretation tool automatically elaborates CPET tests and provides interpretation including text strings and numerical results based on latest scientific guidelines¹.

Built-in Rest and Exercise Protocol editor to design and save any type of protocol.



Wide list of Ergometers can be automatically controlled: (COSMED Bike/Treadmill, Ergoline, HPCosmos, Monark, Trackmaster) and with the optional Ergometer module (LODE, CSafe Treadmill, Cyclus 2, Technogym, Imbramed, Woodway and many others).

Export data in .pdf, .xml, and xls formats. Import data in .xml format.

SQL Database allowing virtually unlimited records and data safety.

Multi-users access rights management (Principal Investigator, Physician, Technician, Administrator...) with event logging.

Compatible with Win 7, 8, 8.1, 10 (32/64). Mac OS compatibility when installed in Virtual PC OS (Parallel, VMware).



Subjects/Visits/Tests database in OMNIA

Measured Parameters with OMNIA

Resting Energy Expenditure (REE)

Resting Energy Expenditure	REE, RMR (Kcal/day)	●
Respiratory Quotient	RQ	●
Substrates	%FAT, %PRO, %CHO	●

Cardio Pulmonary Exercise Test (CPET)

VO_2 max	VO_2 /Kg, RQ, Dyspnea@Max, HR@Max	○
Thresholds	AT, RCP	○
VE Response	VE/VCO_2 slope, VE/VCO_2 intercept, OUES	○
VO_2 /WR	VO_2 /WR slope, O_2 Pulse R2	○
EFVL	Flow/volume loop events	○

Spirometry

Forced Vital Capacity	FVC, FEV1, FEV1/FVC%, PEF	○
Slow Vital Capacity	VC, IVC, EVC	○
Maximum Voluntary Ventilation	MMV, MRf, MVT	○
Broncho-challenge		○

1 ATS/ACCP 2001, ESC 2009, EACPR/AHA 2012, AHA 2010

Networking

OMNIA Network allows to share a single database in either a small network (LAN) or a large network (WAN) environment.

OMNIA Network is based on a Client/Server architecture and allows to run different COSMED devices through simultaneous access of data and run tests via a virtually unlimited number of COSMED products.

The network license includes five clients (simultaneous access) and can be extended with the purchase of additional single licenses.

A user management system allows to define users (Physician, Technician, Administrator, etc.) and roles (which specific feature can a user access).

OMNIA can exchange data with Hospital Information Systems (HIS) via HL7, GDT and with a proprietary Protocol (OCP).

With the optional HL7 module (either standalone or network) OMNIA allows to get data from an HL7 worklist and send results back to Electronic Medical Records (EMR) and Hospital Information Systems (HIS).

Access and security compliant according to US HipAA, ISO 27799:2008, EU 95/46/CE and 2002/58/CE.

OMNIA Network runs on Windows Server 2008 (SP2, R2 SP1) and 2012.

Based on standard SQL database (Express or Standard) to store data securely.

COSMED
The Metabolic Company
Via dei Piani di Monte Savello, 37 - 00041 Albano Laziale - Rome
http://www.cosmed.com

Test date: **25/09/2015**
Report date: **04/11/2015**

Name: SUBJECT DEMO		ID: ---	Gender: Male	Age: 48	Weight (kg): 97.00	Height (cm): 195
Company: COSMED		D.O.B.: 04/03/1967	SSN: ---	BMI (kg/m ²): 25.5	Smoke: ---	Cig/Day: ---
Occupation: demo subject	Technician: Mr. Q	Physician: ---	Dr. House			
Ethnic: Caucasian	Room: ---	Set: ---	Wasserman extended (Exercise Testing (Clinical))			

Exercise (9P Panel) @ 11:26

Test Duration: 13:14		Test Position: Sitting	Reason for test: Exercise Capacity
Exercise Duration: 09:21	Ergometer: COSMED Bike	Reason for stopping test: Leg Fatigue	
Subject Type: Clinical	Protocol Name: 35 Watt Ramp	Test Purpose: Educational	
Test type: Maximal	Max Effort Confirmed: Yes		

W2 mL/min

296 3610 496

AK L

Interpretation:

At peak of exercise respiratory exchange ratio is 1.17, heart rate is 98% of predicted and maximum Rate of Perceived Exertion is Somewhat strong for dyspnoea and Strong for leg pain. Subject achieved maximal effort.

Maximum oxygen uptake is 3610 mL/min (117% of the predicted value) indicating a normal Aerobic Exercise Capacity.

Anaerobic threshold has been identified at 2478 mL/min and at 80% of the VO2max predicted (69% of measured VO2max), which is normal.

VE/VO2 slope is 22.2 and VE/CO2 slope is 26.2, which are both within normal range. This indicates a normal ventilatory efficiency.

Breathing Reserve at peak exercise is 39.4%, which is normal. This indicates an absence of exercise ventilatory limitation.

Variation of IC is -0.03 L. Dynamic hyperinflation is confirmed.

Heart rate reserve is 4 bpm, which is normal.

Heart rate recovery after 1 minute is 29 bpm, which is normal.

Oxygen pulse at peak exercise is 21.5 mL/beat, corresponding to 120% of predicted, which is normal.

VO2/WR slope is 9.46 mL/min/Watt. This indicates a normal aerobic work efficiency.

Maximal systolic blood pressure achieved during exercise is 130 mmHg, which is normal.

SpO2 during exercise did not significantly decrease.

Confirm Report

Customisable header and patient information

Test information

Editable interpretation string

Spirometry		Pre	% Pred	Normal	Class
FVC	L	6.07	97	> 4.32	
FEV1	L	4.80	98		
MVV	L/min	192.0			

Protocol		Meas.	Rest	Warm Up	AT	RC	Max	Normal	Class
Power	Watt	0	203	280	329	> 270			
Revolution	RPM	78	76	76	88				


Metabolic		Meas.	Rest	Warm Up	AT	RC	Max	Normal	Class
VO2	mL/min	646	897	2478	3112	3610	> 2595		Normal
VO2/Kg	mL/min/Kg	6.7	9.2	25.6	32.1	37.2	> 26.8		Normal
METS	---	1.9	2.6	7.3	9.2	10.6	> 7.6		Normal
RQ	---	0.78	0.85	0.93	1.06	1.17	> 1.10		Maximal

Ventilatory		Meas.	Rest	Warm Up	AT	RC	Max	Normal	Class
VE/VO2 slope	---	22.1					< 29.1		Normal
OLUES	mL/min/L/min	3870					> 3239		
VE	L/min	19.8	27.8	60.4	93.2	118.3	> 15.0		Normal
BR	%			68.6	51.5	38.4	> 15.0		Normal

Cardiovascular		Meas.	Rest	Warm Up	AT	RC	Max	Normal	Class
HR	bpm	76	88	132	151	168	> 155		Normal
HRR	bpm	4					< 15		Normal
HRR_1_minute	bpm	29					> 12		Normal
VO2/WR Slope	mL/min/Watt	8.92					> 8.40		Normal, Continual Rise
VO2/HR	mL/beat	8.4	10.2	18.8	20.6	21.5	> 14.4		Normal, Continual Rise Throughout Exercise
SBP	mmHg	0	0	130	130	< 210			Normal
DBP	mmHg	0	0	85	85	< 90			Normal

Gas Exchange		Meas.	Rest	Warm Up	AT	RC	Max	Normal	Class
VO2/AT	---						> 1236		

For each test, users can define parameters and columns to display



Select and edit graphs

t	VO2/Kg	VO2	VO2	R	VE	VE	Vt	VO2/HR	VE/VO2	VE/CO2	PetO2	PetCO2	SpO2	Power	HR	Phase
hh:mm:ss	mL/min/Kg	mL/min	mL/min	---	L/min	L/min	L/breath	mL/beat	---	---	mmHg	mmHg	---	Watt	bpm	---
00:03	7.1	369	389	0.78	9.2	13.1	0.998	4.5	24.9	31.7	105	38	97	0	82	Rest
00:08	7.1	369	287	0.78	9.1	11.2	0.812	4.4	24.5	31.5	106	38	97	0	84	Rest
00:13	6.5	339	258	0.76	8.2	13.7	0.600	4.0	24.3	31.9	105	38	97	0	84	Rest
00:16	7.7	403	390	0.75	9.7	16.1	0.983	4.9	25.5	29.8	105	39	97	0	86	Rest
00:20	9.2	476	363	0.76	11.1	14.6	0.764	5.4	23.4	30.6	104	39	97	0	88	Rest
00:24	8.7	454	349	0.77	10.6	16.2	0.652	5.1	23.3	30.3	106	37	97	0	86	Rest
00:29	16.9	800	734	0.83	23.1	10.7	2.168	10.1	26.5	31.5	111	36	98	0	87	Rest
00:35	7.3	380	343	0.90	10.7	11.3	0.948	4.5	28.2	31.3	110	37	98	0	85	Rest
00:39	9.4	491	425	0.87	13.3	15.5	0.870	5.9	27.5	31.8	111	37	98	0	84	Rest
00:43	7.9	412	325	0.86	11.8	13.9	0.818	4.9	27.5	31.8	110	37	98	0	84	Rest
00:48	0.0	468	407	0.87	12.1	12.0	1.013	5.4	25.9	29.8	108	39	98	0	86	Rest
00:52	8.1	422	356	0.84	11.5	14.5	0.789	4.9	27.2	32.2	109	37	98	0	86	Rest
00:57	8.9	463	392	0.85	11.7	13.1	0.883	5.8	25.5	29.8	108	39	98	0	86	Rest
01:04	7.4	382	327	0.86	9.6	8.3	1.159	4.4	25.1	29.4	107	39	98	0	86	Warmup
01:08	6.1	315	254	0.81	8.4	14.6	0.578	3.6	26.7	33.2	109	36	98	0	88	Warmup
01:14	5.1	248	221	0.84	6.9	16.5	0.658	3.8	26.2	31.8	106	38	98	0	87	Warmup
01:18	6.0	360	285	0.79	9.0	14.2	0.634	4.2	25.0	31.6	105	38	98	0	86	Warmup
01:22	11.1	577	453	0.79	14.4	17.3	0.823	6.7	25.0	31.8	105	38	98	0	86	Warmup

Customisable tabular data

Custom "CPET" printout report

Validation articles

REE

- Ashcraft C.M. et al. "A Test of Validity of a New Open-Circuit Indirect Calorimeter." *J Parenter Enteral Nutr.* 2014 Mar 10
- Sundström M et al. "Indirect calorimetry in mechanically ventilated patients. A systematic comparison of three instruments." *J. Clin Nutr.* 2013 Feb;32(1):118-21
- Blond E. et al. "A new indirect calorimeter is accurate and reliable for measuring basal energy expenditure, thermic effect of food and substrate oxidation in obese and health" *e-SPEN e-Journal of Clinical Nutrition and Metabolism* 6 (2011) e7ee15

CPET

- Gullstrand L., et al. *Validation of the Quark CPET Respiratory gas analyser in the BBB mode.* 2013; *Elite Sport Centre, Bosön*
- Gullstrand L., et al. *Validation of the Quark CPET Respiratory gas analyser (Mixing Chamber).* 2013; *Elite Sport Centre, Bosön*
- Nieman DC, et al. *Validity of COSMED's Quark CPET mixing chamber system in evaluating energy metabolism during aerobic exercise in healthy male adults.* *Res Sports Med.* 2013;21 (2):136-45
- More scientific studies on www.cosmed.com/bibliography

Technical Specifications

Product	Description
Quark RMR	Indirect Calorimetry Laboratory (Ref. C09074-01-99)
Standard packaging	Quark RMR unit, canopy hood (with pipes and backpack), turbine flowmeter w/ sampling line, calibration syringe (3 liters), HR monitor (receiver and transmitter), OMNIA PC software, adapters, pipes, cables, probes and user manual
Standard tests	
Indirect Calorimetry	Resting Energy Expenditure (REE, RMR), w/ face masks or mouthpieces. Respiratory Quotient (RQ) & Substrates Analysis
Optional tests	
ICU kit for vent Patients	Allowing Quark RMR measurement of REE in patients undergoing mechanically assisted ventilation.
Cardio Pulmonary Exercise Test (CPET)	Pulmonary Gas Exchange (VO_2 , VCO_2), VO_{2max} , Sub-max VO_2 , Thresholds (AT, RCP), EFVL, Heart Rate
Spirometry	Forced Vital Capacity (FVC) Pre/Post, Slow Vital Capacity (SVC) Pre/Post, Maximum Voluntary Ventilation (MVV), Broncho-challenge - Bronchial Dilator/Constrictor test
Flowmeters	
Turbine 0-18 (Standard)	Multiuse digital turbine for REE Mask/Canopy Test: Flow range 0-8 l/s; Accuracy $\pm 2\%$ or 20 ml/s (flow) $\pm 2\%$ or 100 ml/min (vent.); Resistance <0.7 cmH ₂ O l/s @ 3l/s; Ventilation range 0.04-50 l/min
Flow-REE (ICU Option)	Disposable PNT (Lilly) for REE ICU Test: Flow range 0-1,7 l/s; Accuracy $\pm 2\%$; Resistance $<2,35$ cmH ₂ O/l/s @ 1 l/s
Turbine 0-28 (CPET Option)	Multiuse digital turbine for CPET Mask/Mix test: Flow range 0.08-20 l/s; Accuracy $\pm 2\%$ or 20 ml/s (flow) $\pm 2\%$ or 200 ml/min (vent.); Resistance <0.6 cmH ₂ O /l/s @ 14l/s; Ventilation range 0.08-300 l/min
Gas Analyzers	
O ₂	Paramagnetic sensor. Range: 0-100% (Standard 0-30% - ICU 0-70% - or user defined); Accuracy: $\pm 0.1\%$; Response time: 120 ms
CO ₂	Nondispersive infrared sensor (NDIR). Range: 0-10%; Accuracy: $\pm 0.02\%$; Response time: 100 ms
Hardware	
Dimensions & weight	Unit: 17 x 30 x 45 cm/8 Kg Canopy: 32 x 50 x 30 cm/0.6 Kg
Interface ports	USB A-B, RS-232, HR-TTL, SpO2
Electrical requirements	100-240V $\pm 10\%$ 50/60 Hz
Internal emergency battery	12V; 1,2 Ah
Environmental conditions	Temperature 0-50 °C (32 - 122 °F); Barometer 400-800 mmHg; Humidity 0-100%
Software	
Available languages	Italian, English, Spanish, French, German, Portuguese, Greek, Dutch, Turkish, Russian, Chinese (Traditional & Simplified), Korean, Romanian, Polish, Czech, Norwegian, Hebrew
PC Configuration	I3 or higher processor speed. Compatible with Windows 7, 8, 8.1, 10 (32 or 64 bit). RAM 4GB (8GB recommended). HD with 4GB of free space (plus tools)
Safety & Quality Standards	
MDD (93/42 EEC); FDA 510(k); EN 60601-1 (safety) / EN 60601-1-2 (EMC)	

CE
0476



COSMED Srl

Via dei Piani di Monte Savello 37
Albano Laziale - Rome 00041, Italy

+39 (06) 931-5492 Phone
+39 (06) 931-4580 Fax

info@cosmed.com | cosmed.com

To know more:

