

# MC

series

**MC-1000**

1100 g × 0.0001 g

**MC-6100**

6100 g × 0.001 g

**MC-10K**


10.1 kg × 0.001 g

**MC-30K**

31 kg × 0.01 g

**MC-100KS**

101 kg × 0.1 g



**MC-10K**



**MC-6100**



**MC-100KS**



**New model  
added!**



# Revealed—Here Are What A&D Use To Manage Their Masses!

Are you finding it difficult to find high performance, quality mass comparators that meet your budget? Discover how A&D, an ISO17025 accredited calibration laboratory as well as being a weighing equipment manufacturer, has the solution for you!

## 1. Enhanced stabilizing filter

High resolution balances like the MC series can respond to even the slightest environmental disturbances. For smooth and accurate mass comparison, the MC series is equipped with a special stabilizing function that copes with the effect of drafts or vibrations.

## 2. Simplified calibration of the balance


Using the internal mass, you can ensure the balance accuracy with one key press when changes to the environment (such as temperature or location) take place. The internal mass value can be corrected as necessary.

## 3. Various options to minimize possible measurement errors

Even under a perfectly controlled weighing environment, there still remain errors that can arise from operation. The following tools are effective in managing such errors:

### **Auto-centering pan (available for MC-1000/6100/10K/30K)**

Corrects eccentric loading automatically by bringing the mass to the center.



MC Series  
with auto-centering pans



### **Remote controller (AD-8922A)**

Enables operations without directly touching and consequently shaking the balance.



### **Static eliminator (AD-1683)**

Instantly neutralizes charged objects such as the operator's clothes to ensure stable weighing.



# Recommended Measurement Ranges

Model	BM-20					BM-252				BM-500				MC-1000				MC-6100				MC-10K				MC-30K				MC-100KS	
	Class	E2	F1	F2	M1	M2	F1	F2	M1	M2	F1	F2	M1	M2	F1	F2	M1	M2	F1	F2	M1	M2	F1	F2	M1	M2	M1	M2			
1 mg																													1 mg		
2 mg																													2 mg		
5 mg																													5 mg		
10 mg																													10 mg		
20 mg																													20 mg		
50 mg																													50 mg		
100 mg																													100 mg		
200 mg																													200 mg		
500 mg																													500 mg		
1 g																													1 g		
2 g																													2 g		
5 g																													5 g		
10 g																													10 g		
20 g																													20 g		
50 g																													50 g		
100 g																													100 g		
200 g																													200 g		
500 g																													500 g		
1 kg																													1 kg		
2 kg																													2 kg		
5 kg																													5 kg		
10 kg																													10 kg		
20 kg																													20 kg		
50 kg																													50 kg		
100 kg																													100 kg		

\*Please refer to the BM series of micro analytical balances for the BM-20, BM-252 and BM-500.

## Wait, there is more! If you are in the industrial sector, how about using the MC series as **Precision Balances With Extended Resolution**

The MC series displays the weighing results with an extra decimal place not normally found in balances of the same capacity. Such extended resolution enables measurements of minute weight changes with heavy objects or containers. Example applications include:

### Management of Liquids

- Resist inks on flat panels (glass)
- Coating materials on silicon wafers
- Moisture contents of catalysts
- Solder pastes and sealing materials on electronic parts
- Resins used for LED lenses
- Surface materials on optical lenses
- Grease in miniature bearings for hard disks

### Management of Powders

- Medicines in pharmaceutical manufacturing – from full flow to minute flow
- Discharges (losses in weight) from feeders
- Wear volumes of metallic materials
- Cutting dust (contaminants) left on mechanical parts
- Particles caught by filters

### Management of Gases

- Fill ration of hydrogen gas
- Gasified components

### Production Line Weighing

- As precision check-weighing sensors with wider ranges



Optional



(Optional for MC-100KS)

MC-10K/30K/100KS

MC-1000/6100

## Specifications

**New!**

Models	MC-1000	MC-6100	MC-10K	MC-30K	MC-100KS
Weighing capacity	1100 g	6100 g	10.1 kg	31 kg	101 kg
Minimum weighing value	0.0001 g	0.001 g		0.01 g	0.1 g
Repeatability (standard deviation) *1	0.0005 g (for 1 kg) 0.0004 g (for 500 g)	0.004 g (for 5 kg) 0.0015 g (for 2 kg)	0.005 g (for 10 kg) 0.004 g (for 5 kg)	0.015 g (for 20 kg)	0.2 g (for 100 kg) 0.1 g (for 60 kg)
Linearity	±0.003 g	±0.03 g		±0.2 g	±2 g
Sensitivity drift (10 to 30°C / 50 to 86°F)	±2 ppm / °C			±3 ppm / °C	±6 ppm / °C
Accuracy immediately after internal calibration *2	±0.010 g	±0.15 g		±1.5 g	±10.0 g
Operating environment	5°C to 40°C (41 °F to 104 °F), 85%RH or less (no condensation)				
Display refresh rate	5 times/second or 10 times/second				
Display mode *3	Gram, kilogram (for MC-100KS only), ounce, pound, pound-ounce, troy ounce, metric carat, momme, pennyweight, grain, counting mode, percent mode, density mode, and a user-programmable unit				
Standard interface	RS-232C				
Weighing pan size	128 x 128 mm	165 x 165 mm	270 x 210 mm		386 x 346 mm
External dimensions	210 (W) x 317 (D) x 86 (H) mm		300 (W) x 355 (D) x 111 (H) mm		Display unit: 260 (W) x 164 (D) x 202 (H) mm Platform: 346 (W) x 443 (D) x 130 (H) mm Cable length: 3 m
Net weight	Approx. 4.6 kg	Approx. 5.1 kg	Approx. 9.3 kg		Approx. 18 kg
Power consumption	Approx. 11 VA (supplied to the AC adapter)				Approx. 15 VA (supplied to the AC adapter)

\* 1 Repeatability when the auto-centering pan or an automatic weight loader is used under good environmental conditions.

\* 2 Accuracy immediately after internal calibration performed under good environmental conditions. The values are for the weighing capacities.

\* 3 Either tael (Singapore/HK jewelry/Taiwan) or tola can be added upon request.

## Options

### For MC-1000/6100

- GX-02 \*4 Quick USB interface with cable
- GX-04 \*4 \*5 Comparator output with a buzzer/RS-232C/Current loop output
- GX-06 \*4 \*5 Analog output/Current loop output
- GX-12 Animal weighing pan
- GX-13 Density determination kit (for MC-1000 only)

\* 4 Only one of GX-02, GX-04, GX-06 or the standard RS-232C interface can be used.

\* 5 Factory-installed options

### For MC-10K/30K

- GX-04K \*6 Comparator output with a buzzer/RS-232C/Current loop output
- GX-06K \*6 Analog output/Current loop output
- GX-07K Waterproof RS-232C cable (25P-9P, 5 m)
- GXK-012 Animal weighing pan
- GXK-015 Carrying case

\* 6 Only one of GX-04K, GX-06K or the standard RS-232C can be used.  
The MC-10K/30K are not IP65 compliant with these options.

### For MC-100KS

- GP-04 \*7 Comparator output with a buzzer/RS-232C/Current loop output
- GP-06 \*7 Analog output/Current loop output
- GP-07 \*8 5-meter extension cable
- GP-12 Animal weighing pan
- GP-21 Underhook
- GP-22 AD-8121B printer support

\* 7 Only one of GP-04, GP-06 or the standard RS-232C can be used.

\* 8 Factory-installed option

## Accessories

- AD-1602-E2 \*9
  - AD-1602-F1 \*9
  - AD-1682
  - AD-1683
  - AD-1684
  - AD-1687
  - AD-1688
  - AD-1689
  - AD-1691
  - AD-8121B
  - AD-8526
  - AD-8527
  - AD-8920A
  - AD-8922A
  - AX-MC1000PAN
  - AX-MC6100PAN
  - AX-MC10K/30KPAN
  - AX-SW128
  - AX-USB-25P
- \* 9 Selectable from 50 g to 20 kg.
- Calibration masses, OIML class E2
  - Calibration masses, OIML class F1
  - Rechargeable battery unit
  - Static eliminator
  - Electrostatic fieldmeter
  - Weighing environment logger
  - Weighing data logger
  - Tweezers for calibration mass
  - Weighing environment analyzer
  - Compact printer
  - Serial/Ethernet converter
  - Quick USB adapter
  - Remote display
  - Remote controller
  - Auto-centering pan (MC-1000)
  - Auto-centering pan (MC-6100)
  - Auto-centering pan (MC-10K/30K)
  - Foot switch
  - Serial/USB converter



Auto-centering pan



AD-1602

**AND**  
A&D Australasia  
"Excellence in Measurement"

ACN 26 007 556 809

Head Office  
32 Dew Street  
THEBARTON South Australia 5031  
Telephone (08) 8301 8100  
Facsimile (08) 8352 7409

Victorian Office  
39 Bakehouse Road  
KENSINGTON Victoria 3031  
Telephone (03) 9371 1555  
Facsimile (03) 9372 1193

New South Wales Office  
Unit 4, 14 Abbott Road  
SEVEN HILLS New South Wales 2147  
Telephone (02) 9674 5466  
Facsimile (02) 9674 2544

sales@andweighing.com.au  
www.andweighing.com.au

Your Dealer



**ASIS**  
SCIENTIFIC

W asisscientific.com.au  
E sales@asisscientific.com.au  
P 08 8340 8444  
F 08 8340 8377  
A 474 Port Road, Hindmarsh, SA, 5007