

# touch<sup>bio</sup>

TouchBio RSV, FLU A/B & Covid-19 Rapid Antigen Combo Test (Nasal)

REF: VMD71 EN

An Antigen Rapid Test for the detection of SARS-CoV-2, Flu A/B and RSV in nasal swab. **For Self-Testing use.**

**IMPORTANT** Read the instructions carefully before taking the test.


Australia Sponsor & Distributor  
Touch Biotechnology Pty Ltd

Customer Support Number: 1300 166 282  
Hours: 9am-7pm (AEST), or 9am-3pm (AEDT), 7 days per week  
Website: [www.touchaustralia.com.au](http://www.touchaustralia.com.au)  
Email: [touch@touchaustralia.com.au](mailto:touch@touchaustralia.com.au)  
Address: 119 Willoughby Road, Crows Nest, NSW 2065

Watch "How to Use Video"  
Scan the QR code for information on how to use the test.



## COMPONENTS PROVIDED



### STEP-1 Wash your hands

Wash or clean your hands and make sure they are dry before starting the test.

After washing your hands, open the box, and check the components before use.

### STEP-2 Read Instructions for use

Read instructions for use carefully before using the test.

### STEP-3 Place the buffer tube into the holder

Carefully place extraction tube into tube holder or tube stand.

Open the lid.

**IMPORTANT** DO NOT DRINK the extraction buffer liquid. If you accidentally drink it immediately consult your healthcare professionals.

DO NOT SPILL any of the extraction buffer liquid. If you spill it, sterilize the area, and repeat the test by using new sampling swab and extraction solution tube.

### STEP-4 Take the sterilized swab

Pull open the swab packaging at the marked point and remove the swab.

**IMPORTANT** DO NOT TOUCH soft head of swab. DO NOT OPEN the swab until you are going to use it immediately.

### STEP-5 Sample Collection

Tilt your head back slightly

Gently insert the swab about 2cm into the **left nostril**. At least with the entire soft swab.

Rotate 5 times

Gently rotate the swab at least 5 times against the nasal wall.

Do the same for the **right nostril**. Gently insert the swab about 2cm. At least with the entire soft swab.

Gently rotate the swab at least 5 times against the nasal wall.

Remove the swab from the second nostril.

**IMPORTANT** IF YOU FEEL DISCOMFORT, STOP IMMEDIATELY. If the swab stick breaks during the sample collection, please use a new swab. Do not insert the swab deeper if you feel strong resistance or pain.

### STEP-6 Insert the swab

Insert the sampled swab into the extraction buffer tube, and dip the tip into the tube.

Rotate the swab tip 10 times along the inner wall of the buffer tube.

And squeeze the tip of the swab 5 times along the inner wall of the tube to keep as much liquid in the bottle as possible.

## TEST PROCEDURE

### STEP-7 Take out the swab

Remove the swab from the tube by squeezing the sides of the tube to release the liquid from the swab.

**IMPORTANT** If squeezing of tube is not done correctly, sample swab absorbs much more liquid from the extraction buffer and that will yield wrong results.

Discard the swab in the biohazard specimen bag.

### STEP-8 Close and Mix the tube

Screw on and tighten the lid on the extraction tube.

And then shake the extraction tube vigorously to mix the specimen and the sample extraction buffer.

**IMPORTANT** Ensure the lid is screwed on properly. Do not spill any of the sample extraction liquid.

### STEP-9 Take out the cassette

Open the foil pouch and take out the test cassette.

Place it on a flat and clean surface.

Perform the test within 15 minutes after the foil pouch is opened.

### STEP-10 Unscrew the white cap

Unscrew the white cap of the tube. This allows drop-wise dispensing of the liquid.

### STEP-11 Test Operation

Add 3 drops of the extraction buffer tube to the FLU A/B sample well marked "S" on the test cassette.

Do the same for the SARS-CoV-2 sample well marked "S", add 3 drops.

Do the same for the RSV sample well marked "S", add 3 drops.

Ensure to add at least 3 drops of the liquid from the specimen tube into the each sample well. If adding less than 3 drops, that will yield wrong result.

### STEP-12 Wait for result

Set timer and wait for 15 minutes.

Read the result at 15-20 minutes.

**IMPORTANT** DO NOT READ the result beforehand or after 20 minutes, even if a line has already appeared at the region "C"

### STEP-13 Read your results

To read your test results, please go to the interpretation of the results section provided below.

### STEP-14 Disposal

Please dispose all parts of the test kits and place them in the biohazard bag that can be disposed in the household waste or rubbish bin. If there are local regulations, please follow them.

### STEP-15 Wash your hands


Wash your hands thoroughly after test completion.

Watch "how to use" video:  
[touchaustralia.com.au/pages/ifu-covid-flu-rsv](http://touchaustralia.com.au/pages/ifu-covid-flu-rsv)

## INTERPRETATION OF THE RESULTS

EACH TEST WINDOW (FLU A/B, SARS-CoV-2 AND RSV ) MUST BE READ INDEPENDENTLY FROM EACH OTHER.


### + POSITIVE RESULTS



THE SHADE OF LINES MAY VARY, BUT EVEN IF A FAINT/WEAK LINE APPEARS, IT SHOULD BE CONSIDERED POSITIVE.

If you have a POSITIVE result, follow the guidance from your local State or Territory Health Department for guidance on confirmation testing if necessary, and if unwell seek medical assistance for SARS-CoV-2 and individuals with a positive result or who are unwell must consult a medical practitioner for follow-up clinical care for Influenza and RSV.


### - NEGATIVE RESULTS



ONLY RED LINES APPEAR IN THE CONTROL REGIONS (C), AND NO LINE IN THE REGION (A), (B), AND (T).

The negative result indicates that there are no Flu A/B, RSV and SARS-CoV-2 particles in the sample or the number of viral particles is below the detectable range. Even if you get a negative result, you still need to follow all public health advice on limiting the spread Covid-19, Flu A/B and RSV. If symptoms persist, repeat testing and consult a medical practitioner for follow-up clinical care.

### ? INVALID RESULTS



NO RED LINE APPEARS IN THE CONTROL REGION (C) FOR FLU A/B. The test is invalid even if there is a line on the region (A), (B) or (A) and (B).

NO RED LINE APPEARS IN THE CONTROL REGION (C) FOR SARS-CoV-2. The test is invalid even if there is a line on the region (T).

NO RED LINE APPEARS IN THE CONTROL REGION (C) FOR RSV. The test is invalid even if there is a line on the region (T).

WHEN NO RED LINE APPEARS IN THE CONTROL REGION (C) FOR ALL THREE TESTS, THE TEST IS INVALID FOR ALL THREE FLU A/B, SARS-CoV-2 AND RSV.

Insufficient sample volume or incorrect procedural techniques are the most likely reasons for the control line (C) failure. Review the test procedure and repeat the test using a new test device. If invalid result continues after repeating, please contact Touch Biotechnology on the provided contact number or email for assistance.

