Optiflow™+ Nasal Cannula

REF MYOPT9**SMALL**

REF MYOPT9**MEDIUM**

REF MYOPT9**LARGE**


Thank you for choosing F&P myAirvo™

Please read this entire guide before you first use your nasal cannula kit. Following the steps in this guide, and referring to the images and diagrams will assist you with daily use of your F&P myAirvo device and its consumables.

LANGUAGES

PAGE	LANGUAGE	PAGE	LANGUAGE
1	English (en)	14	Türkçe (tr)
2	Italiano (it)	15	العربية (ar)
3	Deutsch (de)	16	עברית (he)
4	Français (fr)	17	繁體中文 (zht)
5	Español (es)	18	Polski (pl)
6	ภาษาไทย (th)	19	Română (ro)
7	Português brasileiro (pt-br)	20	Ελληνικά (el)
8	Português (pt)	21	Slovenčina (sk)
9	Nederlands (nl)	22	Čeština (cs)
10	Dansk (da)	23	한국어 (ko)
11	Suomi (fi)	24	Русский (ru)
12	Norsk (no)	25	日本語 (ja)
13	Svenska (sv)		

Single patient multiple use

CE 0123 Conforms with medical device directive 93/42/EEC

Rx only Prescription only

F&P, Optiflow and myAirvo are trademarks of Fisher & Paykel Healthcare Limited.

For patent information, see www.fphcare.com/ip

Manufacturer

Fisher & Paykel Healthcare Ltd,
15 Maurice Paykel Place,
East Tamaki, Auckland 2013

PO Box 14 348 Panmure,
Auckland 1741,
New Zealand

Tel: +64 9 574 0100
Fax: +64 9 574 0158
Email: info@fphcare.co.nz
Web: www.fphcare.com

Importer/Distributor

Australia (Sponsor)
Fisher & Paykel Healthcare Pty Ltd,
19-31 King Street, Nunawading,
Melbourne, Victoria 3131

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Tel: +90 312 354 34 12
Fax: +90 312 354 31 01

UK

Fisher & Paykel Healthcare Ltd,
Unit 16, Cordwallis Park, Clivemont Road,
Maidenhead, Berkshire SL6 7BU, UK

Tel: 0800 132 189
Fax: +44 1628 626 146

USA/Canada

Tel: 1800 446 3908
or +1 949 453 4000
Fax: +1 949 453 4001

myAirvo Optiflow+ Nasal Cannula

REF MYOPT9SMALL

myAirvo Small Nasal Cannula 2-Pack

REF MYOPT9MEDIUM

myAirvo Medium Nasal Cannula 2-Pack

REF MYOPT9LARGE

myAirvo Large Nasal Cannula 2-Pack

INDICATIONS FOR USE

Nasal cannula patient interface for delivery of humidified respiratory gases.

For use with the myAirvo Series humidifiers.

SETUP

- myAirvo Series humidifiers with MYAIRVOKIT1/ MYAIRSPIRAL/900PT560/900PT560E series tubes or kits.
- Flow range: MYOPT9SMALL 10 – 50 L/min
MYOPT9MEDIUM 10 – 60 L/min
MYOPT9LARGE 10 – 60 L/min

HOME AND LONG-TERM CARE FACILITY USE

This product is intended to be used for a maximum of 30 days providing daily and weekly cleaning instructions are followed (see the myAirvo Series User Manual).

i SETUP INSTRUCTIONS

Step 1

Select appropriate size. Prongs must not create a seal in the nares. A clear gap must be visible around each prong.

Step 2

Adjust headstrap to fit. Do not over-tighten.

Step 3

- Ensure headstrap clip is attached to the headstrap, to prevent cannula from being pulled out of the nares.
- Cannula can become unattached if not used with the headstrap clip.

Step 4

Attach breathing tube clip to a secure location (e.g. clothing or bedding) to prevent cannula from pulling off face.

Warnings:

- Appropriate patient monitoring must be used at all times. Failure to monitor the patient may result in loss of therapy, serious injury or death.
- This device is not intended for life support. Do not use on patients who cannot tolerate a brief interruption of therapy, to avoid serious injury or death.
- Nasal delivery of respiratory gases generates flow-dependent positive airway pressure (PAP). This must be taken into account where PAP could have adverse effects on a patient.
- Do not use near a naked flame, to avoid fires.
- Do not crush or stretch tube, to prevent loss of therapy.
- For single patient use only. Reuse may result in transmission of infectious substances. Attempting to reprocess will result in degradation of materials and render the product defective.

Cautions:

- Failure to use the set-up described above can compromise performance and affect patient safety.
- Before connecting the interface, check for adequate gas flow and ensure that the system has warmed up.
- Check for condensate regularly. Drain as required.
- Do not use if individual interface packaging is not sealed.

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