



OMS FRACTURE WALKER

The Fracture Walker is designed to immobilise the foot and ankle to decrease pain and swelling, promoting faster healing following injuries or surgical procedures.



Designed to stabilise the ankle at 90° to promote healing!

KEY FEATURES

- Rocker sole for added comfort and ease of ambulation, as well as a non-slip sole to protect from slips, providing additional security in each step
- Uses Velcro closures, they can accommodate changes in swelling following an injury, providing a more custom fit
- Patented metal interlocker for uprights, provide increased strength and durability
- Breathable moisture wicking liner and EVA foot bed increase comfort

BENEFITS

- To immobilise the foot and ankle following stable foot or ankle fracture, severe ankle sprain, post-surgery application, soft tissue injuries and more.

SIZE	Women (US)	Men(US)	Footbed length (cm)
Small	6 – 8	4.5 – 7	25
Medium	8.5 – 11.5	7.5 – 10.5	28
Large	12 – 13.5	11 – 12.5	30

The fracture boot heel is reasonably thick, and if you walk around lop-sided for several weeks it's common to develop some hip and back pain.

ITEM CODE	
LL7401	Fracture Walker (Moonboot) Length 43cm - Small
LL7402	Fracture Walker (Moonboot) Length 43cm - Medium
LL7403	Fracture Walker (Moonboot) Length 43cm - Large

