



Enbio^s + Enbio PRO World's fastest autoclaves. Class B.



Another dimension. Same performance.




The new definition of sterilization speed.

Time is a relative value. It flows differently under different circumstances. But when you need to have your tools ready, you need them, well, FAST. While other autoclaves might test your patience,



Typical autoclave





Class B Vacuum Pouches+drying

Enbio^s. Speedmaster.


Enbio S is compact to fit even a cramped space. Unbelievably short process times will save your procedure in emergency, or allow sterilization of the instruments as you go, between patients.





Wrapped

Unwrapped



Unit weight 15 kg

Capacity

8 Max. packages amount (circa) Max. recommended batch weight 0,5 kg Process chamber volume 2,71

Pressure & temperature

Working pressure 2,1 Bar Max. pressure 2,3 Bar Max. process tempreature 137°C

Dimensions

Width 25 cm Height 16 cm Depth 56 cm

Electrical

Power supply Installed power Max. power consumption 15 A

Class B

Vacuum Pouches+ drying



Medical Certified

Medical device class II b EN-13060 93/42/EEC Directive



reddot design award winner

230V/50Hz 2250 W

Misc

Ingress protection level IP20 Data archiving Water conductivity*


USB stick < 15 µS/cm

* - demineralized water only

Enbio PRO. Small foot. Big space.

Enbio PRO is an ideal solution as a basic workhorse in a small practice. Superfast sterilization and room to hold most of the instruments, including popular implant casettes.





Unit weight 20 kg

Capacity

16 Max. packages amount (circa) Max. recommended batch weight 1,0 kg Process chamber volume 5,31

Pressure & temperature

Working pressure 2,1 Bar Max. pressure 2,3 Bar Max. process tempreature 137°C

Dimensions

Width 25 cm Height 20 cm Depth 56 cm

Electrical

Power supply Installed power Max. power consumption 15 A

Class B

Vacuum Pouches+ drying



Medical Certified

Medical device class II b EN-13060 93/42/EEC Directive





230V/50Hz 2250 W

Misc

Ingress protection Data archiving Water conductivity* IP20 USB stick < 15 µS/cm

* - demineralized water only

Implantology

enalo

Enbio PRO accommodates all popular implant cassettes that fit in the 29×19×8 space.



Enbio user manual.



1. Plug in 2. Choose program 3. Start working

Installation and operation of Enbio autoclaves does not require any special training.

Over 20 000 autoclaves sold in 5 years.

In just 5 years, Enbio's worldwide community has grown to over 20,000 users, who sterilize instruments in our autoclaves, for the complete safety of their customers and patients.





Designed in Switzerland

GLG

Worldwide presence.

OLG

 \bigcirc

515151616

while a main a second

0

0

Enbio autoclaves are manufactured in the EU and distributed globally through logistics centers covering Europe, North America and the Middle East.



enoio.





Enbio Group AG Eichengasse 3, 4702 Oensingen Switzerland

enbio.com