



# COVERS WHEELCHAIR RELATED








POLYMAILLE®

- Coating of polyurethane on Jersey polyester
- Integral cover over 6 sides
- Non-slip base
- Adapted for home care and medical institutions use
- Hypoallergenic
- Anti-bacterial, anti-germ, anti-fungal
- bi-extensible
- Hyper-breathable (liquid repellent material and permeable to steam), adapted for incontinence. Reduce maceration effects through easing the skin respiration.
- Zip on 3 sides (except P311C)
- Integrated handle



POLYMAILLE® HD

- Coating of polyurethane on Jersey polyamide
- Integral cover over 6 sides
- Adapted for heavy care and medical institutions use
- Hypoallergenic
- Anti-bacterial, anti-germ, anti-fungal
- bi-extensible
- Hyper-breathable (liquid repellent material and permeable to steam), adapted for incontinence. Reduce maceration effects through easing the skin respiration.
- Welded and sewed
- Zip on 2 sides (except P311C)
- With screen printing

	TYPE DE PROTECTION	
	 POLYMAILLE®	 POLYMAILLE® HD
<b>Composition</b>	Coating of Polyurethane on Jersey polyester	Coating of Polyurethane on Jersey polyamide with printed instructions
<b>Version</b>	Integral cover	Integral cover
<b>MATERIAL CHARACTERISTICS</b>		
<b>Anti-allergenic</b>	✓	✓
<b>Antibacterial</b>	✓	✓
<b>Antimicrobial</b>	✓	✓
<b>Antifungal</b>	✓	✓
<b>Bi-extensibility</b>	+	++
<b>Respirability</b>	+++	+++
<b>TAILORING CHARACTERISTICS (Cushions only)</b>		
<b>Kind of assembly</b>	Sewing	Welded seams
<b>Anti-slip basis</b>	✓	
<b>Opening</b>	3 sides (Except P311C: lower face)	2 sides
<b>Integrated handle</b>	✓	
<b>Serigraphy</b>		Printed washing and disinfection instructions
<b>MAINTENANCE &amp; DISINFECTION</b>		
<b>Wash</b>		
<b>Disinfection</b>	 <b>Disinfect using cold-application spray products (Quaternary ammonium - glutaraldéhydes...)</b>	