Hydra-dermabrasion

Facial EMS + Microcurrent



Ace of Face



CE F©





Hydra-dermabrasion

Beginning of skin care and build skin immunity

Hydra-dermabrasion is the first step in skin care

No matter how great the skin care products you apply, acne and inflammation can appear on the skin. The reason for this is the **Demodex folliculorum** and **Body waste**. Demodex follicullorum is an enemy of the skin that destroys the skin balance by destroying collagen as well as producing acne and causing inflammation. Leaving Demodex folliculorum and body waste in your skin is a shortcut to destroying the underlying skin balance.

Hydra-dermabrasion removes skin waste and at the same time, penetrates the skin with a solution suitable for all skin types. It smoothes the skin texture, and makes you experience a greater effect when undergoing other skincare procedures.



How EVE works

As the cleansing solution comes out from the tip, it rotates with a string radial principle, and by sucking waste materials, it is possible to remove waste from the pores without scratching the skin.



Step 1 Cleansing

Removing demodex folliculorum and skin waste.



Step 2 Extractions

Painless, automated extractions use vacuum suction to draw out impurities.



Step 3 Peeling

AHA & PHA & LHA & Enzyme helps the skin loosen dirt and debris from pores without irritation.



Step 4 Hydration

Solution comes out from the tip to nourish and protect the skin.

Features

Solution efficiency

Increased solution efficiency with the specially developed tips that minimises the disadvantages of the existing vortex tips. It doesn't break down easily.

Powerful suction & adjustable levels

Minimises skin irritation and cleanses the skin using a strong vacuum and penetrating power. Adjust the level of the vacuum power that is suitable to the user.

Ergonomic hand piece

Developed with an ergonomic design. It allows for a more comfortable treatment.

Ease of use

- LCD touch screen
- LED guide on tips
- Easy to plug and unplug solutions and sediment bottles.

Available for all Skin types

Cleansing solution that contains LHA (4th generation), a newly developed ingredient from AHA, BHA and PHA, and can be used on all skin types.

What you can expect

- Cleansing
- Peeling
- Hydration
- Brightening
- Exfoliation (blackhead, whitehead, demodex folliculorum)
- Toning
- Acne
- Sebum control
- Skin Immunity

Various Tips

Large Small Cleaning Diamond











Scalp

Cleaning

Diamond



Microcurrent / EMS

The beginning and end of Anti-aging

Microcurrent therapy is where Dr. Erwin Neher and Dr. Bert Sakmann won the 1991 Nobel Prize in Physiology and proved that there is a possibility of recovery when microcurrent is applied to damaged tissue. EVE has developed a microcurrent technique based on this theory. We pursue fundamental skin beauty.

How EVE works

Microcurrent electrotherapy technology increases the production of adenosine triphosphate (ATP) and protein, which plays an important role in cellular energy metabolism, using a low current of μA, which is hardly felt by the human body, thereby promoting skin cell metabolism and lymphatic massage. Pain relief and skin elasticity induction technology.



Features

- Non-invasive, painless non-surgical way to reduce wrinkles and tighten skin in a safe, trustworthy and effective way.
- Quick & effortless easy to learn, simple to operate, fast in delivering results.
- NO down time
- Immediate & long lasting effects Eve provides fast and continuous anti-aging results.
- · Variations of pre-set protocol programs Easy to use even for beginners as the program is set for the purpose of using the device for skin problems such as anti-aging (by age), acne, etc.

What you can expect

- Lifting
- Enhance definition
- Tighten and firm skin
- Stimulate collagen growth
- Drained out lymphatics
- Reduces puffiness
- Improves muscle tone
- Increases cellular activity
- Improved oxygenation of tissue and skin
- Pain relief
- Diminished pigmentation

Various hand pieces

The micro size of the ball intensifies the current and increases effectiveness.

To lift broader areas such as the cheeckbones, periorbital wrinkles and jaw line.

Comb type

To energise hair follicles, promote blood circulation, diminish hair loss, reduce dryness and dandruff.

Gloves

Ideal way to provide a more detail oriented microcurrent treatment.

Electrode patch (Option)

Easy, hands-free treatment targeting the fine lines and wrinkles.















BENEFIT of QUQL

• Storage space



Solution

- New ingredient LHA & Enzyme included

Sensor

- Preventing overflow of sediment bottle
- Preventing malfunction

• Easy to plug and unplug

- No more wasting time to plug solution and sediment bottle.



Detachable hose

- Easy maintenance
- Preventing leakage and break down



Various Applicators



LED GUIDE

- BLUE LED LIGHTS (UV) : Targets and disrupts propionibacterium acnes, improves oily and congested skin
- WHITE LED LIGHTS : Guidelines for users



• GUI (optimising your treatment to the fullest)

- The advanced user interface is not only user friendly, but easy to adopt



• EASY Pre-set protocol programs

- Total of 10 different pre-set protocols that allow you to easily provide customised treatment based on specific skin conditions or concerns in a controlled time setting.



Pre-set protocol programs

EVE provides a skin enhancement program just for you. Enjoy an explicit anti-aging effect with the EVE exclusive programs

Multi-Step Programs

- Anti-Aging group by age
- Facial Acne / seborrhea
- Rosacea
- Restoration post-destruction

Single-step programs

• Balancing Skin Metabolism

Balancing the mitochondrial membrane, thus improving skin cell metabolism

Alleviating Pain

Alleviating pain from facial myalgias by electrical nerve stimulation

Lymphatic Drainage

Improving blood and lymph flow removing waste from the skin cells

Myo-stimulation

Stimulating facial muscles (lifting of sagging skin) / relaxing muscles in hyper tonus

EMS Programs

- Complete comprehensive anti-aging care through EMS Special session.
- Stimulate and tone muscle to tighten the skin.
- Focus on deep skin areas where skin elasticity needs to be improved.





"COMFORT YOU, AND WILL GO UP WITH YOU."

Specifications

Size	20 x 40 x 30 cm (Device)
(Net)	36 x 40 x 105 cm (Cart)
Size	35 x 55 x 50 cm (Device)
(Packaging)	48 x 18 x 76 cm (Trolley)
Weight	13 kg (Device)
(Packaging)	8.5 kg (Cart)
Suction power (Hydra)	650mmHg
Current (Microcurrent)	2~500 UA
Power	100~240V, 50/60Hz,
Consumption	Max 100VA (24V/ 5A)

Contact



- (07) 5577 9881 | 0414 970 019
- mike@beautyequipment.com.au
- www.mmmedicalaesthetics.com.au