

# Clearway 2

by Breas



## ADVANCED TITRATION TECHNOLOGY

- Unique to the Clearway 2 – TreatRepeat®
- Intuitive touchscreen user interface
- Hand-held remote manual switch with start/stop function



## MULTIPLE MODE OPTIONS

- Programmable automatic modes
- Dedicated IPPB mode
- Manual mode



## ENHANCED PATIENT SYNCHRONY

- TreatRepeat® to help improve accuracy and speed of titration and set-up
- iTrigger® – new flow triggering technology
- Synchrony Beep® to aid patient-device synchrony



# Clearway 2

by Breas



## Versatile Airway Clearance

### Breas Medical

Founded in 1991 in Gothenburg, Sweden, Breas has become one of the global leaders in Home Mechanical Ventilation, Airway Clearance and Sleep therapy. Our vision is to provide Better breathing and quality of life to respiratory patients around the world and to help patients transition from Hospital to Home.

Breas is headquartered in Sweden and has subsidiaries in the UK, Germany, Spain, USA and China. Together with network of highly specialized distributors, we serve over 40 countries worldwide. The Breas brands Vivo, NIPPY, Clearway and Z1 are widely recognized and appreciated for their performance, ease of use and attractive design.

NIPPY and Clearway are trademarks of Breas Medical Ltd. © 2020 Breas Medical – All rights reserved. Breas Medical reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your Breas representative for the most current information. Breas and the Breas logo are trademarks of Breas Medical AB. COR-SAL-MAR-000476-Rev0

BREAS.COM

Unit A2, The Bridge Business Centre  
Timothy's Bridge Road  
Stratford-Upon-Avon Warwickshire CV37 9HW. UK

For orders and general enquiries  
orders@nippyventilator.com  
Tel: +44 (0)1789 293 460

**BREAS**







## Clearway 2 – Up close



## Versatile Airway Clearance

Introducing the new Clearway 2 with TreatRepeat®. The second-generation MI-E device by Breas incorporates the latest technology to bring effective and accurate treatment for patients requiring airway clearance and cough assistance.

### Advanced Titration Technology

TreatRepeat function allows clinicians to manually deliver different numbers of insufflations and exsufflations and then save the treatment they have just delivered to the device to be repeated automatically. Helping to improve accuracy and speed of titration and set-up.

### Multiple mode options

Up to four treatment profiles can be saved to the device. The settings used in the manual and programmed modes are also saved to the device. This allows for up to 7 MI-E prescriptions to be saved. The treatment settings in the IPPB and NIV modes are also saved.

The programmable mode enables the delivery of multiple insufflations as part of a cough cycle.

### Enhanced Patient Synchrony

TreatRepeat helps to improve patient device synchrony through optimised manual titration of in/exsufflation and pause delivery times, which are then saved in to a profile.

The new iTrigger on the Clearway 2 provides energy triggering technology based on the class leading eSync feature available in Breas ventilators, to help improve synchronisation.

The Synchrony Beep provides an audible alert to patients just before the device cycles to exsufflation, in automatic MI-E modes, allowing even patients who do not have visual contact with the Clearway 2 to synchronise with the device during MI-E therapy.

For more information or to request a trial contact your local distributor or representative.

