### Sit-u-Up

M BILITY CARE

Since 1992

Phone: (03) 9568 8383

39 Cambro Road, Clayton VIC 3168

info@mobilitycare.net.au

info@mobilitycare.net.au www.mobilitycare.net.au



The Sit-u-Up elevates a person smoothly

from lying down in bed to sitting up.



### **Technical Specifications**

Max user weight	150 kg / 24 stone
Min height	Flat
Total length	75 cm / 29 ½"
Width	70 cm / 27 ½ "
Total Product weight	4.8 kg / 10 ½ lbs
Back/headrest length	73 cm / 28 ¾"
Backrest angle adjustment	Flat to 65 degrees
Airflo Compressor	Mk3

#### **Product code: MPCA130200**

## Easy to fit

- Sits on top of the mattress
- With no fixtures or fittings, it can be used on any size bed
- One airline connects the Airflo to the Pillowlift. Press one button to inflate

# Improves wellbeing

- Stable and comfortable, a person can choose to be lifted on their back or side
- The Sit-u-Up sits on top of the mattress enabling it to be fitted to any bed and causes no disturbance to a second person in double bed
- Minimises moving and handling risks so reducing the risk of musculoskeletal injury for the carer
- Its portable and lightweight which means it can be rolled up for travelling

#### Safe to use

• Includes a structured backrest ideal for those who need extra support

backrest gives extra support

- Stable and compact
- Quilt provides comfort and support
- Inflation and deflation can be stopped at any time, allowing freedom to choose the most comfortable position
- Anti-slip pillow holder

## Simple to maintain

- Lightweight and portable, the Sit-u-Up separates into two parts
- For a routine clean, sponge clean with a detergent solution and dry thoroughly
- The covers can be machined washed
- Spare covers can be purchased

