



DERMALIFT 5D HIFU


Skin Tightening - Face Lifting - Body Slimming - Vaginal Tightening

5D HIFU

HIFU Machine - Multi-Functional




How it works




HIFU-powered transducers deliver concentrated waves of ultrasound into multiple layers of the skin and body tissues without making direct contact with surrounding areas.

Transferring stabilized temperatures of 65~75°C into the targeted zones, the Dermalift aims to induce a process of neo collagenesis for facial tightening and body contouring purposes.

FOR LIFTING & TIGHTENING




1 The dermal layer of the skin naturally synthesizes less collagen over the course of time.



3 The dermal layer of the skin naturally synthesizes less collagen over the course of time.



2 Focused ultrasound waves deliver energy directly within the skin layer causing a shrinking effect.



4 Through the regeneration of collagen, wrinkles are eliminated in the process of both firming and tightening the skin.


FOR SLIMMING & CONTOURING



1 The dermal layer of the skin naturally synthesizes less collagen over the course of time.



3 As a natural healing response of the body the cells regenerate joining together to fill the gaps.




2 The SMASFOCUS V energy destroys the fat cells that have been targeted.




4 The layer of fat naturally exits the body delivering a smooth, contoured body line.

FOR VAGINAL TIGHTENING



Introduction





Vmax single dot handle



1.5mm, 3.0mm, 4.5mm



2D hifu handle



1.5mm, 3.0mm, 4.5mm



3D hifu handle



1.5mm, 3.0mm, 4.5mm



6.0mm, 8.0mm, 10.0mm, 13.0mm, 16.0mm



Liposonix handle



8.0mm, 13mm



Vaginal handle



3.0mm, 4.5mm

Eyebrows, Wrinkles around eyes,
Wrinkles around month

Check, Forehead, Double chin,
Neck wrinkles, V-line forming

Check, Forehead, Double chin,
Neck wrinkles, V-line forming

Abdomen, Inner thighs, Outer
thighs, Hips & Buttocks, Inner
arms, Breast

Body slimming and fat reduction

Vaginal tightening

Forehead Wrinkle

Wrinkle Lifting


Cheek

Neck Wrinkles




Easy to upgrade

Each machine is pre-loaded with a software system with all handle functions added for easy operator use. When you need to add another handpiece to the machine, the software system can automatically identify the new handle to upgrade the existing machine.



Vmax



Liposonix



Vaginal

Summary of Advantages

L Procedure Time

30 - 45 Minutes

♥ Full recovery
Immediate

+ Number of Treatments

1 every 12 - 18 months

✓ Quick To Work

Immediate

📄 Software Perloaded
All functional
software perloaded
⚠️

Risks & Complications

Possible mild redness,
swelling and tenderness
post treatment

✂️ Multifunctional

5 handles optional

📈 Upgrade
Insert new handles

i

Duration Of Results

12 - 18 months



Before and after



Before

After



Before

After



Before

After



Before

After



Before

After



Before

After

Specifications

| | |
|------------|------------------|
| Power | 1200 W |
| Frequency | 4 ~ 7MHz |
| Net Weight | 25kg |
| Dimension | 400*350*1060mm |
| Voltage AC | 100-240v 50/60Hz |

Contact Us

Distributor: Dermal Solutions Australia PL
Medical Equipment Supplier
PO Box 1488 Noosaville QLD Australia
Number:+61753242266
Website:www.dermalsolutions.com.au
E-mail:info@dermalsolutions.com.au

