MINISPIR LIGHT

Does the Minispir Light use the same software as all other MIR products?

No, the Minispir Light has it's own software Winspiro Light.

What are the main differences between the Minispir Light and the Minispir 2?

The Minispir light is a more cost effective version of the Minispir 2 and measures a smaller range of parameters for those who don't need the full feature parameters of the Minispir 2.

FAQ

What parameters does this device measure?

FVC, FEV1, FEV1%, FEV6, PEF, FEF25–75, FIVC, Lung Age, VC, IVC

Does the Minispir Light software feature pre and post test comparison?

Yes, the Winspiro Light software does have the ability for comparison against Pre and Post test.

Can I save the report into my patient files or notes?

Yes, you are able to export the file to a Word Document or PDF and save a copy of the test results to your PC files.

Can I claim the medicare rebate?

Use of this device may allow MBS item claim (#11505 & #11506) Contact Medicare for more information Medicare Rebate Number is 11505 & 11506 <u>View page here 11505</u> <u>View page here 11506</u>

🔊 winspiroLIGHT 3.1 - [Spirometry Standard Mode] - [10]						
File Configuration Privacy Tools ?						
<u>8</u>		A	<u>.</u>	Ŵ		0
Patient	FVC	VC	POST BD	SpO2	Print Export.dog Export.pdf	Help
MARTINI DAVID - Male - Age 28 - 180 cm - 76 kg - BMI 23,46 - Caucasian						

Can this device operate without the PC?

No the minispir device is PC operated only

Where is our device serial number?

It is on the label on the rear of the device. If the label is missing please use the Winspiro PRO software to obtain this information.





Ph: 1300 009 663 Fax: 1300 099 300 Email: info@zonemedical.com.au 7/22 Mavis Court, Yatala, Qld, 4207 www.zonemedical.com.au