


DentiumUSA
Developed by Clinicians for Clinicians

Dentium Osteotome Kit


Osteotomes compress the bone laterally or move it apically, providing a denser bone interface and/or an apical protective "cushion."



Osteotome | Final drill type | Scale 1 : 0.4 / mm

XOF 20 A1 XOF 20 B1		Type A 	Type B 	
XOF 34 A1 XOF 34 B1				
XOF 38 A1 XOF 38 B1				
XOF 43 A1 XOF 43 B1				
XOF 48 A1 XOF 48 B1				

Type	REF	ØA	ØB
XOFK Type A (Convex)	XOF 20 A1	Ø1.7	Ø2.8
	XOF 34 A1	Ø2.3	Ø2.8
	XOF 38 A1	Ø2.7	Ø3.2
	XOF 43 A1	Ø2.8	Ø3.8
	XOF 48 A1	Ø3.0	Ø4.3
XOFBK Type B (Concave)	XOF 20 B1	Ø1.7	Ø2.8
	XOF 34 B1	Ø2.3	Ø2.8
	XOF 38 B1	Ø2.7	Ø3.2
	XOF 43 B1	Ø2.8	Ø3.8
	XOF 48 B1	Ø3.0	Ø4.3



XOFK (Type A)

XOFBK (Type B)

DASK

Dentium Advanced Sinus Kit

DentiumUSA
Developed by Clinicians for Clinicians

DentiumUSA

6731 Katella Ave. Cypress, CA 90630 Tel. 1-877-304-6752 Fax. 1-714-226-0019

Website / E-mail

www.dentiumusa.com, cs@dentiumusa.com

PN2434-6