

Dropside Padding

Provide extra protection for recovering patients with a selection of padded bumper systems for trolley dropsides.

Key Benefits

- Stiff board backed bumper system
- Lift and roll feature allowing use without removal
- Quick install pocket system
- 25mm high density padding
- Welded cover materials for easy cleaning
- Sizes to fit any dropside style

Product Codes

- Padded Bumpers (pair); MPT-CA-PDB-P
- Padded Pockets (pair); MPT-CA-PDB-D

Why Buy

A recovering or restless patient can cause themselves serious harm on the metal risers and rails of standard dropsides. We have developed some ways for the carer to quickly and easily increase the level of protection and provide a soft, padded surface for patient's limbs. Our board system can be flipped to the outside of the stretcher allowing the dropsides to be raised and lowered without having to remove the bumpers from the trolley. Pockets can slide on over the dropside which are quick to use, however must be removed to lower the sides.



Padded Bumpers



Unique access



Padded Pockets



Removable design



Sizing options

v-20012023

"Confidence & Care, every day"

selectpatientcare.com.au

select