



Molini™ Mouth Opener

Great Help for Challenging Situations.

The Molini™ Mouth Opener is a great option for oral access when patients are uncooperative, bedridden, or have poor muscle control. The tapered end acts as a wedge to help get the mouth open when teeth are clenched, and the trim design allows greater visibility of the oral environment. The handle has a bulbous end with a secure grip surface to prevent slipping.

This latest mouth support in our line of products is made from a hard autoclavable plastic and has a tapered end for better insertion. It is ideal when the patient will not or cannot open the mouth. The Molini Mouth Opener was developed in conjunction with a technique where the clinician simultaneously holds the patient's head and the Molini tool, as shown in the photo, with the appropriate points of support intra and extra orally. However, it can also be used for cheek/tongue retraction and as a standard bite block.



Molini™ Mouth Opener	Code
Molini™ Mouth Opener	MMO-1



Creating dental products that help you provide safe, comfortable care.

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

www.mobilemedical.com.au

Open Wide® Reusable Mouth Prop

This Mouth Prop Gives You Control for Safe and Effective Care.

Do you use a mouth prop for patients who have difficulty keeping the mouth open reliably? Those same patients are the ones most likely to attempt to dislodge the prop. Our unique **Wraparound** mouth prop allows the clinician to maintain control with a handle that wraps to the outside of the cheek so it will not block oral access. The mouth prop is designed with a perfectly placed hole to keep a saliva ejector in position. The prop is made of plastic, over-molded with a crush-resistant latex-free* material, for oral access that is safe and comfortable.



Latex-Free*

* Not made with natural rubber latex

Tip: Both styles can be autoclaved repeatedly without getting "gummy".



Open Wide® Reusable Mouth Prop with Wraparound Handle



Open Wide® Reusable Mouth Prop with Tab Handle

Shown are sizes Large (green), Medium (blue) and Small (purple) with **Wraparound** handle.



Open Wide® Mouth Prop with Wraparound Handle

OWR-2S	2 Small Mouth Props/WRAPAROUND
OWR-2M	2 Medium Mouth Props/WRAPAROUND
OWR-2L	2 Large Mouth Props/WRAPAROUND
OWR-S/M	1 each Small & Medium Mouth Prop/WRAPAROUND
OWR-M/L	1 each Medium & Large Mouth Prop/WRAPAROUND
OWR/3 SZS	1 each Small, Medium & Large Mouth Prop/WRAPAROUND

The **Tab** mouth prop is streamlined for procedures that need greater clearance. It has a hole in the tab to allow floss to be threaded for quick retrieval of the mouth prop.

Shown are sizes Large (green), Medium (blue) and Small (purple) with **Tab** handle.



Open Wide® Mouth Prop with Tab Handle

OWT-2S	2 Small Mouth Props/TAB
OWT-2M	2 Medium Mouth Props/TAB
OWT-2L	2 Large Mouth Props/TAB
OWT-S/M	1 each Small & Medium Mouth Prop/TAB
OWT-M/L	1 each Medium & Large Mouth Prop/TAB
OWT/3 SZS	1 each Small, Medium & Large Mouth Prop/TAB

Creating dental products that help you provide safe, comfortable care.


www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au


Open Wide[®] Mouth Rest

A Reliably Open Mouth is Important for In-Office and At-Home Care.

The **Open Wide[®] Mouth Rest** is helpful for those who lack the cooperation or motor control to keep the mouth reliably open. It is made from special high-density foam with a rigid core to add stability. It is sturdy enough to resist biting pressure, yet is comfortable for the person receiving care. The mouth rest has a long handle for safe, easy control. It is considered a disposable item in the office but can also be cleaned and given to a trained caregiver for repeated use at home.



In-Office



At-Home

U.S. Patent No. 5,735,691

Open Wide[®] Mouth Rest

Description

OWD-50
OWD-100

50 pack (Regular)
100 pack (Regular)



Tip: The Thick Mouth Rest has 40% more foam to handle aggressive biting pressure.

Open Wide[®] Mouth Rest Thick

Description

OWDT-50
OWDT-100

50 pack (Thick)
100 pack (Thick)

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand any applicable regulations.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Surround[®] Toothbrush

This Unique Toothbrush Cleans All Surfaces of the Teeth at the Same Time.



Now with a beefier handle!



The Surround[®] Toothbrush is a great choice for patients with low tolerance for daily toothbrushing. It has three bristle pads to clean all surfaces of the teeth at the same time. Whether the user is brushing her own teeth or someone else is doing the brushing, the Surround Toothbrush can make the process easier because more tooth surface can be cleaned with each stroke.



A soft, rubber-like material covers the length of the toothbrush to help protect the oral environment if the person bites down unexpectedly, or falls while using the brush. The soft, end-rounded bristles and heart-shaped, compact head make brushing safe and comfortable.

** Not made with natural rubber latex*



A reusable stay-clean cover comes with every toothbrush!

Regular - Domestic Professional Pricing **Code**

12 Pack	STB-12
72 Pack - Professional Sale Only	STB-72
144 Pack - Professional Sale Only	STB-144

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Surround[®] Toddler Toothbrush

This Toothbrush Helps Teach Young Patients How to Brush.



The Surround[®] Toddler Toothbrush has a unique design that forms a “u” shape to hug the lingual, occlusal and facial surfaces, encouraging the child to move the toothbrush in contact with both teeth and gums. As the child learns, it is natural for the degree of cooperation to vary. With 3 bristle pads, more tooth surface will be cleaned within whatever amount of time the child tolerates toothbrushing.

New! The addition of a “smiley” face on the back of the toothbrush head is our latest development, and it is sure to please the under-6 crowd. The neck has also been made a little beefier, for little ones who have a tendency to bite down on the toothbrush.

** Not made with natural rubber latex*

Toddler - Domestic Professional Pricing

Code

12 Pack
72 Pack - Professional Sale Only
144 Pack - Professional Sale Only

STBC-12
STBC-72
STBC-144

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Stay N Place® Chair Cushion

This Dental Chair Cushion Provides Physical and Emotional Comfort.



Each chair cushion comes with Roll-N-Store Straps for easy carrying and storage.

When body shape and chair shape are not a good match, or when patients become fearful at the sight of a dental chair, the **Stay N Place® Chair Cushion** can help. It presents a welcoming image to help calm patient anxiety. As the patient settles in, the soft foam shifts inside the pillow, filling gaps in the chair and providing a secure, comfy feeling. The Stay N Place Chair Cushion curves in at the neck, allowing your patient's head to tilt back and the chest to be lifted, while leaving you room to maneuver.

A non-slip bottom helps the cushion stay where you place it on the dental chair. The chair cushion is covered in a soft, yet durable vinyl with a choice of dark teal or cream color.

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand any applicable regulations. Never leave patient unattended.

Stay N Place® Chair Cushion

Stay N Place® Chair Cushion

Dark Teal

SPCCdt

Cream

SPCCc

Stay N Place® Chair Cushion For Toddlers

The Infant Chair Cushion Makes it Easier to Accommodate Toddlers.

An over-sized dental chair can feel insecure for toddlers, and it can cause you to bend and twist in order to provide treatment. The **Stay N Place® Chair Cushion for Toddlers** creates a safe, cozy lift for the patient, presents a professional image to parents, and helps you maintain better posture.

The chair cushion is like a bean bag, with soft foam inside. Air escapes through side vents as the child is placed upon it. This settles the cushion into the curves of the dental chair while securely nestling the youngster. It has a non-slip bottom, and comes with a safety belt for added security.

Designed for dentistry, the Stay N Place Chair Cushion has a soft, durable vinyl cover. Clean with a disinfecting agent and wipe dry with a soft cloth.



Before: The dentist is bending and twisting to do her job.



After: The patient is at a more ergonomic height for dental treatment.

Stay N Place® Toddler Chair Cushion

Stay N Place® Chair Cushion

Dark Teal

SPCC INFdt

Cream

SPCC INFc

The Stay N Place Chair Cushion for Toddlers brings the patient to a more ergonomically correct height.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Stay N Place® Pillows

Support in the Right Places Facilitates Good Clinical Outcomes.

No chair can accommodate all patients equally well. If the patient falls short of the chair's headrest, he or she may not be properly supported. **Stay N Place® Pillows** provide the support you need to do your job effectively and the comfort the patient needs to sit without fidgeting.



Neck Pillow

The Neck Pillow supports the neck, helps keep the head tilted back, and can be used as a shoulder roll.

Pediatric - 3.5" H x 3.5" W x 10.5"L
Adult - 4.5" H x 4.5" W x 12.5"L



Knee Pillow

The Knee Pillow helps to center the patient and prevents sliding.

Pediatric - 4.5" H x 4.5" W x 12.5"L
Adult - 6" H x 6" W x 20"L



Plastic Sheaths

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au



Head Cradle

The sides of our Head Cradle push up into a cradling position when the patient puts pressure on the center section.

Pediatric - 3" H x 3.5" W x 9.5"L
Adult - 4" H x 4" W x 12"L

A. Stay N Place® Neck Pillow	Dark Teal	Cream
Neck Pillow - Pediatric	SPNP-Pdt	SPNP-PC
Neck Pillow - Adult	SPNP-Adt	SPNP-Ac
B. Stay N Place® Head Cradle	Dark Teal	Cream
Head Cradle - Pediatric	SPHC-Pdt	SPHC-Pc
Head Cradle - Adult	SPHC-Adt	SPHC-Ac
C. Stay N Place® Head Pillow	Dark Teal	Cream
Head Pillow	SPHdt	SPHc
D. Stay N Place® Conforming Pillow	Dark Teal	Cream
Conforming Pillow	SPCPdt	SPCPc
E. Stay N Place® Knee Pillow	Dark Teal	Cream
Knee Pillow - Pediatric	SPKP-Pdt	SPKP-Pc
Knee Pillow - Adult	SPKP-Adt	SPKP-Ac
F. Stay N Place® Plastic Sheaths	Code	
Sheaths for Pillows (100)	SH-100	
G. Stay N Place® Head & Neck Rest	Dark Teal	Cream
Stay N Place® Head & Neck Rest	SPHNRdt	SPHNRc

Stay N Place® Pillows

Specially Designed Pillows Support the “Chin Up” Position.

All our pillows have straps that attach behind the dental chair to keep the pillow in place, even when the patient sits forward. The vinyl pillow is soft to the touch, resists fluids and is easily cleaned using a surface disinfectant. Pillows are available in dark teal or cream.



The Head Pillow is designed for comfort and gives your patient a sense of security.

One size - 4" H x 11" W x 14" L



The Conforming Pillow takes the shape of the patient's head and neck and has a center cutaway to accommodate pulled-back hairstyles.

One size - 4" H x 7" W x 9" L


Stay N Place® Head & Neck Rest

New, Trimmer Backboard Supports the Patient's Head and Neck

Our Stay N Place® Head & Neck Rest allows the patient to be treated in his own chair - no need to transfer. New! The head area has been narrowed to improve oral access and permits the pillows to be better located to facilitate the “chin up” position.

This low-tech device is a real multi-tasker. It can be used for in-office as well as off-site care. As the patient settles in, the pillows can easily be moved to individualize support. Designed with dentistry in mind, the Head and Neck Rest is easily cleaned with the same protocol you use for dental chairs.

The Stay N Place Head and Neck Rest makes it possible to treat your patient without transferring.



The Head and Neck Rest attaches to the patient's chair.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Stay N Place® Booster Seat

The Booster Seat Saves Your Back and Provides Comfort for Young Patients.



Do you “scoot” patients to the top of the chair, only to have them slide back down when the chair is brought upright? That repeated maneuver is unsettling for the patient, causes back problems for you and lost time for the practice.



Removable lumbar support

The Stay N Place® Booster Seat positions your patient 5” higher in the chair, providing better access and increasing efficiency. This simple addition to the dental chair is more ergonomically correct for you and more comfortable for your patient. It now comes with a removable lumbar support pillow.

Important: Use appropriate monitoring, and never leave patient unattended.

The Stay N Place Booster Seat can be used with or without the lumbar support.

Stay N Place® Booster Seat	Dark Teal	Cream
Stay N Place® Booster Seat	SPBSdt	SPBSc

Rainbow® Toe Guard

This Cleverly Designed Pillow Keeps Patient from Interfering with Monitors.



Left: The patient can dislodge the monitoring device.

Right: The Rainbow Toe Guard is an effective yet comfortable way to solve the problem.



The Rainbow® Toe Guard prevents the patient from using one foot to kick a monitoring device off the other foot. The specially designed pillow is placed between the feet and secured with ankle straps. The belt is used around the pillow and stabilizing board or around the dental chair. This is a reliable way to keep the patient monitored and comfortable.

Rainbow® Toe Guard	Dark Teal	Cream
Rainbow® Toe Guard - One size	RTGdt	RTGc

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au


Stay N Place® Lap Board

Our Lap Board Provides Support and Cleanliness for Knee-to-Knee Exams.

For professionals who use the knee-to-knee technique for infant oral exams, our **Stay N Place® Lap Board** creates a more secure experience for the baby, and a more comfortable and hygienic experience for the adults. The Stay N Place Lap Board is made with a rigid platform, layered with foam and upholstered in soft, durable vinyl. Raised sides help keep baby safely on the Lap Board.



The Stay N Place Lap Board provides a safe, rigid surface for the child while allowing distance between the adults.



Stay N Place® Lap Board

Stay N Place® Lap Board
100 Protective Sheaths

24"L x 16"W
12" x 8" x 30" gusseted

Dark Teal


SP-LAP BDdt
12" x 8" x 30" gusseted

Cream

SP-LAP BDC
SH-100/LAP

Stay N Place® Lap Cushion


The Lap Cushion Provides Another Option for Infant Oral Exams.



The **Stay N Place® Lap Cushion** is a molded foam product that provides support during infant oral exams. Like the Lap Board, it has side rails for safety but also has the "cushiness" that parents appreciate.

Stay N Place® Lap Cushion

Stay N Place® Lap Cushion SP-FLEX 24"L x 17"W
100 Protective Sheaths SH-100/LAP 12" x 8" x 30" gusseted



The Stay N Place Lap Cushion can be held behind the child as he is tilted back. Mom is in sight, which is reassuring to the baby.

Important: Use appropriate monitoring, and never leave patient unattended.

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com


International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Stay N Place® Airway Positioner - Foam

Our Airway Positioner Helps When “Chin Up” Isn’t Happening.



The Stay N Place® Airway Positioner maintains proper “chin up” posture by supporting the patient’s upper torso. Customers use it for sedation and for routine dentistry. With comfortable yet sturdy construction, you can rely upon consistent torso lift throughout the procedure.



The SP-AP FM is 2.5” high at the chest area, 18” long and 10” wide. Our model is 9 years old.

Important: Use appropriate monitoring, and never leave patient unattended.

Stay N Place® Airway Positioner - Foam

Stay N Place® Airway Positioner - Foam	SP-AP FM
100 Protective Sheaths	SH-100/AP

Stay N Place® Airway Positioner - Vinyl

Our Other “Chin Up” Option Has Three Adjustable Pillows.

The Stay N Place® Airway Positioner in Vinyl has two head and neck cradles and a wedge-shaped pillow to lift the chest area. These vinyl-covered pillows can be repositioned to accommodate a range of patients. It can be used with or without the head strap. Available in dark teal or cream vinyl.



The SP-AP VIN is 3” high at the chest, 19” long and 10.5” wide.

Stay N Place® Airway Positioner - Vinyl

Stay N Place® Airway Positioner - Vinyl
100 Protective Sheaths

Dark Teal

SP-AP VINdt
SH-100/AP

Cream

SP-AP VINc
SH-100/AP

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Rainbow[®] Stabilizing System

Protective Stabilization Helps Keep Everyone Safer

You can't always predict which patients will become suddenly agitated when you or a staff member is using a sharp instrument. If there is a risk of sudden, unsafe movement, you owe it to your patients and to your practice to use safety equipment specifically designed for dentistry.

The **Rainbow[®] Stabilizing System** consists of a colorful, breathable mesh wrap that prevents patient overheating, and a padded vinyl-covered board for stability. Together, they help prevent kicking, grabbing and thrashing that can pose a danger for all involved.

See pages 14-17 for more details on our Rainbow[®] Stabilizing System and Accessories.



Remember to use appropriate informed consent protocols specifically for stabilization (as distinct from the dental procedure itself.)

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Rainbow® Board

The Rainbow® Board Helps You Put Safety First During Dental Procedures.

The **Rainbow® Board** provides the foundation for the stabilizing system. It is available as a **Hinged Board** or a **Straight Board**. The Hinged Board bends with the dental chair and allows the patient to be moved to an upright position quickly, which is helpful should the patient have a medical emergency during a procedure. We offer the Straight Board for practitioners whose equipment is better suited to a rigid board, or for those who simply prefer it. Rainbow Boards are used in conjunction with Rainbow Wraps (next page).

The helpful **Rainbow Head Tilt** keeps the upper torso lifted and the head in the “chin up” position.

The **Rainbow Knee Support** relieves pressure on the lower back. It can be used with or without the stabilizing system.



Rainbow® Head Tilt & Knee Support	Length/Width/Height	Dark Teal	Cream
Rainbow® Head Tilt - Small/Medium	9.25" x 9.5" x 2.5"	RHT-S/Mdt	RHT-S/Mc
Rainbow® Head Tilt - Large	12.0" x 11.0" x 3.0"	RHT-Ldt	RHT-Lc
Rainbow® Knee Support - One Size	13.0" x 6.5" x 3.0"	RKSdt	RKSc

The vinyl covered **Rainbow Board** is padded for patient comfort.



The Large Straight Rainbow Board articulates slightly at the bottom to conform to the contours of the chair.



The Hinged Boards fold for easy storage.



The Large Hinged Board articulates in two places to conform to the curves of the chair.

Rainbow® Boards

Height	Width	Rainbow® Board	Dark Teal	Cream
2 Feet	10"	Infant Board	RBldt	RBlc
3 Feet	10"	Small Board	RBSdt	RBSc
4 Feet	11"	Medium Hinged Board	RBHMdt	RBHMc
4 Feet	11"	Medium Straight Board	RBMdt	RBMc
5 Feet	15"	Large Hinged Board	RBHLdt	RBHLc
5 Feet	15"	Large Straight Board	RBLdt	RBLc

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Rainbow® Wrap

Our Unique Wrap is Designed Specifically for Dentistry.

The **Rainbow® Wrap** keeps arms and legs out of harm's way while sharp instruments are in use. Everything about it has been designed for the needs of a dental office. The breathable mesh avoids overheating that often accompanies patient agitation. An opening in the wrap allows you to conveniently use monitoring devices. The wrap has a friendly look, so it is better accepted by patients and parents. But it doesn't stop there. The wrap is made of a moisture-resistant, washable fabric - so even laundering is made simple.



Rainbow® Wristlets	Code	Length/Width
Wristlets Pair/Infant	RWC Infant/PR	(7"L x 1.5" W)
Wristlets Pair/Small-Med	RWC SM-Med/PR	(9"L x 1.75" W)
Wristlets Pair/Large	RWC Lg/PR	(11"L x 2" W)

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand all applicable regulations regarding the use of protective stabilization. Use appropriate monitoring, and never leave patient unattended.



The wristlets are cushioned for comfort, and come in 3 sizes.



Our new Knee Belt has a layer of comfy foam on the inside, and an easily cleaned outer layer.



The RWSB is made of a hook and loop, and is long enough to go around the dental chair.

Rainbow® Knee Belt	Code
Medium Knee Belt	RKB-Mp (6"W by 40"L)
Large Knee Belt	RKB-Lp (9"W by 43"L)
X-Large Knee Belt	RKB-XLp (9"W by 58"L)

Rainbow® Wide Safety Belt	Code
Rainbow® Wide Safety Belt	RWSB (4"W by 60"L)

Rainbow® Wraps with easy-clean back panel

Code	Rainbow® Wrap	Length/Width	Wrap Comes With.....
RWlp	Infant Wrap	20"/43"	Wristlets, Laundry Panel, Storage Bag*
RWSp	Small Wrap	25"/35"	Wristlets, Laundry Panel, Storage Bag*
RWmp	Medium Wrap	32"/44"	Knee Belt, Wristlets, Laundry Panel, Storage Bag
RWLp	Large Wrap	34"/49"	Knee Belt, Wristlets, Laundry Panel, Storage Bag
RWXLp	X-Large Wrap	44"/55"	Knee Belt, Wristlets, Laundry Panel, Storage Bag
RWXXLp	XX-Large Wrap	44"/70"	Knee Belt, Wristlets, Laundry Panel, Storage Bag

* Knee Belt Not Required for Small or Infant Wraps

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Accessories for Your Rainbow® System

Accessories Help Customize the Rainbow® System for Your Practice.

The Head Stabilizer is attached with screws for reliable stability.



The Rainbow® Head Stabilizer is a comprehensive head control device that minimizes both up-and-down and side-to-side movement. It attaches to the Rainbow Board with screws for reliable stability. Adjustable side panels tilt in at a slight angle to improve effectiveness. The thick cushions on the side panels allow for a snug fit while maintaining patient comfort. *Note: The Head Stabilizer must be used in conjunction with the Rainbow Board.

Rainbow® Head Stabilizer	Dark Teal	Cream
Rainbow® Head Stabilizer for Infant Boards	RHlIdt	RHlIc
Rainbow® Head Stabilizer for S/M Boards	RHlCdt	RHlCc
Rainbow® Head Stabilizer for Large Boards	RHlAdt	RHlAc
Replacement Cushions - Infant	RHSC - INFdt	RHSC - INFc
Replacement Cushions - Sm/Med	RHSC - S/Mdt	RHSC - S/Mc
Replacement Cushions - Large	RHSC - Ldt	RHSC - Lc
Replacement Head Strap - RHS - 2		

Rainbow® System Storage Bags

Rainbow® Storage Bags are designed specifically for the Rainbow Stabilizing System. They are made of a tough canvas-like fabric to extend the life of your investment. They will accommodate your board, wrap, and belts and even a few extra items.

Rainbow® Storage Bags - Canvas-like Fabric


- RST-S Rainbow® Storage Bag for Small or Infant System
- RST-M Rainbow® Storage Bag for Medium System
- RST-L Rainbow® Storage Bag for Large System
- RST-H Rainbow® Storage Bag for Hinged System



Shown is the Storage Bag for the Medium System (RST-M).

Rainbow® Triangle Pillows

The Rainbow® Triangle Pillow fills the gap in the dental chair where the seat meets the back. It is covered with a non-slip material, so it stays where you put it. It is particularly helpful when using a straight board. Our customers report it also makes a handy knee pillow, or comfort pillow to fill the gap in the chair - even when a board is not in use.



Rainbow® Triangle Pillow	Dark Teal	Cream
Rainbow® Triangle Pillow	RTPdt	RTPc

Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Robot Arm™ Stabilizers

Help Teach Young Patients How To Sit Safely for Dentistry.



The **Robot Arm™ Stabilizers** are a perfect tool to teach young patients how to keep hands out of harm's way. They start out as flat sheets and are formed into cylinders around arms or legs to prevent grabbing, kicking and bending movements. The Velcro® material makes the arm stabilizers easy to apply and remove during the appointment, as a patient learns. They are colorful and fun, appealing to parents and patients alike.

The Robot Arm™ Stabilizers are easy to clean. See the back cover for a QR code that takes you to a page on our website with **cleaning instructions** for all of our products.

Important: It is the clinician's responsibility to assure the appropriateness of any medical device for a given patient. Be sure to understand all applicable regulations regarding the use of protective stabilization. Never leave patient unattended.

Velcro® is a registered trademark of Velcro BVBA.

Code	Size (Circumference x Height)	Usage Guide
RAS-SM - Pair	10" x 8" (up to 9" arms/legs)	Ages 1-4
RAS-MED - Pair	12" x 10" (up to 11" arms/legs)	Ages 5-8
RAS-LG - Pair	14" x 12" (up to 13" arms/legs)	Ages 9-12
RAS-XL - Pair	16" x 18" (up to 15" arms/legs)	Teens
RAS-XXL - Pair	20" x 18" (up to 19" arms/legs)	Adults

Patient Education on Stabilization

These Educational Tools Help You Explain Stabilization.



The **Protective Stabilization Model** is a 11" stabilizing system complete with "patient". It helps you demonstrate to caregivers how a patient is safely swaddled during treatment. It will become a key part of your informed consent procedure for parents, especially when there is a language barrier. The model can also be used to help patients know what to expect.

Riley and the Hugging Blanket is a beautifully illustrated book about a patient who gets upset at the dental office, despite trying to be brave. The "Hugging Blanket" is explained to Riley, who accepts its use. Riley learns to relax, and soon the Blanket is no longer needed. Riley is very proud! The teaching book is available in English and in Spanish.

Educational Tools - Please specify English or Spanish for the Books

- PSM Protective Stabilization Model
- HBK-1 Riley and the Hugging Blanket - Single Copy
- HBK-3 Riley and the Hugging Blanket - 3 Pack
- HB-PKG Educational Package (1 PSM and 3 Books)



Creating dental products that help you provide safe, comfortable care.

www.specializedcare.com

International Distributor: Mobile Medical Systems International /02 9252 4410/info@mobilemedical.com.au

Size Guidance



What Size To Buy?

Up to this approx. weight/age:	Recommendation
26 pounds(Infants).....	Infant Board - Infant Wrap
35 pounds(1-4 years old).....	Small Board - Small Wrap
60 pounds(5-8 years old).....	Medium Board - Medium Wrap
145 pounds(9-12 years old).....	Large Board - Large Wrap
185 pounds(Teen/Adults).....	Large Board - X-Large Wrap
250 pounds(Teen/Adults).....	Large Board - XX-Large Wrap

Solutions for Your Everyday Challenges

Infant patient requires knee-to-knee exam	Lap Board.....Page 11 Lap Cushion.....Page 11
Patient lacks cooperation skills	Stabilizing System..... Pages 13-17 Protective Stabilization Model.....Page 18 Robot Arm™ Stabilizers.....Page 18
Patient does not fit well in the dental chair	Chair Cushions.....Page 7 Dental Chair Pillows.....Pages 8 & 9 Head & Neck Rest.....Page 9 Booster Seat.....Page 10
Patient won't keep the mouth open	Mouth Opener..... Cover Mouth Props.....Page 2 Mouth Rests.....Page 3
Staff Training for patients with special needs	Protective Stabilization DVD.....Page 19 D-Termined Program (Patients with Autism).....Page 19
Patient does not cooperate for daily oral hygiene	Mouth Rests.....Page 3 Toothbrushes.....Pages 4 & 5
Patient tries to kick pulse oximeter off toe	Toe Guard.....Page 10
Caregiver and/or patient requires education	Oral Care Kits DVD, Mouth Rests, Toothbrushes.....Page 6 DVD: How to Brush the Teeth of Another Person.....Page 6 Oral Care Programs for Group Living Situations.....Page 6 Protective Stabilization Model.....Page 18
Patient does not maintain "chin up" position	Dental Chair Pillows.....Pages 8 & 9 Airway Positioners.....Page 12 Head Tilt and Knee Support.....Page 14