

# Stretcher Dropside Options



## Why Buy?

There are a range of dropside options available to customise any stretcher to suit your special requirements. We understand that there are different ways that patients must be transported and procedures performed. Safety comes first, which is why we've developed a selection of side rail options, to ensure you can have the right outcome for your department's requirements.

## Key Benefits

- Multi-riser dropside options for closer spacings
- Split fold option for safer patient access
- Horizontal rail designs for longitudinal support
- Flush folding option for mattress transfer or x-ray
- Mid position lock option for easier patient access

v-04122025

# MODSEL

[modsel.com](https://modsel.com)

# Stretcher Dropside Options



**12 riser support dropside option**  
Can be installed folding towards the head or foot end  
Fits all stretcher models except the eye stretcher VSM-10V



**Split dropside option**  
Each section folds down outwards towards the head and foot ends  
Fits all stretcher models except the eye stretcher VSM-10V



**Mid-position catch option**  
Stops the dropside halfway at 45 degrees for partial support  
2 presses to lower and then ratchet lifts back to highest/fully upright position



**VSM-10V split dropside option**  
Improved access without compromising safety.  
Split-section dropsides provide greater flexibility when positioning patients for other specialist procedures such as Endoscopy. Each section can be lowered independently for close access to the head and torso, while maintaining protection where needed.

## Product Codes

- 12 riser support dropsides - VSM-107
- Split dropsides (both sides) - VSM-106
- Split dropsides (Patients Left Hand Side) - VSM-106L
- Split dropsides (Patients Right Hand Side) - VSM-106R
- Mid-Position Dropside Catch (Pair) - VSM-111

v-04122025