



Demonstration Detail Agreement Form: Autochair Hoist Range

Installation item (mm) measurements when stripped of removable attachments

Folded  
Height \_\_\_\_\_ Width \_\_\_\_\_  
Length \_\_\_\_\_ Weight \_\_\_\_\_ Kg


Car Boot and boom Measurements (mm) -

Height \_\_\_\_\_  
Width \_\_\_\_\_ Depth \_\_\_\_\_  
Radius from post base to mid scooter \_\_\_\_\_  
See more detailed measurements below

DATE OF DEMO

<b>Customer</b>	<b>Address</b>	<b>Phone</b>	<b>Email</b>
<b>Therapist</b>			
<b>Autochair Hoist Chosen</b>			
<b>Item to be hoisted:</b>			<b>Notes</b>
<b>Vehicle Make:</b>	<b>VIN:</b> _____	<b>Rego:</b>	<b>Year:</b>
	<b>Model:</b>		<b>Kms:</b>
<b>Any special requirements</b>			

A                      B                      C                      D                      E                      G                      F



Dealers, Get these measurements from <http://autochaironline.co.uk/> Or measure if needed

**Note for Dealer** A-B = 120mm approx, so please allow this clearance when calculating fitting. Deduct 100mm from the total width of the rear opening to calculate G. Please visit [autochaironline.co.uk](http://autochaironline.co.uk) and sign in with your username and password to see the rear measurements for 300 vehicles and the colour picture fitting instructions. Lost your username and password – just call Mobilitycare on 03 95688383 or email [info@mobilitycare.net.au](mailto:info@mobilitycare.net.au)

**DECLARATION:** I understand that some trim and carpet may either need to be removed, cut or drilled to facilitate the installation of the Autochair-Hoist. Some rear seats may not be useable when hoist is operated I agree that the Autochair Hoist has been demonstrated to me and that I am happy that I can satisfactorily operate it in accordance with the above conditions:

Signed by Client \_\_\_\_\_ Date \_\_\_\_\_

Signed by Therapist \_\_\_\_\_ Date \_\_\_\_\_

For Dealer \_\_\_\_\_ Date \_\_\_\_\_