

# SV-A Series



*2ml Sample Measurement  
Now Possible!*

- **2ml Sample Measurement (SV-1A)**
- **X-Y-Z Stage for Fine Position Adjustment**
- **Anti-corrosive, Titanium Sensor Plates**
- **Standard WinCT-Viscosity, Serial-USB Converter**



- *The sensor unit can be detached for field research purposes.*
- *A portable carrying case is provided as standard.*



**AND** ...Clearly a Better Value  
A&D Company, Limited  
<http://www.aandd.jp>

# SV-A Series



## Specifications

	SV-1A	SV-10A
Measurement Method	Sine-wave Vibro Viscometer using Tuning Fork Vibration Method (30Hz)	
Viscosity Measurement Range	0.3~1000mPa.s	0.3~10000mPa.s
Accuracy	1% Repeatability (S.D., 20~30 °C, No condensation)	
Minimum Sample Amount	2ml~	10ml~
Temperature Measurement	0~160 °C / 0.1 °C (32~320 °F / 0.1 °F)	
Power Supply	AC Adaptor (Approx. 14VA)	

## Options

- AX-SV-51 Stand Set with X-Y-Z stage
- AX-SV-52 X-Y-Z stage
- AX-SV-53 Software Package with serial USB converter
- AX-SV-54 Sample Container Set  
45ml X 5pcs., 10ml with cover X 5 pcs., Polycarbonate  
13ml X 2 pcs., Glass  
Glass Container Holder X 1 pc.  
Circulation Water Jacket X 1 pc.
- AX-SV-55\* 2ml with lid X 10 pcs., Polycarbonate  
2ml X 10 pcs., Glass  
Sample Container 45ml X 5 pcs  
Container Holder for 2ml X 5 pcs., Polycarbonate  
Container Table for 2ml container X 1 pc.  
Circulation Water Jacket X 1 pc.
- AX-SV-56\* Container Holder for 2ml X 5 pcs., Transp./Black
- AX-SV-57\* Container Table for 2ml container X 2 pcs.
- AX-SV-58\* 2ml with lid X 100 pcs., Polycarbonate
- AX-SV-59\* 2ml with lid X 5 pcs., Glass  
Container Table for 2ml container X 1 pc.

\* for the SV-1A only



AX-SV-52



AX-SV-54



AX-SV-55

# A&D

...Clearly a Better Value

**A&D Company, Limited**  
3-23-14 Higashi-Ikebukuro, Toshima-ku, Tokyo 170-0013 JAPAN  
Telephone:[81](3) 5391-6132 Fax:[81](3) 5391-6148  
<http://www.aandd.jp>

**A&D ENGINEERING, INC.**  
1756 Automation Parkway, San Jose, CA 95131 U.S.A.  
Telephone:[1](408) 263-5333 Fax:[1](408) 263-0119

**A&D MERCURY PTY. LTD.**  
32 Dew Street, Thebarton, South Australia 5031 AUSTRALIA  
Telephone:[61](8) 8301-8100 Fax:[61](8) 8352-7409

**A&D INSTRUMENTS LTD.**  
Unit 24/26 Blacklands Way Abingdon Business Park,  
Abingdon, Oxon OX14 1DY UNITED KINGDOM  
Telephone:[44](1235) 550420 Fax:[44](1235) 550485

**<German Sales Office>**  
Große Straße 13 b 22926 Ahrensburg GERMANY  
Telephone:[49](0) 4102 459230 Fax:[49](0) 4102 459231

**A&D KOREA Limited**  
Manhattan Bldg. 8F, 36-2 Yoido-dong, Youngdeungpo-gu, Seoul, KOREA  
Telephone:[82](2) 780-4101 Fax:[82](2) 782-4280

**A&D RUS CO., LTD.**  
Vereyskaya str.112 Kuntsevo Block, 121357,  
Moscow, Russia  
Telephone: [7] (495) 937-33-44 Fax: [7] (495) 937-55-66