

# CentroGlide VL



## Do more in this maneuverable wheelchair.

With a vertical lift system to help users reach hard to get places, the CentroGlide VL has all the industry-leading features of the CentroGlide, and more.

### Communicate on your level.

Powered seating on the CentroGlide VL rises to an additional 30cm. This allows one to be at eye level while communicating in social settings, giving self-confidence and added visibility.

### Easy transfer.

Being able to raise or lower the seat helps transferring from the wheelchair to another surface like a bed, or sit to stand position, much easier. Adjusting the height makes the transfer simple, quick and smooth.

### Live life on your own terms.

The CentroGlide VL's powered lift system gives you added independence by reducing the requirements for home modifications. Tasks like cooking at the stove or reaching for that microwave is possible. You now have the ability to access those hard-to-reach high shelves without straining or needing to ask for assistance.

## Designed for comfort and convenience.

This industry-leading power wheelchair has the additional convenience of powered seat raising making it even more comfortable and easy to use.



### Customised cushioning.

For ultimate comfort, custom cushioning can be fitted by Glide or by the distributor.



### Active suspension.

Our patented 6-wheel APT suspension uses power traction to smoothly traverse kerbs and prevent nosedive when descending.



### Wheels.

The CentroGlide VL has 14" drive wheels and 7" castor wheels as standard. There is a range of options such as wider off-road tyres and larger or smaller castors.



### Colour choices.

Choose from an extensive colour palette for the ultimate customisation.



### Vertical Lift.

Replacing the standard seat tilt mechanism, users of the Centro VL can raise their seat by 30cm.



### Posterior Tilt

For comfort and convenience, the seat has a 50° posterior tilt.



### Battery life.

The powerful battery will give you the freedom of all-day operation with a range of up to 40km on a single charge.



# CentroGlide VL sets the benchmark.

## Features.

- Power base with removable seating module
- 50-degree seat tilt standard
- Available in Adult and Paediatric versions
- High quality German 2-4 pole high torque motors
- Vehicle tie down points
- Low maintenance with accessible electronics
- Turning circle between walls 1145mm
- Active Pozi-Track (APT) kerb climbing ability
- 80mm ground clearance
- Max ramp incline/decline 12° (as per AS3695 testing)
- Optional anterior tilt
- 30cm vertical lift seat with 46cm seat height
- Electronic brakes with manual park
- R-net LCD Std controls (dynamic optional)

## Safety.

- Full compliance AS 3695 175kg
- Crash tested to AS 3696:19 to max user weight 159kg

## Arrange a test drive.

Chat to a Glide Wheelchair specialist and experience the extraordinary features of the CentroGlide VL in person. Find your nearest distributor at [glide.com.au/nearest](http://glide.com.au/nearest)

## Specifications.



**Recommended User Weight**  
Up to 175kg AS 3695 full compliance



**Overall Length**  
900mm (foot plate removed)



**Seat Widths**  
34 to 67cm (custom width optional)



**Overall Width**  
620mm (wheel base)



**Seat Depths**  
31 to 64cm



**Seat Height**  
46cm (top of rail)



**Seat Raise Height**  
300mm



**Arm Rests**  
Flip up (adjustable)



**Overall Height**  
1030mm Top of push handle (custom height available)



**Leg Rests**  
Swing away detachable



**Brakes**  
Electronic and manual park



**Seat Tilt**  
0-50 degrees std



**Frame Colour**  
Wide choice



**Battery**  
70-80 amp/hr (260Lx170Wx220H) Sealed



[www.glide.com.au](http://www.glide.com.au)

(08) 9345 3400 | [sales2@glide.com.au](mailto:sales2@glide.com.au)

25 Ledger Road, Balcatta, Western Australia 6021

@GlideProductsAustralia

Glide Products