# Transfer Chair | VSM-20



# Why Buy?

## Save space, gain flexibility

Make better use of your valuable space with Modsel's compact VSM-20 stretcher chairs. Every aspect of design has been tested and refined to deliver ultimate comfort and efficiency. With an impressive 400kg/880lb lifting capacity and superb manoeuvrability, Modsel stretcher chairs optimize operator and patient positioning for a wide range of procedures

- Dual column design for exceptional strength and stability
- Pop up push handles can be used with patient sitting or lying flat
- Multi-layer mattress provides contoured comfort without adding bulk

# **Application**

Ideal for Day Surgery, Recovery, Transit Lounge, Fast Track ED & Cardiac

# **Key Benefits**

- Greater throughput: Patients can remain in the same chair for admission, procedure and recovery, reducing risk and enhancing efficiency.
- Leading weight capacity: Modsel stretcher chairs are built with a super strong subframe that can hold and lift a class-leading 400kg/880lb.
- Supercharged: The long-life, cordless, fastcharge battery system keeps things moving all day (and night).
- Clean and hygienic: The smooth, gap-free platform makes cleaning easy and prevents fluid ingress to meet the highest infection control standards.

v-05122025



# Transfer Chair | VSM-20



### Low and safe

Excellent 480mm/18.9" low height with fold-under dropsides.

The ultra-low height is combined with dropsides that glide down vertically under the chair top to provide safe access with no protrusions to catch on the patients legs or clothes.



#### Lies fla

The backrest and leg rest move independantly or together from chair to bed.

You now have the flexibility of sitting a patient up to almost 90 degrees or lying down completely flat.



#### Head and foot tilt

A generous 15 degrees top tilt both directions.

Position a patient how you need them or quickly lower their head down (and feet up) in an emergency.



#### Simple transfer

Conveniently stored handles unfold for staff friendly transportation.

Unfold the push handles and move forward effortlessly on the large castors. Our unique patented handle design correctly positions the carer while preventing unwanted forces being applied directly to the patient.



#### Forward tilt for safer transitions

Positioning support when entering or exiting the chair.

The chair's forward tilt angle helps patients move into a seated or standing position with greater stability. Ideal for recovery and discharge phases, this function supports safe, staff-assisted transfers while reducing strain on both patient and clinician.



## Effortless cornering

Perfectly centered steer wheel for ease of movement.

Lock in the steer system to ensure the chair tracks straight along hallways and then navigates the tightest turning circles possible.



### Multi-layer comfort mattress

Engineered to deliver support without unnecessary bulk.

The multi-layer design combines soft memory foam, contouring immersion layers and a protective base, topped with a 4-way stretch cover. Together, these materials reduce pressure points and enhance patient comfort during longer procedures or recovery periods.





## **Specifications**

- Overall Size: 780mm/30.7" wide x 1130mm/44.5" long (1800mm/70.9" lying flat). 100mm/4" wider option available
- Mattress Construction: 100mm/4? thick, premium 4-layer foam, with welded pressure relieving cover
- Mattress size: 1800mm/70.9" long x 650mm/25.6" wide with 500mm/19.7" deep seat. 100mm/4" wider option available
- Height Range: 480mm/18.9" 860mm/33.9"
  with lie-flat, forward and reverse tilt
- **SWL:** 400kg/880lbs
- Backrest: 0 to 85° with CPR emergency quick release
- Castors: Precision 150mm/6" twin wheel with 4wheel braking and center steering
- Trendelenburg: 15° both head and foot down
- Wheelbase: 590mm/23.2" wide overall x 750mm/29.5" long
- Undercarriage clearance: 140mm/5.5"

## **Product Codes**

- Standard Chair: VSM-20
- 100mm Wide Chair Option: VSM-201

v-05122025

