



HP-30^{with} PCA

Patient-Controlled Analgesia

- Drive lock, prevent unauthorized syringe replacement.
- Screen lock, prevent unauthorized change in parameter.
- More than 2000 historical events will be recorded for review.
- Clear visibility of important information with medical grade large touchscreen.

As a member of HP family , all advantages of HP technology are inherited.



PCA Controller




Syringe Anti Removal Cap



PCA Key

Specifications

Dimensions	258(W) × 75(H) × 152(D)mm
Weight	1.70kg (including the battery)
Accuracy	≤ 2%
Battery Operating	Time Min. 10h at 5ml/h
Initial Bolus	0 or 0.1-9999.99ml
Basal Rate	0 or 0.1-2000ml/h
Lockout Time	00:00:01-99:59:59 (hh:mm:ss)
PCA Bolus Amount	0.1-9999.99ml
Bolus Rate	0-2000ml/h
PCA Modes	P mode; CP mode; LCP mode



MEDCAPTAIN MEDICAL TECHNOLOGY CO., LTD.



Address: 12th Floor, Baiwang Research Building, No.5158 Shahe West Road, Xili, Nanshan, 518055 Shenzhen, Guangdong, PEOPLE'S REPUBLIC OF CHINA

Tel: +86 755 26953369

E-mail: info@medcaptain.com

Fax: +86 755 26001651

Website: www.medcaptain.com

MEDCAPTAIN is a trade mark of Medcaptain Medical Technology Co., Ltd.

Specifications are subject to change without prior notice.

©2019 Medcaptain Medical Technology Co., Ltd. All rights reserved.

EN-HP-30 with PCA-2P-Version 2.2

DISTRIBUTOR:

