

Braun ThermoScan[®] PRO 6000

Attribute	Description
ExacTemp™ & PerfectTemp™ Technology	Insures accuracy, even when placement of the thermometer in the ear canal is not ideal. ExacTemp™ helps support reliability and accuracy by detecting the stability of probe placement during measurement. PerfectTemp™ supports accuracy by detecting direction and depth of probe placement in ear and adjusts temperature calculation.
Pre-Heated Tip and Sensor	Activated when thermometer is removed from cradle allowing the 3 – 4 second warm up period to occur while loading the probe cover and inserting into the ear so thermometer is ready when user is ready to push measurement button.
Recharging Option	Charging Station provides 700 measurements on a fully charged Rechargeable Battery; 350,000 measurements over the life of Rechargeable Battery Pack.
Probe Cover Storage	20 covers (small cradle), 40 covers (large cradle), or 200 covers (charging base) opt ons.
Security Options	6ft and 9ft tether attaches to cradle with hex screw to secure attachment. Large cradle option can be mounted to wall. Electronic count down time security with option of 1 – 12 hours of operation until unit locks.
Back lit Display	Improves ability to read display in low room light conditions such as the PACU or night rounds.
Temperature Mode	The Braun PRO 6000 provides an oral equivalent reading to allow facilities to use a variety of Welch Allyn devices on patients (or other oral thermometers in conjunction with the PRO 6000) without having to convert between an oral and ear reading.
Unadjusted Ear Mode	User enabled. Easily identified with a “U” in the bottom right corner of the LCD display
Product Availability	The PRO 6000 is a professional thermometry device, sold through authorized Welch Allyn medical distributors. It is labeled “This thermometer is for professional use only. ”, and is sold with dedicated probe covers for the device which are marked as “For use with the PRO 6000 series”.

