

# MEDI GUARD

MEDICAL REFRIGERATORS AND FREEZERS

## Medi Guard Chest Freezer MG-86W50

# MEDI GUARD



### Temperature Control.

Ensure products are kept at the best temperature possible with a high precision microcomputer temperature control system.



### Quiet Operation.

Operates at 60dB(A) when compressor system is fully running. Quiet operation means less disturbance in the lab where the freezer is running.



### Security System.

Excellent alarm system equipped with a number of alarm functions - see specs for full list.

## MEDI GUARD

Leading specialists in vaccine fridges and freezers.

[eurochill.com.au](http://eurochill.com.au)

 1/13-15 Kevlar Cl, Braeside, VIC 3195

 (03) 8372 0435

 [sales@eurochill.com.au](mailto:sales@eurochill.com.au)

## Product structure data

<b>Temperature Range</b>	-40°C to -86°C
<b>Fridge Storage Volume(L)</b>	50L
<b>Net/Gross Weight (kgs)</b>	83/125
<b>Product Dimensions (mm) (wxdxh)</b>	695x740x985
<b>Interior Dimensions (mm) (wxdxh)</b>	430x400x350
<b>Packing Dimensions (mm) (wxdxh)</b>	810x830x1180

## Technical data

<b>Cooling Type</b>	Direct
<b>Defrost type</b>	Manual
<b>Climate Class</b>	N
<b>Rated Voltage/Frequency (V/Hz)</b>	220V/50Hz
<b>Rated Power (W)</b>	540
<b>Rated Current (A)</b>	2.5
<b>Decibel Rating (dBA)</b>	60
<b>Refrigerant/Amount (g)</b>	Mixed CFC free
<b>Sensor</b>	PT1000

## Control

<b>Controller</b>	Microprocessor
<b>Display</b>	LED display

## Alarms

<b>Yes</b>	Power failure, sensor failure, low battery, door ajar, high temperature, low temperature.
------------	---

## Accessories

<b>Castors/type</b>	4 castors
---------------------	-----------

**MEDI GUARD**

Leading specialists in vaccine fridges and freezers.

[eurochill.com.au](http://eurochill.com.au)

 1/13-15 Kevlar Cl, Braeside, VIC 3195

 (03) 8372 0435

 [sales@eurochill.com.au](mailto:sales@eurochill.com.au)