

Dental Materials · Laboratory Materials · Equipment



# **Excellent Hydrophilic Impression Material**

[Type 3]

Vonflex S Light XLV

Vonflex S Light XV

Vonflex Light XV

Vonflex S Light XV

Vonflex S Light XV

Vonflex S Light XV

Vonflex S Light XV

Vonflex Light XV

# **FEATURES AND BENEFITS**

- Vonflex S Light is an addition polymerization silicone composed of hydrophilic vinyl polysiloxane which helps to create precise shape of oral tissue. which meets the requirements of ISO 4823.
- Vonflex S Light is convenient to use and has an outstanding stability, helps to make precise impression taking.

## **INDICATIONS**

- Taking impression for crown and bridge
- Taking impression for inlay and onlay

# TECHNICAL DATA

- Excellent detail reproduction by optimal flowability
- Procedure of underwater test.
   (Place a coin under the water and apply the Vonflex S Light to the coin.)
- Vonflex S Light underwater test result





## **PACKAGE**

ALL PRODUCTS			PRODUCT CODE	KFDA Reg. CODE	
		50ml x 2cartridges, w/mixing tips	VF310-SLN0	Vonflex S Light 0102	
	Normal	50ml x 12 cartridges	VF360-SLN0	Vonflex S Light 0112	
Vonflex S Light		50ml x 48 cartridges	VF324-SLN0	Vonflex S Light 0148	
	Fast	50ml x 2cartridges, w/mixing tips	VF310-SLF0	Vonflex S Light Q 0102	
		50ml x 12 cartridges	VF360-SLF0	Vonflex S Light Q 0112	
		50ml x 48 cartridges	VF324-SLF0	Vonflex S Light Q 0148	

ALL PRODUCTS			PRODUCT CODE	KFDA Reg. CODE	
Vonflex S Light XLV	Normal	50ml x 2cartridges, w/mixing tips	VF310-SLNX	Vonflex S Light XLV 1402	
		50ml x 12cartridges	VF360-SLNX	Vonflex S Light XLV 1412	
		50ml x 48cartridges	VF324-SLNX	Vonflex S Light XLV 1448	

VERILCOM
Variage to Mode foods

Dental Materials | VPS Impressions



# **Excellent Fixation Characteristic Material**



## FEATURES AND BENEFITS


- Vonflex S Supermix is an addition polymerization silicone composed of hydrophilic vinyl polysiloxane which helps to create precise shape of oral tissue. which meets the requirements of ISO 4823.
- Void-free and uniform mixing when used with the automatic mixer
- High tear resistance
- Superior initial hydrophilicity
- Superior dimensional stability
- Easy to take off the impression body from the patient's oral because of its excellent elastic recovery
- Extended working time

### **INDICATIONS**

- Taking impression for crown and bridge
- Taking impression for inlay and onlay


# TECHNICAL DATA

Contact angle




Good flowability and Accurate impression on the oral surface where saliva is present.





Snap set technology(optimal working time and short setting time) helps to make better impression results.

## • Tear strength



## **PACKAGE**

ALL PRODUCTS			PRODUCT CODE (KFDA Reg. CODE)
Vonflex S Supermix		380ml x 1 cartridge, Mixing tip 10ea, Fixation ring 1ea	VF3SK-SSNE
	Normal (Heavy type)	380ml x 2 cartridges	VF338-SSNE
		380ml x 4 cartridges	VF34B-SSNE

34 www.vericom.co.kr



# **Excellent Fixation Characteristic Material**



## FEATURES AND BENEFITS


- Vonflex S Supermix MONOPHASE is an addition polymerization silicone composed of hydrophilic vinyl polysiloxane which helps to create precise shape of oral tissue. which meets the requirements of ISO 4823.
- Void-free and uniform mixing when used with the automatic mixer
- High tear resistance
- Superior initial hydrophilicity
- Superior dimensional stability
- Easy to take off the impression body from the patient's oral because of its excellent elastic recovery
- Extended working time

### **INDICATIONS**

- Taking impression for crown and bridge
- Taking impression for inlay and onlay


# TECHNICAL DATA






Good flowability and Accurate impression on the oral surface where saliva is present.





 $\label{thm:constraint} Snap\,set\,technology (optimal\,working\,time\,and\,short\,setting\,time) \\ helps\,to\,make \quad \begin{array}{c} better \quad impression\,results. \end{array}$ 

#### • Tear strength



## **PACKAGE**

ALL PRODUCTS			PRODUCT CODE (KFDA Reg. CODE)
Vonflex S Supermix MONOPHASE		380ml x 1cartridge, Mixing tip 10ea, Fixation ring 1ea	VF3SK-SMNI
	Normal (Medium type)	380ml x 2cartridges	VF32B-SMNI
		380ml x 4cartridges	VF34B-SMNI

VERILCOM
35
36

Dental Materials | VPS Impressions



# **Excellent Fixation Characteristic Material**

[Type 1]

Tube Type

Vonfex

Von

# FEATURES AND BENEFITS

- Vonflex S Heavy is an addition polymerization silicone composed of hydrophilic vinyl polysiloxane which helps to create precise shape of oral tissue. which meets the requirements of ISO 4823.
- Vonflex S Heavy is convenient to use and has low deformation, helps to make precise impression taking.
- Dual phase impression technique
- Precise duplication of models
- All impression techniques that operators need heavy viscosity materials

# **INDICATIONS**

- Taking impression for crown and bridge
- Taking impression for inlay and onlay

# TECHNICAL DATA

• Easy handling property





# **PACKAGE**

ALL PRODUCTS			RODUCTS	PRODUCT CODE	KFDA Reg. CODE
	Normal	Injection	50ml x 2cartridges, w/mixing tips	VF310-SHN0	Vonflex S Heavy 1602
		Injection	50ml x 12 cartridges	VF360-SHN0	Vonflex S Heavy 1612
		Injection	50ml x 48cartridges	VF324-SHN0	Vonflex S Heavy 1648
Vonflex S		Tube	Base 150ml & Catalyst 150ml x 2ea	VF330-SHNT	Vonflex S Heavy 16T
Heavy	Fast	Injection	50ml x 2cartridges, w/mixing tips	VF310-SHF0	Vonflex S Heavy Q 1602
		Injection	50ml x 12 cartridges	VF360-SHF0	Vonflex S Heavy Q 1612
		Injection	50ml x 48cartridges	VF324-SHF0	Vonflex S Heavy Q 1648
		Tube	Base 150ml & Catalyst 150ml x 2ea	VF330-SHFT	Vonflex S Heavy Q 16T

36 www.vericom.co.kr



# Proven One Step Material Based on A-silicone

[Type 2] Medium-bodied consistency



### FEATURES AND BENEFITS

- Vonflex S Medium is an addition polymerization silicone composed of hydrophilic vinyl polysiloxane which helps to create precise shape of oral tissue. which meets the requirements of ISO 4823.
- Vonflex S Medium is convenient to use and has ideal consistency for monophase technique, helps to make precise impression taking.

#### INDICATIONS

- Taking impression for crown and bridge
- Taking impression for inlay and onlay

## **PACKAGE**

ALL PRODUCTS			PRODUCT CODE	KFDA Reg. CODE
Vonflex S Medium	Normal	50ml x 2cartridges, w/mixing tips	VF310-SMNO	Vonflex S Medium 2502

Dental Materials | Bite Registration



Thixotropic Bite Registration Material



# FEATURES AND BENEFITS

- Vonflex S Bite, as the additional polymerization silicone type, is composed of vinyl polysiloxane impression materials designed for making precise occlusal registrations.
- Vonflex S Bite provides excellent detail for accurate impressions.
- Standard bite registration in the end bite position
- Portrayal of antagonists
- Coating of bite forks
- Key material for needle point registration
- Production of small model segments

#### **INDICATIONS**

Bite registration

### **PACKAGE**

	ALL PRODUCTS	PRODUCT CODE	KFDA Reg. CODE
Vonflex S Bite	50ml x 2cartridges, w/mixing tips	VF310-SBNB	Vonflex S Bite 0202

VERICON
37
Wang to Roof total