

# **BioSonic**® UC300 / UC300R Ultrasonic Cleaning System

### Unmatched flexibility

Multiple basket and beaker arrangements
Adjustable basket dividers for instrument sorting
Accommodates most racks to fit popular cassette systems

### Unique hinge action cover

Allows cover to stay fully open Reduces noise levels for a quiet operation Traps aerosols within and directs condensation back to the tank

### Recessed model available

Optional remote digital timer for more user friendly operation





# BioSonic® UC300/UC300R Ultrasonic Cleaning System

Tank Capacity: 3.8 gallons (14.4 liters) Use Capacity: 3.0 gallons (11.4 liters) Overall Dimensions: 19 x 12 x 12.25 inches (48.26 x 30.48 x 31.12 cm)

14 x 9 x 6.75 inches

Tank Dimensions: (35.6 x 22.86 x 17.15 cm)

Output: 300 watts Operating Frequency: 40 kHz



Above countertop surface: 17.125 x 12.125 x 2.25 inches (43.5 x 30.8 x 5.7 cm) Below countertop surface: 14.25 x 10.5 x 14 inches (36.2 x 26.7 x 35.6 cm)

For maximum performance and longevity of your BioSonic Ultrasonic Cleaner, we recommend the use of genuine BioSonic solutions and accessories.







## **BioSonic** Solution Concentrates in MeterDose® bottles











Solution: IIC30 IIC34 11C38/11C42 11039 BioSonic BioSonic BioSonic **BioSonic** BioSonic **Enzymatic Ultrasonic Plaster and Stone** Germicidal Ultrasonic **Super Concentrate Ultrasonic Cement General Purpose Cleaning Concentrate** Remover Powder **Cleaning Concentrate** Remover Concentrate Cleaning Solution Concentrate 8 oz. (236 ml) **Contents:** 16 oz. (473 ml) 8 oz. (236 ml) 1.85 lb. (840 g) UC38 16 oz. (473 ml) 16 oz. (473 ml) UC42 **Diluted Yield:** 8 gal. (30 liters) 8 gal. (30 liters) 3 gal. (11.4 liters) 8 gal. (30 liters) UC38 .4 gal. (1.5 liters) 16 gal. (60 liters) UC42 **Dilution Ratio:** 1:128 1:135 **Contaminants:** Blood, debris, saliva and Protein soils (dried blood, Plaster and stone Blood, debris, saliva and Temporary cements general cleaning tissue, saliva) and gengeneral cleaning. Cleans and permanent cements eral cleaning. Effectively instruments as it disin-(including zinc-phosphate, breaks down and loosens fects the ultrasonic bath. polycarboxylate) hardened debris. Items: Instruments, crowns, Crowns, veneers, brack-Instruments, crowns, Instruments, crowns, Dentures, partials, castings, denture flasks, bridges, dentures and bridges, dentures and bridges (no plastic parts) ets and instruments ultrasonic scaler tips ultrasonic scaler tips mixing bowls and nonaluminum articulators

### **Order Information**

UC300	BioSonic Ultrasonic Cleaner Accessories not included	UC315115	Remote Digital Timer for UC300R
	Accessories not included	UC325	Cassette Basket for UC300 / UC300R
UC300R	BioSonic Recessed Ultrasonic Cleaner Accessories not included	UC330	Two-Beaker Holder
UC310	Large FingerGuard™ Instrument Basket	UC53	Borosilicate Beaker 600 ml with cover and positioning ring
UC311	Twin FingerGuard™ Instrument Basket	LICORE	Basket Section Inserts
UCB2	Small Parts Basket	UC335	

### Coltène/Whaledent Inc.

235 Ascot Parkway Cuyahoga Falls, OH 44223 / USA Tel. USA & Canada +1 800 221 3046 +1 330 645 8704

Fax +1 330 916 7077 info.us@coltene.com

#### Coltène/Whaledent AG

Feldwiesenstrasse 20 9450 Altstätten / Switzerland Tel. +41 (0)71 757 53 00 Fax +41 (0)71 757 53 01 info.ch@coltene.com

#### EC REP

Coltène/Whaledent GmbH + Co. KG Raiffeisenstrasse 30

89129 Langenau / Germany Tel. +49 (0)7345 805 0 Fax +49 (0)7345 805 201 info.de@coltene.com





coltene.com mybiosonic.com