

## VEHICLE ACCESS & MODIFICATION FOR DISABLED & AGED PERSONS

**VEHICLE ACCESS  
FOR PERSONS**



**TRANSFERRING YOUR  
WHEELCHAIR OR SCOOTER**







**9%** of Australians  
have a disability with  
Moderate to Profound core  
activity limitations.\*

**12%** of Australians are  
involved in caring for a  
person with a disability  
or an aged person.\*

Almost 1 million carers  
are primary carers.\*

\*ABS; Disability, ageing and carers 4430.0 2009



Accessing a vehicle helps to restore  
the independence of a disabled  
or elderly person.

World renowned manufacturers in  
vehicle access solutions:

**autochair**  
Innovating Vehicle Adaptations since 1983

**Autochair** is the leading  
supplier "in vehicle" transfer  
of persons, wheelchair &  
scooters with over 45,000  
satisfied users worldwide.  
Manufactured in UK.



**AMF-Bruns** is the leading  
supplier of vehicle  
modifications. Over 4,000  
vehicles are modified  
annually from their factory in  
Apen, Germany.



*Compliant with International & Australian Safety Standards*



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



# SMART TRANSFER

“ The guys who fitted the Smart Transfer person lift were absolutely fantastic. I could not have wished for better or friendlier guys. ”

Ms O'Malley, Penicuik

Simple 2 button operation

## Why choose The Smart Transfer?



Move straight from chair into vehicle



Works with standard car seats in any make or model



Simple to install, assemble and operate



No extra lifting required



Can lift up to 150kgs



Easily transferred to your next vehicle



Standard installation takes just 3-4 hours



Manufactured in UK

**45,000**  
satisfied  
users  
worldwide



**autochair**  
Innovating Vehicle Adaptations since 1983



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



# SMARTLIFTER series

## Get yourself mobile with The Smart Lifter

Designed to lift scooters and wheelchairs weighing up to 200kg. Straightforward to fit, suitable for a wide range of vehicles and easily detachable when not in use, the Smart Lifter makes light work of loading and unloading.

**autochair**  
Innovating Vehicle Adaptations since 1983



Manufactured in the UK



Easy to use



Lifts all types of scooters and wheelchairs



Lifts up to 40, 80, 100, 150, 200kgs



Over 300 vehicle specific installation kits



Transfer easily to your next vehicle



Leaves space in boot when not in use



## Lift & loads in just 60 seconds!

**SMARTLIFTER LM**



**45,000**  
satisfied  
users  
worldwide



Strap reaches and lifts low loads



Boom arm can fold



Maximum Safe Working load of  
40kg & 80kg



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



## SMARTLIFTER LC



Quick release  
removable  
spreader bar



48% narrower boom



Increase lifting range by 200mm



Quick release spreader bar



Maximum Safe Working load of 80kg & 100kg

The NEW  
Compact LC is  
ultra-light and slim

## SMARTLIFTER LP



Simple 4  
button  
operation

The strongest  
Smart Lifter for loads  
of up to 200kg




Maximum Safe Working load of  
150kg & 200kg



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



# AUTOCHAIR EXCLUSIVE LIFTING ATTACHMENTS



**A.** Quick release removable spreader bar



**B.** Spreader Bar with keyhole attachments or karibiners.


Suitable for scooters with seat posts.



**C.** Halo seat post lifter attachment (75kg capacity).

Suitable for lifting a scooter with removable battery.

“ A lifter is useless unless you have the correct attachments to lift your wheelchair, powerchair or scooter. Autochair has them! ”



**D.** Wheelchair lifter attachment (50kg capacity) for smart lifter lc compact

Spreader bar with 4 karibiners to lift a powerchair.

Available options for 80kg - 200kg



The seatbelt restraint ensures your wheelchair or scooter is safely stowed



MANY MORE OPTIONS AVAILABLE



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



# SMARTBASE

“ After struggling for ages to manoeuvre my husband from his chair to his wheelchair, we discovered the SmartBase... a really big help around the home ”



## Why choose The SmartBase?

-  Adjustable leg width to suit space available
-  Simple to assemble and operate
-  No extra lifting required
-  Can lift up to 150kgs
-  Easily manoeuvred in and around the home
-  Folds flat for ease of use at home, work or on holiday
-  Manufactured in UK

“ The SmartBase works with our best-selling person lift, the Smart Transfer. ”

**autochair**  
Innovating Vehicle Adaptations since 1983



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



# PANORAMA LINEAR LIFT



## Why choose The Panorama Linear Lift?



Patented, newly developed aluminium arms make the AL1 the lightest lift on the market



See-through platform design ensures excellent rear view.



The fully automatic lift for rear and side entry



Roll-off stop opens automatically when touching ground even on uneven surfaces



900mm lifting height



Other Linear Lift design options available



Manufactured in Germany



Lifting capacity of 400kg



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



# CASSETTE LIFT K90



Ultra slim view when not in use



Large lifting platform with barriers

## Why choose The Cassette Lift K90?



Ultra-slim design guarantees ultimate ground clearance



Large lifting platform, suitable for large electric wheelchairs up to 300kg



Vehicle interior remains unchanged



No changes to the car body necessary



Manual emergency operation



Manufactured in Germany



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



# FUTURESAFE HEAD & BACKREST



## “ Safe vehicle restraints solutions

### Why choose FutureSafe?



Only AMF Bruns possesses patented 20G front impact & 10G rear impact head and back supports



Shoulder belt increases safety in case of front impact



Height adjustable to your individual needs



Easily swivelled manually



Floor or sidewall mounted



Manufactured in Germany



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)





# for person and wheelchair.

## PROTEKTOR SYSTEM RESTRAINTS - SILVER SERIES

Meets ISO 10542-1 + 2, DIN 75078-2, CSA Z 605, FMVSS/CNVSS, **AS 2596**, SAE J 2249, 2001/85/EC.



Find Us [www.mobilitycare.net.au](http://www.mobilitycare.net.au)



# VEHICLE ACCESS CHECKLIST

- ☐ Do you require vehicle mobility for a disabled or elderly person?
- ☐ Do you need to transport your wheelchair, powerchair or scooter?
- ☐ Does the carer need to transfer a person safely and easily?
- ☐ Is patient comfort and safety when driving a priority?
- ☐ Is ease of access vital at the beginning and end of the journey?
- ☐ Is a discrete and dignified solution important to you?
- ☐ Do you require a comprehensive vehicle modification?
- ☐ Is value for money and worlds best quality important?
- ☐ Would you prefer a solution that can easily be transfered to another vehicle?

If you TICKED any of the above questions please contact us about a free demonstration and quote.

For more information and video demonstrations visit our website.

We provide quotes for you and your therapist in a format that may assist you with state or federal funding.



YOUR LOCAL DEALER:

CONTACT US

**MOBILITYCARE**  
*Since 1992*

Phone: (03) 9568 8383  
39 Cambro Road, Clayton VIC 3168  
[info@mobilitycare.net.au](mailto:info@mobilitycare.net.au)  
[www.mobilitycare.net.au](http://www.mobilitycare.net.au)