



INTUBATION IS EASY — UNTIL IT ISN'T.

McGRATH® MAC
Enhanced Direct Laryngoscope

Medtronic


IMPROVE YOUR VIEW WITH THE McGRATH® MAC VIDEO LARYNGOSCOPE

THE ANTERIOR CAMERA IS DESIGNED TO

- Improve grade of view compared with direct laryngoscopy
- Provide better visual confirmation of tube placement
- Reduce tunnel vision experienced with direct laryngoscopy

MINIMAL CHANGE IN INTUBATION TECHNIQUE

- Macintosh type blade
- Used like a traditional laryngoscope



ENHANCE YOUR CHANCE OF FIRST-ATTEMPT SUCCESS.¹

You prepare for every eventuality before you intubate.

So when your patient is on the table, when those critical moments are passing, you are better prepared for whatever happens.

The McGRATH[®] MAC enhanced direct laryngoscope (EDL) is a single handheld device that combines the benefits of:

- Direct laryngoscopy (DL) — the familiar technique of a Macintosh laryngoscope
- Video laryngoscopy (VL) — video technology

THE McGRATH[®] MAC ENHANCED DIRECT LARYNGOSCOPE

INTUITIVE

Uses standard DL technique combined with VL technology

PATIENT-FRIENDLY

Offers improved access to the patient's airway with a slim blade profile²

ON-SCREEN COUNTDOWN

Includes a 250-minute battery with a minute-by-minute on-screen charge countdown that confirms device battery life

DURABLE

Drop tested to two meters, the McGRATH[®] MAC EDL can endure rough handling and still perform in a fast-paced medical environment

EASY TO CLEAN

Provides sterile integrity with disposable blades and is fully immersible. Both screen and handle can be treated with high-level disinfection



ORDERING INFORMATION

CODE	DESCRIPTION	SIZE	QUANTITY
300-000-000	Handle	N/A	1
350-072-000	McGRATH® MAC blade size 1	1	50/case
350-017-000	McGRATH® MAC blade size 2	2	50/case
350-005-000	McGRATH® MAC blade size 3	3	50/case
350-013-000	McGRATH® MAC blade size 4	4	50/case
X3-003-000	McGRATH® MAC X blade™ size 3	3	10/case
340-000-000	250-minute battery pack	N/A	1

MIXED BLADE KITS:

AUSTRALIA ORDERING CODE	NEW ZEALAND ORDERING CODE	KIT DESCRIPTION	QTY SIZE 1 BLADES	QTY SIZE 2 BLADES	QTY SIZE 3 BLADES	QTY SIZE 4 BLADES	QTY SIZE X BLADES
K000MB	K800MB	Adult Blade Kit	0	0	10	10	5
K001MB	K801MB	General Blade Kit	5	5	10	10	5
K002MB	K802MB	Paediatric Blade Kit	10	10	0	0	0

TECHNICAL SPECIFICATIONS

McGRATH® MAC Handle

Size: 180 mm x 68 mm x 110 mm

Weight: 200 g

Power: Proprietary 3.6V lithium battery pack (c. 250 minutes)

Light source: High-intensity LED

Display: 2.5" LCD color display

Camera: CMOS

Materials: Durable medical-grade thermoplastics with reinforced structural alloy core. The device and packaging are latex-free.

McGRATH® MAC Disposable Blades

Material: Fog-free medical-grade optical polymer

Packaging: Packaged sterile for single use

For a demonstration or quote, please contact your local Medtronic representative.

REFERENCES

1. Su YC, Chen CC, Lee YK, Lee JY, Lin KJ. Comparison of video laryngoscopes with direct laryngoscopy for tracheal intubation: a meta-analysis of randomised trials. *Eur J Anaesthesiol.* 2011;28(11):788-795.
2. Levitan R, Heitz J, Sweeney M, Cooper R. The complexities of tracheal intubation with direct laryngoscopy and alternative intubation devices. *Ann Emerg Med.* 2011;57 (3):240-247.

Medtronic

Medtronic Australasia Pty Ltd
2 Alma Road
Macquarie Park NSW 2113
Australia
Tel: +61 2 9857 9000
Fax: +61 2 9889 5167
Toll Free: 1800 668 670

Medtronic New Zealand Ltd
Level 3 - Building 5, Central Park Corporate Centre
666 Great South Road,
Penrose, Auckland 1051
New Zealand
Fax: +64 9 918 3742
Toll Free: 0800 377 807

medtronic.com.au
medtronic.co.nz

© Medtronic.
2021 All Rights Reserved. Printed in Australia.
Ccare 379-01-21 ANZ. #9131-022021