


Solero

Microwave Tissue Ablation System



SIMPLE | SPEEDY | SCALABLE



SOLERO MICROWAVE TISSUE ABLATION SYSTEM

The Solero* Microwave Tissue Ablation (MTA) System and Accessories are indicated for the ablation of soft tissue during open, laparoscopic, or percutaneous procedures. The Solero MTA System is not intended for cardiac use.

SYSTEM SPECIFICATIONS

- 2.45 GHz Operating Frequency
- Solid State Generator
- Powers up to 140 W
- Intuitive Touch Screen User Interface
- Integrated Peristaltic Pump
- Applicator Coolant Temperature Monitoring
- Reflected Energy Monitoring
- No Grounding/ Dispersive/ Neutral Electrode Required
- Single Applicator System
- Applicator Available in 14 cm, 19 cm, and 29 cm Lengths



SKU	Description
H787127400000	Solero MW Generator
H7877001060010	Solero MW Applicator 14 cm
H7877001060020	Solero MW Applicator 19 cm
H7877001060030	Solero MW Applicator 29 cm

The Solero Microwave Generator and Accessories are not FDA cleared for sale in the United States. Please refer to the Solero Generator Operator's Manual and the Solero Applicator Directions For Use for complete instructions, warnings, and precautions.

DISTRIBUTED BY



AUSTRALIA

GENERAL ENQUIRIES: 1300 886 385

FAX: 1300 668 440

EMAIL: info.au@getzhealthcare.com

NEW ZEALAND

GENERAL ENQUIRIES: 0800 733 333

FAX: +6 49 447 2001

EMAIL: orders.nz@getzhealthcare.com

www.getzhealthcare.com



USA > 14 Plaza Drive, Latham, NY 12110 > tel: 800-772-6446 or 518-798-1215 > fax: 518-798-1360
International > Haaksbergweg 75 (Margrietoren), 1101 BR, Amsterdam Z-O > The Netherlands
tel: +31 (0)20 753 2949 > fax: +31 (0)20 753 2939

www.angiodynamics.com

*AngioDynamics, the AngioDynamics logo, Solero and the Solero logo are trademarks and/or registered trademarks of AngioDynamics, Inc., an affiliate or a subsidiary. © 2016 AngioDynamics, Inc. ANGB 358 GL Rev 01 12/16