

Simple Scan, Clear Result

EdgeFlow UH10

The intelligent bladder scanner




Deep learning technology for simplicity and accuracy

AI

Experience LiveContour™ technology. It automatically identifies the contour of the bladder in real-time ultrasound images.

Just sweep the probe over a patient's abdomen. MaxContour™ allows you to obtain the optimized images from a series acquired during scanning.




LiveContour™

Experience DualActive™ probe, revolutionary imaging technology

DualActive™ all-electronic probe shows you real-time bladder images of sagittal and transverse planes simultaneously in live-imaging.

Powered by DualActive™ probe technology and deep learning technology, EdgeFlow UH10 facilitates seeing into the whole bladder clearly with confidence and makes your exam swift and easy.

• DualActive™ probe



• **Ergonomic Design**

Lightweight and easy-grip design
Intuitive user interface



• **Convenient Cradle Station**

Dual-function of printing & charging

• **DualActive™ probe**

T-shaped array presents images to allow simultaneous observation of the sagittal and transverse planes of the bladder



Lithium-ion battery
Charging time : 3 hours 50 mins
More than 400 scans when fully charged
Dual charging system : C-type adapter/Cradle station



Accuracy : $\pm 7.5\%$ (over 100ml)
Volume measurement range : 0 - 999ml



Weight : Console & probe 419g



Electronic probe without periodic calibration



IPX4 level probe protection

*Accuracy is based on the recommended user guideline with a tissue-equivalent phantom.

Edgecare

Copyright © 2024 Edgecare Inc. All rights reserved. Edgecare, Edgecare logo, EdgeFlow, DualActive, LiveContour, and MaxContour are registered and unregistered trademarks of Edgecare Inc.

www.edgecare.co.kr