



Patient **Support** Systems

MEDICAL • DESIGN • INNOVATION

Air Chair *Active*

Our unique active alternating air cell system provides the ultimate in care and comfort.

Featuring a semi-automatic cushion system with battery backup for up to 8 hours.



Unit C3, 15 Narabang Way Belrose, Sydney NSW 2085
t: 02 9844 5456 f: 02 9844 5445 e: info@patientsupportsystems.com

www.patientsupportsystems.com



Patient **Support** Systems

MEDICAL • DESIGN • INNOVATION

Product Specs



Air Chair Active

Order code: C01-BLU-ALT

Safe working load: 150kg

Overall height: 1300mm

Overall width: 710mm

Overall depth: 560mm

Seat level height: 700mm

Seat width: 510mm

Seat depth: 560mm

Seat back adjustment angle: 97-136 °

Seat tilt adjustment angle: 0-17 °

Leg rest adjustment angle: 0-90 °

Included accessories:

Active Alternating Air Cell System, Carilex Pump,
Meal tray, Headrest, Footrest Pads

Unit weight: 50kg

Package dimensions: 1150x700x1300mm



Features

Air Chair *Active* uses an advanced Carilex active air cell system throughout the seating surface to deliver the best possible pressure care. This system features intelligent pressure sensing to adjust cell pressure to patient movements, as well as customisable comfort settings. This is in addition to our foam system and 4-way stretch material. The chair is highly adjustable, to cater for a range of users and requirements.

- Advanced Active Air Cell System with IPS (Intelligent Pressure Sensing) Technology
- Customisable comfort settings - digitally adjust internal air cell pressure for patient comfort
- Battery back up - up to 8 hours battery backup, for use in patient transportation or in the event of power failure
- Whole body therapy - even in sitting position
- Adjustable therapy modes - offers dynamic and static therapy, with multiple cycle times
- Special 4-way stretch material covers user area and provides additional pressure relief and mitigates skin tear risks
- No-gap seating surface
- Water resistant to prevent liquid ingress, fire retardant
- Easy to remove covers for cleaning and replacement
- Completely adjustable for patient positioning, with ergonomic controls for caregiver, including recline, tilt in space, leg rest adjustment
- Head rest is height adjustable and removeable
- Multi-position fold away arm rests, wings and footrest pads provide ease of access and transfer
- Individual locking castors with single rear directional lock castor

Unit C3, 15 Narabang Way Belrose, Sydney NSW 2085
t: 02 9844 5456 f: 02 9844 5445 e: info@patientsupportsystems.com

www.patientsupportsystems.com