

TheGlobalBeautyGroup Technology for Beautiful Business





Welcome to a new world of medi-aesthetics!





06





80







16

20

22









TWIST: MEDICAL-GRADE MICRONEEDLING FROM GERMAN NEEDLING INNOVATORS

HYDROLUX: ADVANCED SKIN RESURFACING + PRODUCT INFUSION

HAIRLASE: INTUITIVE OPERATION FOR CONSISTENT HAIR REMOVAL RESULTS

MEDILUX: 5 WAVELENGTHS OF MEDICAL-GRADE LED

ULTRALUX + ULTRALUX PRO: 12 IN-DEMAND TREATMENTS IN 1 DEVICE

SUPERLUX + SUPERLUX PRO: AWARD-WINNING MULTIFUNCTIONALITY

POWERLUX: VERSATILE IPL MEETS POWER, PERFORMANCE AND VALUE

ULTRALASE: ADVANCED Q-SWITCHED ND:YAG LASER TECHNOLOGY

SLIMLUX: SLIM, TONE, TIGHTEN AND LIFT WITH ULTRASOUND + RF

TECHNOLOGY FOR BEAUTIFUL BUSINESS

TheGlobalBeautyGroup Technology for Beautiful Business





REGENERATE

ANALYSIS





REBORN: MEDICAL-GRADE FAT REDUCTION IN A WHOLE NEW LIGHT

WITH TECHNOLOGY ENHANCED SKIN

OBSERV 520X: GO BENEATH THE SURFACE





THUZZLE: SOLVE 'THE THERMIC PUZZLE' WITH THE EVOLUTION OF THERMAL ENERGIES

PLEXR PLUS: THE NEXT GENERATION OF MEDICAL-GRADE PLASMA TECHNOLOGY





THE GLOBAL BEAUTY GROUP **EDUCATION HUB**



Innovating the industry since 2010, The Global Beauty Group are driven by a mission to make people feel wonderful.

Our focus on delivering clinics with forward-thinking treatment innovations and an unparalleled level of service has supported over 4,500 beauty businesses succeed with advanced aesthetic and medical-grade technology.











Award-Winning Excellence

The Global Beauty Group and our range of technology have been awarded some of the industry's highest honours, including MyFaceMyBody Supplier of the Year, the Australian Beauty Industry Award for National Wholesaler of the Year and the Dermascope Aesthetician's Choice Award for Favourite Multi-Functional Equipment.

Industry Leading Support

Your support team at The Global Beauty Group includes gualified Dermal Therapists, Bio-Medical Engineers, Customer Service Specialists, Marketing Professionals, Clinical Education Experts and experienced Account Managers. From rapid response technical and clinical support 7 days a week, comprehensive marketing resources, online educational programs with Education HUB and face-to-face education programs, our clients enjoy ongoing support that leads the industry.

Best-In-Class Technology

In addition to delivering our own LUX Series range of award-winning technology, The Global Beauty Group are proud to partner with international aesthetic and medicalgrade leaders to deliver best-in-class innovations from around the globe. With Medical CE, TGA listing, ISO Certification and FDA Approved technologies, the range delivers trusted quality, safety and reliability.



SLIM, TONE, TIGHTEN AND LIFT WITH ULTRASOUND + RF

Take advantage of the rapidly growing demand for face and body enhancement services with the versatile SlimLUX:

1. ULTRASOUND BODY CONTOURING

High-frequency ultrasound waves are utilised to contour and sculpt targeted areas. A course of 5-10 treatments is typically required for optimum results, with occasional maintenance sessions recommended to support long-term body goals and generate ongoing clinic bookings.

2. RADIO FREQUENCY (RF) SKIN TIGHTENING

RF uses controlled heat to target the dermal layer, generate new collagen and improve skin laxity. Services are suitable for almost all areas of the face and body including full face skin tightening, Brazilian butt lifts, eye lifts, arm tone-ups, neck lifting, and postnatal tummy toning.

3. SONOPHORESIS SKIN CARE INFUSION

Ultrasound skin care infusion facials result in greater absorption of topically applied skin care. Treatments plump and hydrate fine lines and wrinkles, as well as assist cellulite management product absorption.















An effective, non-surgical and non-invasive alternative to procedures such as liposuction and skin lifting.

Services can be performed from as little as one session.

Body Confidence for Every Body

demographics.

offer treatments to:

- Contour, shape and tone the body
- Tighten lax skin on the face and body
- Reduce the appearance of cellulite
- Infuse skin care with sonophoresis

Treatments are comfortable, require no downtime and allow clients to get on with daily activities as usual.

SlimLUX delivers a comprehensive range of technology enhanced treatments to sculpt, tone, contour, tighten and lift. Treatments are highly versatile, and target a wide range of client concerns and

With SlimLUX you can comfortably, non-surgically and non-invasively

- Reduce the appearance of fine lines and wrinkles



ADVANCED Q-SWITCHED ND:YAG LASER TECHNOLOGY

UltraLASE is an advanced Q-Switched Nd: YAG laser designed for clinics and practices wanting to meet the growing demand for professional and reliable Laser Tattoo Removal plus Carbon Facial services. Treatments are non-invasive and work via the delivery of laser light energy to the targeted pigment in the treatment area.

KEY FEATURES

- Suitable for a broad range of skin types and tattoo pigmentation colours
- Both 532nm and 1064nm treatment wavelengths are available
- · Perform treatments easily with preprogrammed options plus tailorable manual modes
- User-friendly, with an intuitive user interface and seamless switching between wavelengths
- 6-12 sessions are typically required to achieve optimal results, depending on tattoo quality
- On-trend Carbon Facials add additional treatment revenue options











A 2018 McCrindle Research study showed that while 20% of people have one or more tattoos, 31% of inked individuals regret at least one of their tattoos.



Carbon Facials are an advanced skin rejuvenation treatment also commonly known as the "China Doll Facial" or "Hollywood Facial".

By combining cutting-edge LUX Series technology with a special carbon cream, Carbon Facials work to boost collagen and elastin production. Treatments refine skin imperfections such as enlarged pores, pigmentation, acne scarring and fine lines, while offering an instantaneous glow to your client's complexion.

Carbon Facials are the perfect skin corrective treatment for clients wishing to achieve a smoother, clearer, more radiant skin tone without having to spend too much time in the treatment room. Treatments are comfortable, fast and require no downtime.



IDENTIFY Ink in tattooed areas creates pigment that sits in the dermal layer of the skin.



TARGET Treatments target pigment in unwanted tattoos with Q-Switched Nd: YAG laser energy.



SHATTER

The tattoo then absorbs the energy, fraamentina the ink into fine particles.



REMOVE

Shattered ink is processed by the body's natural filtration (lymphatic) system.





VERSATILE IPL MEETS POWER, **PERFORMANCE AND VALUE**

Ready to power up your business with high-performance IPL services? PowerLUX offers medical-grade IPL (ARTG #295158) that provides power, results and versatility in one tabletop device.

Compact in size, the PowerLUX fits neatly on a table top or treatment trolley and is perfect for even the most spacelimited treatment rooms. As well as a streamlined design, the PowerLUX combines 1,200 watts of power, a variety of spot sizes and an integrated skin cooling system (reaching as low as -4°) to deliver all top 6 IPL treatments.

KEY FEATURES

- 6 in-demand IPL treatments plus optional Radio Frequency upgrade
- Integrated skin contact cooling from -7°c
- Pre-programmed settings based on Fitzpatrick Skin Type
- Manual modes for clinicians wanting greater customisation
- Interchangeable magnetic safety filter system
- 15x35mm, 15x50mm and 12x40mm spot sizes available
- Large 8.4 inch full colour touch screen with intuitive interface



IPL Photo Rejuvenation



IPL Pigmentation Correction

Correction causes sun damage (such as freckles) to react and darken. The targeted pigment pigmentation, sun damage and a more uniform skin tone with clients typically seeing results within 2-4



IPL Vascular Therapy

530nm and 560nm wavelength IPL energy is absorbed into the chromophore (haemoglobin) of targeted redness during LUX Series IPL Vascular Therapy sessions, causing damage to the unwanted vessel wall and leading it to collapse then disappear. Treatments are highly effective, with vascular lesions such as

telangiectasia and cherry angiomas typically resolved with 2-4 sessions.



IPL Acne Management



IPL Hair Removal

target melanin within hair follicles. Energy is absorbed into the bulb of the hair and damages the dermal degrees helps deliver effective yet



LUX Series IPL Technology

The medical-grade, clinically proven range of high-performance LUX Series IPL treatments has been designed by The Global Beauty comfort plus ease-of-use front of mind. LUX Series IPL technology is available in a range of devices intended to make industry leading

From the compact PowerLUX to the multifunctional SuperLUX with UltraLUX PRO, LUX Series IPL technology delivers exceptional results for your clients and clinic.





Dual-Mode Hair Removal

LUX Series Dual-Mode combines 3 advanced technologies: Bipolar Radio Frequency (RF) + Intense Pulsed Light (IPL) + Skin Contact Cooling. RF allows IPL energy to travel deeper into the skin to target hair with less energy, increasing treatment comfort. Unlike some other IPL devices, Dual-Mode technology can achieve great results on even soft and thin hair plus shades of blonde and red hair.



SUPERLUX + SUPERLUX PRO: **AWARD-WINNING MULTIFUNCTIONALITY**

SuperLUX is the next generation of multifunctional, medicalgrade IPL (ARTG #295158) for advanced skin rejuvenation results and faster, pain free permanent hair removal.

The 2017 Dermascope Aesthetician's Choice Award for the best multifunction equipment, SuperLUX boasts state-of-the-art SHR (Super Hair Removal) technology.

In addition to all 6 PowerLUX IPL applications, SuperLUX adds the rapid and comfortable next-generation hair removal technology SHR. Utilising high pulse repetition rates at low energy levels, SHR scans the treatment area to gradually heat and destroy the hair follicle. The result is pain free permanent hair removal treatments in half the time of a traditional IPL hair removal session.



After

THE SCIENCE BEHIND SHR IPL HAIR REMOVAL

Comfortable, permanent hair removal results on larger areas in half the time of standard IPL.



and gradually heats up

7 IPL **APPLICATIONS**

IPL Acne Management **IPL** Pigmentation Correction IPL Vascular Therapy IPL Photo Rejuvenation IPL Hair Removal Dual-Mode Hair Removal SHR IPL Hair Removal

TREATMENTS Maximise your treatment potential with optional LUX

Tightening, or upgrade to SuperLUX PRO with Nd:YAG laser and LUX Series Tattoo

KEY FEATURES

- Powerful and multi-functional device with 2,000 watts of power
- Intuitive graphic user interface with 10.4 inch full colour touch screen
- Pre-programmed settings based on Fitzpatrick Type, plus manual mode options
- Features 4 powerful, in-demand and clinically proven IPL skin correction treatments
- Features 3 advanced hair removal treatments that target a wide range of hair and skin types
- Optional upgrades to Tattoo Removal, Carbon Facials and RF Skin Tightening
- Interchangeable magnetic safety filter system plus 15x35mm, 15x50mm and 12x40mm spot size options
- Integrated skin contact cooling from -7°c for more comfortable treatments





10 LEADING

Series Radio Frequency Skin Removal plus Carbon Facials.

UNLIMITED POTENTIAL

Non-invasive and with minimal discomfort or downtime SuperLUX instantly updates your treatment menu to target new demographics, boost booking rates and increase profit margins.



ULTRALUX + ULTRALUX PRO: 12 IN-DEMAND TREATMENTS IN 1 DEVICE

UltraLUX PRO is the ultimate all-in-one device that delivers unprecedented power and scope for clinics wanting to maximise their service potential. Combining 12 of the most popular LUX Series applications into a streamlined design, this easily manoeuvrable treatment powerhouse expands your treatment menu while maximising competitive edge, booking potential, cash flow and daily productivity levels.

KEY FEATURES

- 2,000 watts of power with 12+ hours of continuous working time
- Intuitive graphic user interface with 12 inch full colour touch screen
- 6 IPL applications, SHR Hair Removal, Ultrasound Body Contouring, RF plus Sonophoresis
- Pre-programmed settings based on Fitzpatrick Type, plus manual mode options
- Optional PRO upgrade to Tattoo Removal and Carbon Facials
- Interchangeable magnetic safety filter system plus 15x35mm, 15x50mm and 12x40mm spot size options
- Integrated skin contact cooling from -7°c for more comfortable treatments

INTUITIVE

The easy-to-use, intuitive interface pre-programmed with adjustable Fitzpatrick skin type and treatment area for faster, more accurate and

ADAPTABLE

EFFECTIVE





- SHR (Super Hair Removal) Pain free, ultra-fast permanent hair removal
- Dual-Mode IPL Hair Removal For lighter pigmented hair, including shades of red and blonde • IPL Hair Removal - For smaller or difficult to treat areas not suitable for SHR
- IPL Photo Rejuvenation Visibly reduce the appearance of pores, fine lines and wrinkles
- IPL Pigmentation Correction Reduce skin blemishes, sunspots and uneven skin tone
- IPL Vascular Therapy Reduce redness conditions, vascular lesions, rosacea and facial flushing • IPL Acne Management - Treat breakouts and reduce acne causing bacteria
- RF Skin Tightening Tighten, smooth and tone lax skin on the face and body
- Ultrasound Body Contouring Sculpt and tone to reduce stubborn areas
- Sonophoresis Skin Care Infusion Superior product penetration for accelerated results
- Professional Tattoo Removal (UltraLUX PRO) Effective, permanent tattoo removal results on a broad range of ink types and colours
- Doll " or "Hollywood Facial"



INNOVATIVE, REVENUE-GENERATING TREATMENTS

• Carbon Facials (UltraLUX PRO) - Advanced skin rejuvenation commonly known as the "China



5 WAVELENGTHS OF MEDICAL-GRADE LED

MediLUX delivers unprecedented service power and treatment scope, allowing clinicians to push their LED service potential higher than ever before.

Included on the Australian Register of Therapeutic Goods (ARTG #317123), the highpowered medical-grade MediLUX delivers blue, green, yellow, red and near-infrared LED wavelengths to target a wide range of common client skin concerns.

An intuitive user interface allows clinicians to tailor treatments with either single strategic wavelengths or craft their own layered services in manual mode. 20-minute pre-programmed services designed for optimum results can also be chosen from, including Acne Management, Skin Tone Correction, Redness Reduction, Age Management and Dermal Repair.















INTUITIVE OPERATION FOR **CONSISTENT HAIR REMOVAL RESULTS**

HairLASE is engineered for skin specialists and clinicians of all skill levels seeking straight-forward, results-driven professional laser hair removal services to drive home better outcomes and higher profits in the treatment room.

With laser hair removal remaining one of the most in-demand aesthetic treatments available, HairLASE is engineered for excellence with simple pre-programmed protocols, large intuitive user interface and a comfortable, ergonomic handpiece to make revenue boosting laser hair removal treatments easier and more effective than ever.

KEY FEATURES

- Ergonomic treatment handpiece enables fast treatments and superior clinician comfort
- Long sweeping, short sweeping and stamping application modes available
- HairLASE recommends the optimal settings for each client based on Fitzpatrick skin type
- Advanced skin cooling reaches temperatures as low as -4°C to provide superior comfort
- Up to 15hz pulse frequency delivers rapid, powerful hair removal treatments
- No consumables means low running costs and high return on investment treatments
- Large 12mm x 12mm spot size is suitable for a wide range of treatment areas
- Intuitive user-friendly system and 10.4 inch colour touch screen
- Streamlined design provides ultimate portability between treatment spaces











HairLASE is designed to effectively treat and remove unwanted hair your client's unique objectives.

 Arms 	• Bikini
 Armpits 	Hairline
 Legs 	• Lip Area
Chest	Beard
 Back 	• Chin



HairLASE combines diode laser technology, medical science and ease of use to deliver consistent and reliable results-driven laser hair removal services.

3 reasons to get up and running with HairLASE • Encourages repeat business (a course of treatments is always recommended) • No consumables keeps running costs low and puts more profits in your pocket • Increases overall customer satisfaction with results they can see, touch and feel



ADVANCED SKIN RESURFACING + PRODUCT INFUSION

Thanks to HydroLUX, the exfoliation and skin renewal of microdermabrasion meets the intense product infusion of hydrodermabrasion for enhanced facial rejuvenation, operator ease, client comfort and treatment efficacy.

Treatments can be performed with a single targeted treatment solution or layered with multiple solutions to create a fully customised treatment experience.

KEY FEATURES

- Extra-large solution containers reduce time spent refilling and emptying
- 3 solution containers mean customising treatments is just 1 button press away
- Adjustable vacuum and flow provide for enhanced treatment efficacy
- Interchangeable tip options marked with tip size plus colour coded with grit intensity
- · Handy and easily accessible tip case incorporated into the design of the device
- Dry microdermabrasion to liquid infusion hydrodermabrasion footswitch
- Full colour, 8 inch display and intuitive graphic user interface make treatments easy
- Cleaning is fast and effortless with an in-built one press 2-minute cleaning cycle





CUSTOMISE TREATMENTS TO UNIQUE CLIENT NEEDS 80 grit, 100 grit, 150 grit, 200 grit and 250 grit tips available.

8MM TREATMENT TIP

Enables targeted decongestion on labial folds and lip lines.

15MM TREATMENT TIP The ideal size of tip for standard facial treatments to quickly manoeuvre around the face comfortably.





25MM TREATMENT TIP

Makes in demand larger body area treatments like back and buttocks "bacials" fast and easy.

TWIST

MEDICAL-GRADE MICRONEEDLING FROM **GERMAN NEEDLING INNOVATORS**

The next generation of microneedling and Collagen Induction Therapy (CIT) has arrived with XCellaris PRO TWIST.

Engineered and manufactured in Germany, XCellaris PRO TWIST is developed by the makers of the original Dermaroller®. Combining medical-grade technology (ARTG #332420, #332421, #332265) with unparalleled comfort, safety and precision XCellaris PRO TWIST delivers advanced innovations that enhance clinical treatment outcomes.

NEEDLES WITH 500X MAGNIFICATION



VERSATILE TREATMENTS

- Generates new collagen to promote skin rejuvenation
- Reduces the appearance of fine lines and wrinkles
- Reduces the appearance of facial and body scarring
- Reduces the appearance of enlarged pores
- Addresses unwanted stretch marks
- Boosts skin elasticity for plumper, younger looking skin
- Improves the effectiveness of topical products
- Refines and smooths uneven skin texture and tone





DISCOVER THE DIFFERENCE

- Innovative Technology -- Medical-Grade Equipment -- German Engineering -

INNOVATIVE MICRONEEDLES

from 0-2.5mm displace the skin during the puncture process. This within minutes of treatment to

REGULATED **TECHNOLOGY**

ENGINEERED FOR SAFETY

ACTIVE NEEDLE RETRACTION

INDUSTRY LEADING SPEEDS

OPTIMISED NEEDLE COUNT

A needling tip with too many needles can lead to unnecessary significant pressure to penetrate the skin. XCellaris PRO TWIST 6 needles to maximise treatment

INDUSTRY-LEADING LASER TECHNOLOGY TO **ERASE, RENEW AND REGENERATE**

TRI-BEAM™ is an industry-leading Q-Switched Nd: YAG laser designed for clinics that want to deliver premium medical-grade (ARTG #392750) services.

From full colour spectrum laser tattoo removal to skin rejuvenation and the ability to address difficult-to-treat pigmentation such as melasma, TRI-BEAM offers point-of-difference features for clinics wanting to maximise their treatment outcomes and competitive advantage.

KEY FEATURES

- Full colour spectrum laser tattoo removal with industry-leading treatment speeds
- 1064nm wavelength handpiece included to effectively target black and darker inks
- 532nm wavelength handpiece included to effectively target red and orange inks
- 585nm and 650nm wavelength handpiece options to target sky blue and green inks
- Rich-PTP™ targets a wide range of difficult-to-treat epidermal and dermal pigment
- Variable spot sizes (2mm to 10mm) with auto-sensor adjustment feature
- Zoom collimated handpiece available for faster, more comfortable treatments
- Multiple operation modes for individualised tattoo and skin corrective services
- User-friendly operation with light weight 360° rotating articulated arm











TARGET DIFFICULT TO TREAT PIGMENTATION

Caith

Rich-PTP™ (photo-acoustic toning pulse) technology is unique to TRI-BEAM and provides the ability to correct a range of difficult to treat epidermal and dermal pigmentation, including melasma. By distributing high peak energy into a rapid double pulse, optimal tissue interaction is achieved. Thermal relaxation time for targeted melanosomes is also optimised for enhanced safety.

ACHIEVE ENHANCED SKIN REJUVENATION RESULTS

has to

Re

Carbon Facials are taken to a whole new level with TRI-BEAM thanks to Gen Technique™ - a Quasi-Long Pulse (pulse width 300 µs) utilising the 1064nm wavelength for deep penetration. Gen Technique™ combines thermal lifting and toning by using the fundamentals of selective photothermolysis for dramatic results in skin revitalisation, firming and rejuvenation treatments.

10mm

2mm

Auto-detect spot size function for rapid treatments with superior clinical outcomes

ENJOY THE CONFIDENCE OF ADVANCED SAFETY

TRI-BEAM Fractional Mode delivers a high peak energy pulse while avoiding trauma to much of the epidermis. By minimising the possibility of over-treating targeted chromophores and surrounding tissue, Fractional Mode enables users to treat conditions such as Café au lait, Nevus of Ota and Solar Lentigines with a high level of confidence and minimal side effects.

OBSERV[®]520_X

GO BENEATH THE SURFACE WITH TECHNOLOGY ENHANCED SKIN ANALYSIS

OBSERV

Equipped with 9 comprehensive observation modes to view and analyse skin health from the deepest layer right through to the surface, OBSERV 520x allows skin professionals to illustrate client skin irregularities in detail and discuss the best course of treatment to return the skin to balance.

This interactive approach, combined with clearly defined visual evidence plus analysis and recommendation software, is the ideal way to build customer trust during a skin consultation. OBSERV 520x empowers skin professionals to guide the patient through confidently and knowledgeably to treatment and skin care recommendations that increase client engagement and investment.

KEY FEATURES

- 9 modes of comprehensive dermal and epidermal analysis
- Supports clinicians in providing accurate recommendations
- Builds client trust, loyalty plus treatment and skin care sales
- Consultation images are taken in less than 10 seconds
- 'Face Positioning System' for aligned before and after images
- Visually supports your expertise as a skin professional
- Intuitive iPad app facilitates comprehensive consultations
- Images and recommendations can be emailed to the client









SURFACE TEXTURE

The face is illuminated with a natural light for overall appearance evaluation of the skin

MODE

MODE Shows early signs of ageing such as fine lines, wrinkles, visible pore structures and bumpy skin







PARALLEL POLARIZATION

Reveals fine surface textures, pore structures and emerging wrinkles

CROSS POLARIZATION Unobstructed view

()

of dermal structures, vascular anomalies and pigmentations



360 LIGHT MODE

Innovatively and accurately illustrates areas of volume on the face. Ideal to showcase results with injectables like dermal fillers or technology enhanced treatments for age management.





PIGMENTATION MODE

Shows skin tone inhomogeneities such as hyperpigmentation, melasma and sunspots







Shows thinning skin structures, active sebaceous glands and the hydrolipid film





VASCULARITY MODE

Vascular anomalies, inflammation and compromised skin barriers become clear







Distinguishable fluorescence reveals hidden skin abnormalities



REB

MEDICAL-GRADE FAT REDUCTION **IN A WHOLE NEW LIGHT**

Discover the next generation of medical-grade (ARTG #367990) body contouring with ReBorn. Developed by Israeli based LightFective, ReBorn showcases forward-thinking engineering for leading fat reduction results. Innovative Power LEDs provide cutting edge non-surgical services for easy, consistent, comfortable treatments and unbeatable clinical results.

940nm IR (infrared) light is utilised during treatments, as it is effectively absorbed by fat cells (5mm-14mm under the skin) while not impacting upon the epidermis. A patented "multisource approach" for homogeneous heat distribution across the entire treatment area elevates the temperature of fat cells to 42-47°C, triggering apoptosis.



KEY FEATURES

- Clinically proven fat reduction results with safe and effective technology
- Matrix array of 63 (9x7) Power LEDs in the 940nm IR (infrared) range
- 7x5cm² spot size with 80° optical angle is x1.45 bigger than competitors
- Targets fat cells (5mm-14mm under the skin) without impacting the epidermis
- Patented "multisource approach" for homogeneous heat distribution
- Proprietary cooling technology keeps patient safe and comfortable
- 4 small, light and manoeuvrable applicators can treat multiple areas at once
- No consumables, disposables, messy gels or hidden charges
- Effortless and pain-free services are ideal as a "lunch-break" procedure
- Hands-free treatments leave clinicians free for other treatments



SUPERIOR RESULTS

Clinically proven to deliver safe, effective and consistent results for fat reduction services.



NON-INVASIVE

Comfortable, noninvasive and suitable for a wide range of client types and treatment areas.



60%-50%

<50%

ZERO DOWNTIME

ReBorn treatments last 25mins, after which clients can immediately resume daily activities.



SAFETY FIRST

Power LED technology avoids the epidermis, while integrated cooling ensures safety and comfort.

MEDICAL-GRADE RESULTS, CLINICALLY PROVEN **TECHNOLOGY**

ReBorn provides hands-free, clinically proven treatments that remove unwanted fat cells for long term fat reduction results.

Treatment protocols consist of 3 comfortable 25-minute treatment sessions, spaced 4-8 weeks apart. As the body flushes out the treated fat cells patients see results within 6-8 weeks, with optimal results visible after 12 weeks.

Studies show ReBorn achieves an average temperature on targeted fat cells that is 2.4°C higher than leading competitive fat-reduction devices. Treatments are highly effective in destroying fat cells, with improved patient comfort compared to other fat reduction devices.













Unlike other light-based fat reduction technologies, which suffer heat disparities of up to 60%, ReBorn delivers homogenous heat for enhanced treatment outcomes and more reliable results than more traditional fat reduction methods.

Power LED powered ReBorn is not associated with the risk of paradoxical adipose hyperplasia (unexpected overgrowth of fatty tissue) that can be experienced with other forms of fat reduction, such as cryolipolysis. Inbuilt contact cooling helps to prevent the complications seen with other forms of highperformance fat reduction for a safer and more comfortable client experience.



SOLVE 'THE THERMIC PUZZLE' WITH THE EVOLUTION OF THERMAL ENERGIES

Thuzzle is developed by GMV, an Italian based technology pioneer at the forefront of aesthetic and medical device development.

Utilising over 10 years of medical technology experience, GMV has developed the medical-grade (ARTG #385416) Thuzzle with a range of effective modalities designed to deliver the ultimate in total face and body reiuvenation.

KEY FEATURES

- Addresses some of the most in-demand medi-aesthetic concerns
- 5 handpieces deliver optimal spot sizes for services on face and body
- In-built temperature and motion control plus skin contact recognition
- Patented technology utilises waves across the entire spectrum
- Italian engineered to CSQ, Medical CE and TGA standards
- 10.1 inch full colour touch screen and advanced graphic user interface
- 40 pin and 64 pin fractional radial tips for versatile treatments

VERSATILE TREATMENTS



BIPOLAR RF SKIN TIGHTENING

CONTROLLED

TRIPLE WAVE

Elastic retraction,

neocollagenesis,

lipolisys and

vascularization

on decolleté,

arms and

knees.

CONTROLLED CROSSLINKED TETRA-WAVE



Elastic retraction, neocollagenesis, vascular activation and deep lipolisys on the abdomen, hips, gluteus and legs.

CONTROLLED DOUBLE WAVE



Elastic retraction and neocollagenesis, lipolisys and vascularization on the periocular and perioral region.









.

TETRA-WAVE FAT REDUCTION

Tetra-Wave technology combines 4 linked electromagnetic waves for nextgeneration circumferential fat reduction

=

FRACTIONAL REJUVENATION



SUPERFICIAL



Retraction and neocollagenesis, specifically for use on smaller treatment areas such as the facial region.

FRACTIONAL RADIAL

Skin resurfacing, lifting and targeting of uneven texture and tone such as acne scars, skin spots and enlarged pores.



PLEXR PLUS - THE NEXT GENERATION OF **MEDICAL-GRADE PLASMA TECHNOLOGY**

Plexr Plus is a medical-grade (ARTG #388450) plasma technology developed by GMV, an Italian based company at the forefront of medical innovation with over 10 years of experience in plasma technology.

An increasingly popular non-surgical cosmetic dermatology technique, plasma treatments utilise a plasma flash (small electrical charge) to create a radius of microplasma. Energy is delivered to the superficial skin cells, transitioning the targeted treatment area from a solid to a gas state without damaging surrounding tissue (a process known as sublimation).

Treatments encompass a range of versatile, non-surgical procedures that smooth, tighten, lift, plump, repair and rejuvenate the face and body.

Advanced touch screen user

interface shows settings and

WIRELESS HANDPIECES

recharge when placed in the

MULTIPLE HANDPIECES

4 hours of operation

video tutorials

docking station

100%	+2.6	70%
of patients saw	average	of patients
improvements	increase in	improved their
in upper eyelid	eyelid ratings on	ratings by
grade	all patients	3 grades

VERSATILE TREATMENTS

- Tighten and smooth fine lines and wrinkles
- · Firm lax, sagging skin on the face and body
- Lift drooping eyelids (non-surgical blepharoplasty)
- Clear superficial and cystic acne
- Treat hypertrophic and atrophic scarring
- Remove unwanted skin tags (acrochordon)
- Remove moles, pigmented lesions and angiomas

MEDICAL APPLICATIONS

- Remove senile warts (seborrheic keratosis)
- Treat keloid scarring and hyperkeratosis
- Remove deposits of cholesterol (xanthelasma)
- Perform gum correction treatments (gingivectomy)
- Perform eyelid correction (ectoprion and entropion)

treatment room

Treatment tips are the

only required consumable

ITALIAN ENGINEERED

of medical professionals to

Remove blood vessel growths (hemangiomas)





SKIN LESION REMOVAL

In-demand skin lesion services (such as moles, skin tags, senile warts and cherry angiomas) provide the opportunity to treat common client concerns and target a wider range of client demographics than traditional skin rejuvenation.

NON-SURGICAL **BLEPHAROPLASTY**

Tighten and lift the appearance of the periorbital area (upper eyelids and lower eyelids) with non-surgical blepharoplasty. Only superficial corneocytes are sublimated to ensure comfortable eyelid movement post treatment.













PLASMA SKIN REJUVENATION

FINE LINE + WRINKLE REDUCTION SKIN TIGHTENING AGE MANAGEMENT

SCAR REDUCTION SERVICES

Target 2 of the most common forms of scarring with hypertrophic (raised) and atrophic (indented) medi-aesthetic services available on Plexr Plus. Treatments for keloid scarring are also available as a medical application.









The Global Beauty Group Education HUB

A BLENDED LEARNING AESTHETIC EDUCATION PATHWAY LIKE NO OTHER

Access to the innovative online learning platform The Global Beauty Group Education HUB is included with your device to set you up for success with your new technology. Enjoy comprehensive training in the theory and practice of performing treatments, interactive learning modules plus treatment and consultation advice from The Global Beauty Group Dermal Education Team.

In addition to training programs, The Global Beauty Group Education HUB features all marketing, clinical and technical resources to help you and your business maxamise your investment. Discover consultation forms, before and afters, user manuals, step by step treatment protocols, social media posts, home care guides, client brochures and posters plus much more.

Efficient Learning

Interactive modules cover successfully and safely performing reatments, as well as treatment theory

96% Pass Rate

Fun, interactive and engaging learning modules plus online e-textbooks, quizzes and podcasts

Annual Membership 1-year membership comes free with device purchase, and

is available to renew

annually at a cost

24/7 access to s knowledge an resources when needed without

Lower Operating Costs

95% of students rate 5 stars

"Very easy to learn and EVERYTHING was very straight forward and easy to follow" - Helai Reshteen, Perth







TheGlobalBeautyGroup Technology for Beautiful Business

The Global Beauty Group p: 1300 006 607 e: info@theglobalbeautygroup.com w: theglobalbeautygroup.com.au