



Air Alternating Mattress

Patented Auto 3D Sensors & Heel Suspension!

The Air-Lo 8 Alternating Mattress Replacement leads the field in affordable prevention of pressure areas and skin breakdowns.

It's compact, powerful yet almost silent control unit can quickly support patient loads of up to 200kg whilst alternately allowing pressure to be released across the rest of the body.

When only the very best prevention, support and therapy is required, the Air-Lo 8 Mattress Replacement is the best value choice for safely caring for patients assessed at high risk.

Static head section eliminates patient dizziness and disorientation

20 transverse air cells alternate to provide maximum prevention and therapy.

They also feature permanently inflated safety sections for total peace of mind

All waterfalls are turned and welded for maximum infection control and zip protection

High quality two way stretch top cover is fully welded and can be high temperature washed for total infection control

Tilt Sensor!
automatically adds more pressure when patient is inclined

Cell disconnection valves!
allow the end four foot cells to be completely deflated to assist in suspending a patient's heels

Medical grade power cord is run inside the mattress to avoid dangerous threading under the bed or on the floor

Patented automatic sensor suite!
No more guessing correct settings!
Just press one button and walk away! After 7 minutes the correct pressure is set for patient weight and position on the mattress.
Three modes of operation with full safety timeouts to therapy mode.
Intermittant operation ensures long life for this super silent unit.

An optional internal Air-Lo Hyper-Stretch Coverlet can easily be fitted under all Air-Lo top covers to ensure maximum patient comfort and insulation



Unit C3, 15 Narabang Way Belrose, Sydney NSW 2085
t: 02 9844 5456 f: 02 9844 5445 e: info@patientsupportsystems.com

www.patientsupportsystems.com



Air-LO 8

Specifications

Control Unit:

Dimensions:	W 27 x H 21 x D 10cm
Weight	2.5kg
Mode of Operation	Non-Continuous
Power Rating	240V 50Hz 12VA
Transport Function	Hose End Cap
Air Flow	10 Litres per minute
Auto Startup	✓
Auto Sensor for Patient Weight Profile and Position	✓
Comfort Override Control	✓
Self Diagnostics	✓
Audible/Visual Alarm	✓
Mute Function	✓
Static Mode	Timeout - 60 minutes
Care Mode	Timeout - 20 minutes
Incline mode with Auto Tilt Sensor	✓
Auto Anti Tamper - 2 Second Unlock Function	✓
CE Certification and C-Tick Listed N27882	✓

Warranty

2 Years Pump
2 Year Soft Goods

ARTG: 175810

Mattress Replacement

Dimensions	L 200 x W 90 x H 23cm
Alternating Cell Type	Safety Cell in Cell
Overall Cells Count	20 - (3 static Head Cells included)
Quick release valves for last foot zone cells	4
Cell Material	Pure TPU
Top Cover	Dahlia Bi-Elastic 2 way stretch - welded
Base/s	Nylon TPU welded
Cell Cycle	1 in 2
Alternating Mode	7 minutes per alternation
Infection Control	Welded Top Cover
Positioning Handles	2 aside
Base Attachment	Side Straps

Optional Hyperstretch Coverlet

Dimensions	L 200 x W 107 x H 1cm
Attachment	Base Surround Zip
Cover Material	High Performance Micromesh
Foam	Dunlop Marathon Hypersoft HS18-35 with Ultrafresh