




Alma X-Lift

The Smartest Lift Yet

Alma **PrimeX**



Alma PrimeX

Alma Prime X is more than a device — it's the driving force behind the Alma X-Lift protocol. Designed to deliver non-invasive facial lifting, contouring with minimal downtime, and high patient satisfaction. Prime X empowers clinics to treat skin laxity and subcutaneous fat in one seamless solution.

Beyond the Trend



Dual Energy, Smarter Outcomes

Prime X combines unipolar RF and patented non-focused ultrasound to target both deep dermal structures and subcutaneous fat. This synergy allows practitioners to achieve visible lifting, tightening, and sculpting in fewer sessions — without compromising comfort or safety*.



Customisability and Maximising ROI

With the ability to add applicators based on client demand, Prime X offers a scalable solution that grows with your business. Its streamlined workflow ensures treatments are fast, repeatable, and cost-effective — empowering clinics to deliver premium results with strong margins and minimal complexity.



From Full Face to Full Potential

While Alma X-Lift is the signature facial protocol, Prime X is engineered for versatility. Its applicators are designed for small, delicate facial zones as well as larger body areas — giving clinics the flexibility to expand treatment offerings and grow with patient demand.

*Designed to maximise safety; however, a thorough assessment by a trained healthcare professional is essential. Clients must disclose all relevant medical history and potential contraindications during consultation.

Alma X-Lift



Seamless Synergy for Facial Lifting

Alma X-Lift combines two clinically proven technologies — functional ultrasound and high-frequency radiofrequency — in a single, harmonised protocol to address facial laxity from multiple angles. This synergy allows practitioners to contour, tighten, lift the face with precision, and comfort making it a powerful alternative to traditional platforms.



Ultrasound – Sculpt with Precision

Alma's patented non-focused ultrasound technology (AlmaWave) targets subcutaneous fat to reduce volume and refine facial contours. This energy creates acoustic pressure that gently disrupts adipose tissue, making it ideal for treating areas like the submental region and lower cheeks. The result is a slimmer, more sculpted appearance, all without surgery or downtime.



Radiofrequency – Rebuild from Within

The unipolar RF technology behind X-Lift delivers deep, controlled dermal heating using a frequency of 40.68MHz. This high-powered energy stimulates fibroblasts, promotes collagen and elastin regeneration, and enhances overall skin firmness. Unlike other devices, the treatment is comfortable, consistent, and leaves skin visibly smoother and tighter.



Contour – SpeedSculpt Technology

The SpeedSculpt applicator is engineered for precision. Using focused ultrasound through a compact sonotrode, it effectively reduces fat in smaller areas like the submental triangle. This step forms the foundation of the X-Lift protocol — debulking and refining the lower face before tightening.



Tighten – TuneFace RF Heating

TuneFace delivers unipolar RF energy combined with vacuum-assisted technology to heat the dermis deeply and uniformly. By stimulating collagen production in the mid to deep dermis, it visibly improves skin texture and elasticity. It's ideal for cheeks, nasolabial folds, and other delicate facial areas.




Lift – The Signature Finish

In the final stage of the X-Lift protocol, TuneFace is used to reinforce lift in the lower face. The mechanical vacuum gently elevates tissue while RF energy stimulates structural tightening, giving patients a natural, youthful contour — with no needles, no pain, and no downtime.



The technology behind the lift



Alma X-Lift

Visible Lift and Redefined Jawline

Noticeable elevation in the lower face and jawline contour after just 3 treatments, indicating effective tissue tightening and structural support.

Reduction in Skin Laxity

Improved firmness throughout the cheek and submental areas, demonstrating collagen stimulation and dermal remodeling via unipolar RF.

Improved Skin Texture and Smoothness

Smoother skin surface and more even tone across the treated area, highlighting the dual effect of ultrasound contouring and RF-based rejuvenation.



Courtesy of Alma Lasers Clinical Department



Courtesy of Alma Lasers Clinical Department

Jawline Redefinition and Submental Tightening


The lower face shows a significant improvement in contour, particularly along the jawline and under the chin — highlighting effective fat reduction and skin tightening in the submental zone.

Enhanced Mid-Face Lift

The directional lift seen through the cheek and lower face shows deeper structural support and volume repositioning, achieved through synergistic RF and ultrasound treatment.

Smoother, Firmer Skin Texture

Skin appears visibly firmer and more uniform, indicating improved dermal density and collagen regeneration after 3 treatments.



Alma Lasers Australia

G05, 241A O'Riordan Street, Mascot NSW 2020

Tel. 1300 511 862

info@alma-lasers.com.au

www.alma-lasers.com.au

© 2025. All rights reserved.

Connect with Alma

