

DRÄGER DRUGCHECK® 3000

Saliva Test Kits



The Dräger DrugCheck® 3000 significantly reduced the issue of a false negative and positive results. The device has two indicators that inform the tester that they can proceed to the next step, confidently knowing that the procedure is correct.

Rapidly tests individuals simultaneously for up to six substances, including Amphetamines, Benzodiazepines, Cannabis (THC), Cocaine, Methamphetamines and Opiates.

Benefits

- Find out within minutes if a person recently consumed certain drugs.
- It does not require electricity and can be used anywhere.
- The compact and quick oral fluid-based drug test yields reliable results affordably and efficiently.
- High sensitivity (THC down to 15ng/ml)
- Operating temperature 5 - 40°C.



(+61) 7 3361 2000
sales@pacdatasys.com.au / info@drugdetectionsolutions.com.au
pacdatasys.com.au / drugdetectionsolutions.com.au

DRÄGER DRUGCHECK® 3000

Saliva Test Kits



Applications

- Workplace health and safety.
- Mining, logistics, manufacturing, agriculture, among others.
- Law enforcement.
- Correctional facilities.
- Hospitals and rehabilitation clinics.
- Personal use.

Features

The disposable DrugCheck® 3000 cannot be manipulated and is hygienic to use.

- **Unambiguous on-the-spot drug detection:** Its compact design makes the DrugCheck® 3000 easy to transport. It can be made ready for use quickly and easily and is suitable for roadside drug checks and drug testing in safety-relevant workplaces. The test kit has no electrical parts, making the drug tester safe to use in hazardous areas and even as point-of-care-testing (POCT).
- **Fast and sensitive analysis for up to six substance classes:** Check individuals for five substance classes simultaneously with Dräger DrugCheck® 3000: cocaine, opiates, amphetamine, methamphetamines/designer drugs (e.g. Ecstasy, MDMA), benzodiazepines (e.g. in prescription medicines) and cannabis (THC).
- **Dräger Drugcheck App and QR Scan:** Before performing a test, the QR on the test can be scanned. User and subject details are entered, and test performed. Once results are visible, the app will scan the result, photograph it, and email it to your records' nominated location. Simple.



(+61) 7 3361 2000
sales@pacdatasys.com.au / info@drugdetectionsolutions.com.au
pacdatasys.com.au / drugdetectionsolutions.com.au