

Probe

Model mode	Applications	Transmit frequency (MHz)	Max. depth	View field	Array radius	Biopsy guide
Linear rectal probe						
L5FVC	Bovine Equine	B mode: 4.0/4.7/5.0/5.5/6.2	150mm	64mm		Invalid
L7FVC	Bovine Equine	B mode: 6.5/7.1/7.5/8.3/9.0	134mm	50mm		Invalid
Convex probe						
C3FC	Bovine Equine Ovine Porcine	B mode: 2.5/3.0/3.5/4.0/5.0	244mm	69°	R60	Invalid
Micro-convex probe						
C5FC	Ovine, Canine Porcine	B mode: 3.5/4.0/5.0/6.5/7.1	150mm	72°	R20	Invalid
Linear probe						
L1FC	Porcine (back fat)	B mode: 2.5/3.0/3.5/4.0/5.0	244mm	150mm		Invalid