2-FOLD MAT DESIGN FOR SAFETY, EASY STORAGE AND MULTI-FUNCTION

Y-CABLE TO REDUCE THE RISK TRIP HAZARD AND REPLACEMENT COST

LOW PROFILE

NON-SLIP MATERIAL

PORTABLE

EASY TO USE

EASY TO MOVE FROM ROOM TO ROOM

WATERPROOF

WASHABLE

INTEGRATES TO MOST NURSE CALL SYSTEMS

PLUG & PLAY OPERATION

1 YEAR WARRANTY



PVC 2-Fold Floor Sensor Mat-Corded



(TGA Registration: 223902 – RCM Approved – AS/NZS61000.6.3.2012)

PRODUCT CODES:

Floor Sensor Mat only: MED-FM-03

Floor Sensor Mat & Y Cable (X = 1, 2, or 3): MED-FM-03-6T(X)

Indigo Care provides a smart Australian designed 2-Fold Premium Floor Sensor Mat that is designed for easier storage and trip-hazard free. It is manufactured from a heavy duty, easy to clean waterproof PVC vinyl with a textured slip reduction surface. It is designed to be placed on the floor at the side of the resident's bed providing coverage for the full length of a standard bed at 1220mm in length. Multifunctional, it can also be used as a chair sensor mat as the floor sensor mat is operable in a folded position and can be placed on a large chair or moved to the sitting area and placed in front of the resident's chair.

Unlike other floor mats, our Premium Floor Sensor Mats are designed to work either with a cable (corded) or completely cordless with our Wireless Receiver and Transmitter. In addition, it has a small storage compartment for loose cables and wireless transmitters.

The "Y" cable is included with the hardwired, corded mat. The Y cable has a spare 6.3mm socket allowing the nurse call bell pendant to plug in as well so that either pressing the call button or standing on the mat will alert/trigger the nurse call system's alarm. The Y Cable is detachable to reduce trip hazard and to provide an economical replacement cost if required.

KEY SPECIFICATIONS:

- Dimensions: 610mm X 1220mm X 6mm (24" X 48" X ¼")
- Electrical Resistance < 100 Ohm
- Black PVC Yoga Pattern
- Bottom: Black non-slip
- Weight Range of 30 kg 200 kg
- 2-Fold Mat for Easy Storage
- Operable in fold position (suitable for chair situation)
- Wired Connection- use Y Cable (included)
- Wireless Connection- use RRU (wireless receiver) and RTU (cordless transmitter). RRU and RTU can be purchased as kit set or individual unit. Refer to Cordless Set-up Procedure for connection and programming.

PRODUCT ACCESSORIES:

Included with the Floor Sensor Mat: 1 x Y Cable: T1, T2, or T3 Wiring Configuration:

- T1-Tip & Ring
- T2-Mono
- T3-Tip & Sleeve

ORDER INFORMATION
Email your purchase orders to sales@indigocare.com.au

INDIGO CARE (Horizon Electronics Company, T/as) 2/93 Rushdale Street, Knoxfield, VIC, 3180 PH +61 3 9763-2988 www.horizone.com.au ABN 36007302801