

## ZOLL AED Plus®

Not all AEDs are created equal

### IMAGE



[Dropbox link to Images](#)

### MESSAGING

**SMART. READY. TRUSTED.**  
Not all AEDs are created equal

The ZOLL AED Plus® was the first AED with Real CPR Help® to help rescuers know when they are performing high-quality CPR. Today, the ZOLL AED Plus® continues to set the benchmark for design, dependability, durability, and performance.

The ZOLL AED Plus® defibrillator supports you with:

- **Real CPR Help®** provides integrated, real-time feedback on CPR compression rate and depth
- **Long-life consumables** - 5-year battery and pad life means fewer replacements, a low cost of ownership, and an AED that's ready for the rescue
- **Designed to be durable**, the AED Plus can be used in various environments

---

#### ZOLL AED Plus® Semi-Automatic

PN: 20200000102011130

---

#### ZOLL AED Plus® Fully Automatic

PN: 22300700502011130

---

#### Inclusions

- LCD screen showing voice prompt messages with step by step directions
- 10 Type 123 Lithium Batteries. Five (5) year shelf-life (8000-0807-01)
- CPR-D padz® one piece electrode pad with Real CPR Help® (8900-0800-01). Five (5) year shelf-life.
- CPR-D accessory kit including CPR barrier mask, scissors, gloves, prep razor, towel and a moist towelette.
- Soft Carry Case
- Operator's Guide and Basic Life Support Poster
- Access to ZOLL Online Learning
- (7) Seven year limited warranty

## ZOLL AED Plus®

Not all AEDs are created equal

### CONSUMABLES



[Dropbox link to Image](#)

#### CPR-D·padz®

PN: 8900-0800-01

CPR-D·padz® one piece electrode pad with Real CPR Help®. Supplied with gloves, barrier mask, scissors, razor, wet wipe and dry wipe. Five (5) year shelf life.

This revolutionary one-piece adult electrode ensures fast and accurate placement, then senses and reports the depth and rate of the chest compressions to the ZOLL AED Plus®. With a five-year shelf life, the CPR-D·padz electrode has the longest shelf life of any disposable electrode.



[Dropbox link to Image](#)

#### pedi·padz® II

PN: 8900-0810-01

pedi·padz® II Paediatric Multi-Function Electrodes

Designed for use with the AED Plus. The AED recognises when pedi·padz II are connected and automatically proceeds with a paediatric ECG and adjusts energy to paediatric levels. Twenty four (24) month shelf-life. One pair.



[Dropbox link to Image](#)

#### Type 123 Lithium Batteries

PN: 8000-0807-01

- Duracell Lithium 123A batteries for AED Plus
- 5 year life within the unit
- Sleeve of 10 batteries

### OPTIONAL EXTRAS



[Dropbox link to Image](#)

#### Standard Wall Cabinet (depth: 9")

PN: 8000-0855

Visibility of your AED Plus is essential so that it is easy to find in an emergency situation. The 9" alarmed wall cabinet holds the ZOLL AED Plus in its carry case with an extra set of electrodes.

- Dimensions: 44.4 cm W x 44.4 cm H x 22.7 cm D
- Includes: alarm and setup instructions

## ZOLL AED Plus®

Not all AEDs are created equal

### OPTIONAL EXTRAS



[Dropbox link to Image](#)

#### Standard Wall Cabinet (depth: 7")

PN: 8000-0817

Visibility of your AED Plus is essential so that it is easy to find in an emergency situation. The 7" alarmed wall cabinet holds the ZOLL AED Plus (on a bracket without the carry case) and one spare set of electrodes.

- Dimensions: 44.1 cm W x 44.1 cm H x 22.6 cm D
- Includes: alarm and setup instructions



[Dropbox link to Image](#)

#### Mounting bracket (use to mount unit on the wall)

PN: 8000-0809-01

- Mounting bracket (use to mount the ZOLL AED Plus on a wall without the carry case.)



[Dropbox link to Image](#)

#### Clear AED Plus Wall Cabinet

PN: 8000-0856

- Dimensions: 30.9 cm W x 56.6 cm H x 16.3 cm D
- Clear AED Plus Wall Cabinet



[Dropbox link to Image](#)

#### ILCOR Flush Wall Sign (AED)

PN: 8000-001260

- ILCOR Flush Wall Sign (AED)



[Dropbox link to Image](#)

#### ILCOR 3D Wall Sign (AED)

PN: 8000-001261

- ILCOR 3D Wall Sign (AED)

## ZOLL AED Plus® Trainer2

Not all AEDs are created equal

### IMAGE



[Dropbox link to Image](#)

### MESSAGING

#### Not all AEDs are created equal

Complete Simulation Control at the Instructor's Fingertips

The ZOLL AED Plus® Trainer2 allows CPR and first aid instructors to have students experience what it's like to save a life with the AED Plus. Choose from four different fixed rescue scenarios or have complete manual control of the entire rescue. A wireless remote lets you control multiple AED Plus Trainer2s.

#### ZOLL AED Plus® Semi Automatic Trainer2

PN: 8008-0050-05

#### ZOLL AED Plus® Fully Automatic Trainer2

PN: 8008-000052-05

#### Inclusions

- Trainer Wireless Remote (1 ea 8008-0007)
- CPR-D padz Training Electrodes (1 pr 8900-0804-0)
- Batteries (2-AA size and 4-D size)
- Trainer 2 Operator's Manual
- 12 months factory warranty

### CONSUMABLES



[Dropbox link to Image](#)

#### CPR-D-padz Training Electrodes

PN: 8900-0804-01

CPR-D-padz Training Electrodes with reusable "Z-design" electrode with CPR hand placement indicator and one (1) pair of disposable adhesive gels.

Note: The disposable gels must be applied to the reusable pad before use

#### Replacement Adhesive gels for CPR-D padz Training Electrode

PN: 8900-0803-01

- 5 pairs per case

## ZOLL AED Plus® Trainer2

Not all AEDs are created equal

### CONSUMABLES



[Dropbox link to Image](#)

#### pedi-padz II Training Electrodes

PN: 8900-000861-01

- 6 pairs per case



[Dropbox link to Image](#)

#### Replacement Trainer2 Wireless Remote Controller

PN: 8008-0007

- with 2AA batteries



[Dropbox link to Image](#)

#### ZOLL AED Plus® Trainer Carry Case

PN: 8000-0375-01

## ZOLL AED Plus®

Not all AEDs are created equal


### VIDEO & WEBSITE LINKS



#### ZOLL AED Plus Promotional Video


- Promotional Video – ZOLL AED Plus features and benefits.
- See why not all AEDs are created equal

<https://info.zoll.com/aed-plus>



#### ZOLL AED Plus Step-by-Step Demo

<https://bit.ly/2uUGISE>



#### Watch the ZOLL AED Plus Training Videos:

- Introducing the ZOLL AED Plus
- What's In the Box
- Demonstration
- Deliver a Shock
- Performing CPR with Real CPR Help
- Using the ZOLL AED Plus on Children & Infants
- Changing Pads and Batteries

<https://bit.ly/2YWU3rt>

#### More Information

To learn more about The ZOLL AED Plus®, please contact your local representative or visit ZOLL online at <https://bit.ly/2uSz9vx>