

# StatSensor®

## Creatinine & eGFR Hospital Monitoring Device



- Simple Test Procedure
- 1.2 µL Sample Size
- Results in 30 Seconds
- eGFR calculated by CKD-EPI/MDRD & Cockcroft-Gault equations
- No User Calibration
- No Interferences
- Designed for Hospital Use
- Accuracy and Reliability
- Connectivity to EMR

For more information on StatSensor Creatinine, please contact us on **1300 845 762** or at [poc@amsl.com.au](mailto:poc@amsl.com.au)

[amsl.com.au](http://amsl.com.au)

PR-100-128 August 2019

