






About Plasma Fibroblast

What is Plasma Fibroblast?

Plasma Fibroblast is an innovative development in aesthetic treatments using the fourth state of matter - Plasma. The Plaxpot device contains a



small metal probe, that when exposed to atmospheric gas (nitrogen) produces a tiny plasma flash between the tip and skin, leaving a tiny superficial thermal trauma. This causes the skin in its vicinity to tighten and lift – giving results noticeable after just one treatment. As new collagen and elastin is produced, over the next 90 days, the tightening effect will continue to progress.

Why is Plasma Fibroblast a revolutionary treatment?

Plasma Fibroblast offers an alternative to invasive surgical procedures, the treatment can be applied with absolute precision to areas that other modalities might not be able to target as effectively.

Additionally, Plasma Fibroblast delivers results that rival surgical procedures, with the advantage of minimal downtime and low risk of infection. The Plaxpot Plasma Pen allows you to treat areas such as eyelids that generally can't be treated with other modalities (e.g - laser, RF or IPL) within a salon or clinic environment. Results can be seen for up to 18 months with as little as two or three treatments, and minimal maintenance is required.

How Does Plasma Fibroblast Benefit Clients?

Plasma Fibroblast can be performed on the face, neck, hands and most parts of the body, and may assist with lifting and rejuvenating the appearance of droopy eyelids, lower eyebags, crow's feet, nasolabial folds, neck wrinkles, skin tags, redness and blemishes, as well as improving overall skin tone and texture by stimulating collagen and elastin production.

When a treatment is performed by a trained/certified Plaxpot practitioner, Plasma Fibroblast is a completely safe procedure that's suitable for all ages and most skin tones. The ideal treatment for anyone who wants to achieve a smoother, more youthful looking complexion with incredibly natural results, it's also perfect for those who want smoother skin without the use of fillers.

Plaxpot Attachments

Unlike other devices on the market, the Plaxpot has additional treatment tips giving therapists the ability to perform additional treatments, other that standard Plasma Fibroblast. These tips can be used to improve product absorption and diminish signs of ageing without causing a thermal trauma, and can be used alongside other modalities such as microdermabrasion, LED and epiblading to achieve superior results for clients and maximise ROI for clinics.

The Plaxpot also features variable energy outputs that allow therapists to control the intensity of the treatment, reducong the risk of scarring, hyper and hypopigmentation, and ensure a faster healing time when compared to other plasma fibroblast devices.



Plasma Arc Tip:

Treatment using a disposable sterile tip or enclosed reusable tip to generate a plasma arc.



Lifting Tip:

Plasma lifting for instantly younger-looking skin, without the post-treatment scabbing.



These tips are used in conjunction to create channels in the skin for better product absorption.

Inclusions

- Plaxpot Device
- Charging dock
- 3 x Attachments
- Plasma Arc Tip, Lifting Tip, Permeating Tip)
- 100 x Disposable Plasma Tips
- Training delivered by TADLI valued over \$2600
- 10 x 15ml DermaHeal Post Care RRP \$87 each
- 2 x N95 masks



Plaxpot Training

Course Content

- Plasma Essentials: Dive into the fundamental knowledge about Plasma and its role in aesthetics.
- Skin Interaction Techniques: Discover how Plasma interacts with various skin layers and the effects it can induce.
- Treatment Protocols & Techniques: Equip yourself with both traditional and innovative approaches for optimum results.
- Transdermal Delivery: Understand the different techniques for permeating the skin using Plasma, such as blemish treatment, lifting, and more.
- Wound Healing: Delve deep into the body's natural healing processes post-treatment and how to accelerate recovery.

Theory: 6 hours - live online, via Zoom or self paced.

Practical: 6 hours - in clinic, face to face.

Certificate of completion (PDF) provided upon assessment.



- Contraindications & Criteria: Recognise when Plasma Fibroblast treatment is not suitable, ensuring safety and integrity in all procedures.
- Enhanced Consultation Process: Master the art of consulting potential clients, setting realistic expectations, and ensuring their readiness.
- Skincare & Supplementation: Learn about topical and oral enhancements to amplify the results of Plasma treatments.
- Clinical Application: Get hands-on training on actual cases, fortifying your practical knowledge.
- Specialised Treatments: Dive into niche areas such as pigment and vascular lesion treatments, understanding specific techniques and intervals.
- Safety Protocols: Ensure client and practitioner safety with a deep understanding of health guidelines and best practices.

ROI & Finance

Plasma Fibroblast offers an unrivalled return on investment. Treatments that take as little as an hour can generate results that clients are willing to pay \$400-\$1200, making this device one of the most profitable treatments on your menu.

With each treatment costing less than \$20 to perform, the consumable costs for Plasma Fibroblast are incredibly low - meaning even more profits for your business! In fact, one trained therapist alone could pay off this device in a few weeks and generate enough income annually to cover their wage, your rent and so much more!

High Return on Investment

Treatments	Return/Week	Return/Year
2/week @ \$700	\$1400	\$72,800
4/week @ \$700	\$2800	\$145,600
8/week @ \$700	\$5600	\$291,200



Area	Pricing	Tx Time	
Upper Eyelids	\$600-1000	30 mins	
Crows Feet	\$200-600	30 mins	
Upper Lip	\$300-600	30 mins	
Lower Lip	\$300-600	30 mins	
Nasolabial	\$400-700	30-45 mins	
Jowls	\$500-800	30-45 mins	
Forehead Lines	\$400-700	30-45 mins	
Nech Bands	\$900-1500	45-60 mins	
Neck Wrinkles	\$1200-1500	45-60 mins	
Body (A5 area)	\$1200-1500	60-90 mins	
Body (A6 area)	\$600-900	60 mins	
Blood Vessels (A7 area)	\$129-149	Dependent on size of area treated	
Spot Treatments (skin tone, small blood vessels, skin tags, milia)	From \$50 and \$10 per lesion after that		



Results with Plaxpot





















About **TADLI**

Raising the standard of skin education & technology and bridging the gap between beauty, dermal & health care practitioners, ensuring an integrated approach.



Established in 2007, The Australian Dermal & Laser Institute (TADLI) offer diverse and comprehensive aesthetic device training for various manufacturers and direct to clinics. We believe that "best practice" is not negotiable and safe, superior results are the minimum standard. This is key to having credible staff members, ensuring greater client satisfaction, and guaranteeing a successful business.

TADLI are proud to offer an array of aesthetic devices that are both cost effective and results driven. Our strong focus on education and maximising treatment results, makes us a trusted partner to ensure clinics receive a comprehensive experience.

Every device includes award-winning training, providing users the expertise to



maximise potential.

Award-Winning Training:

Continuous Clinical Support:

We stand by our partners with ongoing clinical support, providing insights, solutions and advice whenever required.



Trust and Integrity:

Our core values emphasise building enduring relationships based on trust, making us a dependable partner in your clinic's aesthetic journey.



Promotion of Safe Practices:

With our extensive training and continuous support, clinics can ensure they're always delivering treatments that are safe and compliant.



Authentic Female Leadership:

Founded and run by women who've thrived in their own aesthetic clinics, we prioritise genuine care and industry expertise over mere profit-making.

Device Specifications

Size	197 x 40 x 38mm	
Screen	1.5" full-colour OLED	
Output Frequency	70kHz	
Power Level	3 steps (0.7W, 1.0W, 2.0W)	
Output Mode	Pulse (frequency: 10, 20, 40, 100, 200Hz) and continuous	
Power	Internal Li-Polymer Rechargeable Battery 7.4V, 1,000mAh	
Life Time	Fully charged, max 3 hours	
Charge Time	3 hours for full battery	
Battery	General life of recharging battery is 500 x recharging and discharging	
Warranty	l year guarantee	





Distributed by The Australian Dermal & Laser Institute (TADLI) info@dermalandlaser.com.au | (03) 8001 5004 | www.dermalandlaser.com.au