## REDUCE THE RISK OF FALLS & INJURIES

ALLOWS OUR FLOOR, BED OR CHAIR MATS TO TRANSMIT WIRELESS ALARMS TO CALL POINT

NO WIRES AROUND RESIDENT/PATIENT-NO TRIP HAZARDS

PORTABLE

**EASY TO USE** 

EASY TO MOVE FROM ROOM TO ROOM

INTEGRATES TO MOST NURSE CALL SYSTEMS VIA OUR RECEIVERS

PLUG & PLAY OPERATION

**DURABLE & RELIABLE** 

**COST EFFECTIVE** 

# Wireless Transmitter (RTU)



PRODUCT CODE:

### DE: MED-RTU-FM - 6.3mm Jack - Floor Mat MED-RTU-CM - RJ12 - Chair & Bed Mat

Indigo Care's Wireless Transmitter (RTU) is an exclusive ID encoded RF Transmitter. It communicates directly with the Wireless Receiver Box (RRU) or Dual Function Universal Call Point (UCPV5.1). This plug-and-play device converts a sensor mat with hardwired connection to a cordless mat.

The Wireless Transmitter is available in 6.3mm plug and RJ12 connector. It is designed to connect directly to floor mats, chair mats and bed mats.

Simply connect the Floor Mat Transmitter to the floor sensor mat. Then program the Transmitter to RRU or Dual Function Call Point by following the instructions printed on the RRU. The sensor mat will now trigger the call point when activated within the transmitting area and alert your staff via your nurse call system.

#### **KEY SPECIFICATIONS:**

- Box Dimensions: 43 x 41 x 12mm
- Connection: 6.3mm mono jack plug or RJ12 connector
- No External Antenna
- LED low battery indicator
- Transmits up to 30m (line of sight)
- Uses a 3V battery-CR2032
- Frequency: 433.92MHz
- Warranty: 12 months

## INDIGO CARE

INDIGO CARE (Horizon Electronics Company, T/as) 2/93 Rushdale Street, Knoxfield, VIC, 3180 PH +61 3 9763-2988 www.horizone.com.au ABN 36007302801

### **CORDLESS RANGE:**

- Cordless Receiver (RRU)
- Cordless Transmitter (RTU)
- Cordless Under Mattress Sensor Mat
- Cordless Floor Sensor Mat
- Cordless Bed Sensor Mat
- Cordless Cushion Floor Sensor Mat
- Cordless Chair Mat
- Cordless Waterproof Neck Pendant
- Cordless Waterproof Wrist Pendant
- Cordless PIR Motion Sensor
- Cordless Suspense Point
- Cordless Doorway Sensor

ORDER INFORMATION Email your purchase orders to sales@indigocare.com.au