




SecurePlus® AUTO CPAP 3.0 PRO



SecurePlus® AUTO CPAP 3.0 PRO



Professional Care

◦ Intelligent algorithm

The optimized therapy algorithm with additional respiratory-events monitoring guarantees a more comprehensive and comfortable treatment experience.

◦ Ultra-quiet

The newly developed blower and airway design enable ultra-quiet operation of the BreathCare III.

◦ Mask-fit test

A simple mask fitting test before sleep ensures the optimal fit, providing patients with greater confidence as they fall asleep.

Ergonomic Interaction

◦ Convenient touch screen


The complex operations of the CPAP device have been simplified into easy-to-use touch controls, the large and responsive screen makes settings easier than ever.

◦ Auto-all

Automized settings including auto-ramp, auto-humidity, auto-heated tube, auto-light, smart start/stop, altitude compensation and leak compensation, allowing the patients to receive treatment with minimum confusion and improved compliance.

◦ Wireless data system

A seamless cloud-based system enables all patient care network participants to find their most suitable portal to therapy data, with convenient functionality.



INTELLIGENT DATA MANAGEMENT

Managing the therapy data through multiple platforms, seamlessly monitor the treatment process on your mobile phone, office computer or simply a web browser. Meeting various demands of different roles in the patient care network.



WIRELESS TRANSFER

Comprehensive and professional respiratory therapy analysis system based on security-assured cloud server for patient management. WiFi module are required for this function.



QR CODE

User-friendly sleep quality assessment with essential parameters. Unique QR code will be generated after each therapy and can be easily interpreted by scanning with mobile app BreathSecure.



SD CARD/USB

Traditional and reliable data transfer method. Stored data can be interpreted with computer software BreathSecure Station and automatically generate comprehensive report.

SecurePlus® AUTO CPAP 3.0 PRO

Color	● White
Mode	CPAP/Auto CPAP
Noise	27±2 (dBA)
Pressure	4-20 (cmH₂O)

Package	Main device	x1
	Water tank	x1
	Heater tubing	x1
	CPAP Mask	x1
	Power adapter and Cable	x1

Warranty	Device(2 years)
	Water tank and Heater tubing(90 days)



SD card	x1
Portable Bag	x1

Wipes	x1
CPAP Filter	x1

Professional Care

Intelligent Algorithm
Multiple-events Monitoring
Ultra-quiet Noise
Efficient Humidifier
Central Apnea Detection
Mask-fit Test

Auto-all

Auto Ramp
Auto Heated Tube
Auto Light
Auto Start/Stop
Altitude Compensation
Leak Compensation

Wireless Data System

WIFI Module
Free Data-cloud Access
OTA
BreathSecure Station
BreathSecure+ APP

Packaging parameters

Model
BREATHSECURE 3.0
AUTO CPAP

QTN
(PCS/CTN)
1

G.W./N.W.
(kgs)
3.7/2.4

Size
(cm)
41x31x24



Plus Medical Pty Ltd
1A Coronation Ave, Kings Park,
NSW 2148, Australia
Tel: +61-2-9881 0368
Email: sales@plusmedical.com.au
Website: www.plusmedical.com.au
CPAP Division: www.cpapau.com





Australian Government

Department of Health and Aged Care
Therapeutic Goods Administration



Public Summary

Summary for ARTG Entry:	485651	Plus Medical Pty Ltd - CPAP system with Mask - Home CPAP unit	
ARTG entry for	Medical Device Included Class IIa		
Sponsor	Plus Medical Pty Ltd		
Postal Address	1A Coronation Avenue, Kings Park, NSW, 2148 Australia		
ARTG Start Date	9/04/2025		
Product Category	Medical Device Class IIa		
Status	Active		
Approval Area	Medical Devices		
Conditions			
- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information. - Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.			
Manufacturers			
Name	Address		
Plus Medical Pty Ltd	1A Coronation Avenue Kings Park, NSW, 2148 Australia		
Products			
1 . CPAP system with Mask - Home CPAP unit			
Product Type	Medical device system	Effective Date	9/04/2025
GMDN	60711 Home CPAP unit		
Intended Purpose	This system consists of a CPAP unit, CPAP masks that are compatible with the machine, and a pack of wipe for cleaning the masks. It is intended to improve the quality of sleep for individuals by delivering positive airway pressure. It is intended for the treatment of Obstructive Sleep Apnea (OSA) in patients weighing over 66 lb (30 kg). It is intended for home and hospital/institutional use. Manufacturer Details Address Certifi		
Specific Conditions			
No Specific Conditions included on Record			

© Commonwealth of Australia. This work is copyright. You are not permitted to re-transmit, distribute or commercialise the material without obtaining prior written approval from the Commonwealth. Further details can be found at <http://www.tga.gov.au/about/website-copyright.htm>.

Public Summary

SecurePlus[®]

CPAP MASK SERIES

Rebound Silicone Cushion

The mask is made of rebound silicone cushion, which is airtight and soft. It is designed to reduce facial pressure to make it fit and comfortable.



Adaptive Silicone Cushion

The silicone cushion has an adaptive design at the nose bridge, which can be automatically adjusted according to different face shapes, so as to better relieve the pressure on the high nose and improve user experience.

Quiet Vent-Hole

Honeycomb design to reduce airflow erosion intensity. Gentle and quiet.



Quick-Release Connector

Using a quick release connector to easily separate or connect the mask and the device.

Soft and Light Headgear

Selected soft skin-friendly fabrics, light and stable, help patients adapt to the treatment comfortably.



TECHNICAL SPECIFICATIONS ►

Picture					
Model	YF-01	YN-02	YF-02	YN-03	YP-01
Certification	TGA	TGA	TGA	TGA	TGA
Type	Full Face Mask	Nasal Mask	Full Face Mask	Nasal Mask	Pillow Mask
Therapy Pressure Range	4-30 cmH ₂ O	4-30 cmH ₂ O	4-30 cmH ₂ O	4-30 cmH ₂ O	4-20 cmH ₂ O
Size	S/M/L/XL	S/M/L	S/M/L/XL	S/M/L	S/M/L
Tube Connection	φ22mm	φ22mm	φ22mm	φ22mm	φ22mm
Flow Resistance [At 50 L/min]	≤0.5cmH ₂ O	≤0.5cmH ₂ O	≤0.5cmH ₂ O	≤0.5cmH ₂ O	≤0.5cmH ₂ O
Flow Resistance [At 100 L/min]	≤1.5cmH ₂ O	≤1.5cmH ₂ O	≤1.5cmH ₂ O	≤1.5cmH ₂ O	≤1.5cmH ₂ O
Expiratory Resistance with Anti Asphyxia Valve Open [At 50 L/min]	≤2.0cmH ₂ O/ [L/S]	--	≤2.0cmH ₂ O[L/S]	--	--
Inspiratory Resistance with Anti Asphyxia Valve Open [At 50 L/min]	≤2.0cmH ₂ O/ [L/S]	--	≤2.0cmH ₂ O[L/S]	--	--
Anti Asphyxia Valve Open Pressure	≤2.0cmH ₂ O/ [L/S]	--	≤2.0cmH ₂ O[L/S]	--	--
Anti Asphyxia Valve Closed Pressure	≤2.0cmH ₂ O/ [L/S]	--	≤2.0cmH ₂ O[L/S]	--	--
G.W./N.W	147/106g	90/60g	148/107g	129/99g	268/81g
Packing	10PCS/CTN 35x28x28cm	10PCS/CTN 35x24x24cm	10PCS/CTN 35x28x28cm	10PCS/CTN 35x28x28cm	10PCS/CTN 34x22x24.5cm

INFORMATION INTERCONNECTION



SecurePlus® AUTO CPAP 3.0 PRO PARAMETERS

Specification	BreathSecure 3.0
Model	YH-680
Picture	
Color	White/Black
Certification	TGA
Mode	CPAP Auto CPAP
Pressure	4~20cmH ₂ O
Noise	27dBA+2dBA
Screen	3.5-inch
Efficient Humidifier	•
FPS.Tech	•
Auto-all	•
Central Apnea Detection	•
Leak Reminder	•
Data Management	•
Wireless Transfer	•
Cloud Platform	•
G.W./N.W.	3.7/2.4kg
Packing	1SET/CTN 41×31×18cm
Warranty	Device(2 years) Water tank and Heater tubing(90days)



Plus Medical Pty Ltd

1A Coronation Ave, Kings Park,
NSW 2148, Australia
Tel: +61-2-9881 0368
Email: sales@plusmedical.com.au
Website: www.plusmedical.com.au
CPAP Division: www.cpapau.com