



POWER UPGRADE

WITH OVER 33 % MORE POWER

V-MAX mini

JUST PUSH IT





EASY AND FATIGUE-FREE OPERATION

EVEN ON LONGER TRIPS

THE MINI IN DETAIL

Total weight, approved	190 kg ⁽¹⁾ (Person, wheelchair, V-MAX mini)
Mounting possible	Seat width 36 cm (camber 0°)
Max. mounting width	Seat width 60 cm
Speed forward	intuitively adjustable 6 km/h
Range with one charge of the batteries	up to 15 km ⁽²⁾
Maximum possible incline	up to 15 % (8,5°) ⁽³⁾
Batterys	36 V / 6,4 Ah
Direct current motor	36V / 250 W
Weight, power unit	6 kg
Weight, interface (incl. cables)	0,8
Weight, Lithium ion battery pack	1,5 kg
Total weight without battery pack	6,8 kg
Total weight incl. battery pack	8,3 kg

- (1) Please note the maximum load stipulated by your wheelchair manufacturer.
- (2) According to ISO 7176-4 the range depends among other things, on the weight of the person to be transported, the surface of the ground, the topography, etc.
- (3) Climbing capacity at a total weight until 190 kg (person weight approx. 160 kg). The climbing capacity is dependent, among other things, on the weight of the person to be transported, the surface of the ground, the topography, etc.



AAT Alber Antriebstechnik GmbH
Postfach 100560 · D-72426 Albstadt

Phone: +49.7431.1295-0 · Fax: 1295-35
info@aat-online.de · www.aat-online.de

LIGHT IN WEIGHT

STRONG POWER ASSIST

Pushing a wheelchair can often be very labourious for the accompanying person. A lot of strength is needed, especially when going uphill. Even short trips are often exhausting and tiring.

The innovative V-MAX mini power assist therefore supports the attendant and thus enables more mobility. The user-friendly operating unit allows easy handling and ensures reaching all functions

with only one hand. The lightweight lithium-ion battery has a range of up to 15 km and weighs only 1.5 kg.



**THE POWER ASSIST
FOR YOUR WHEELCHAIR**

WITH THE MINI

BENEFIT TO THE MAXIMUM



QUICK AND EFFORT- LESSLY READY TO GO

The low weight and compact design make it quick and easy to attach and detach the drive and save space when stowed away.



SIMPLE AND EASY TO USE

In addition to the simple and clear controls, the small, powerful drive wheel ensures unparallelled legroom for the attendant.



STRONG AND ENDURING

The powerful drive in combination with optimum traction supports you perfectly when pushing, even on rough terrain and longer excursions.



booklet_v-max-mini_gb_12_2023 (1,1)

This brochure can also be displayed in big letters for easier readability. It is available for free download at www.aat-online.de. In addition, you will find product-specific user manuals with more technical data and further information on our homepage.



AAT online



Contact agents



Social Media



The add-on kits of AAT coincide with the guidelines of the European Union 2006/42/EG, 2007/47/EG and 93/42/EG.

In the interest of progress we reserve the right for technical changes

AAT is certified according to DIN EN ISO 9001:2015 and DIN EN ISO 13485:2016.



Engineered in germany



PERFECT GRIP

EVEN ON ROUGH TERRAIN