

We've got your cryogenic and cryosurgical needs covered

To compliment our service of delivered liquid nitrogen across Australia and New Zealand, we also supply storage and application equipment to assist you with a holistic solution for your facility.

BOC supplies a series of cryogenic dewars in various capacities for convenient, economical and safe storage of liquid nitrogen.

Features and benefits:

- High efficiency, super insulated
- Superb vacuum and thermal performance
- Lightweight and durable design
- Two-year general warranty
- Five-year vacuum warranty



Complementary to our liquid nitrogen dewars we offer the Brymill™ range of cryoguns for safe and efficient cryosurgical applications.

Features and benefits:

- Includes four tips with sizes of 0.04 to 0.016 inches
- Includes both a 20-gauge straight and a 20-gauge bent spray
- Patented safety autovent feature
- Streamlined relief valve maintains constant operating pressure
- Two-year warranty




BOC cryogenic storage and therapy equipment



Dewar products	DW03	DW10	DW20	DW30	DW50
Storage volume	3.9L	10.5L	20.5L	31.3L	50L
Static evaporation rate (L/day)	0.135	0.117	0.097	0.28	0.43
Neck opening (mm)	51	51	51	51	51
Height (mm)	403	570	695	560	740
Diameter (mm)	225	260	380	433	433
Net weight (kg)	2.7	5	10.9	11.3	14.8
Gross weight (kg)	6.4	15.2	31.3	42.6	63.3

Complementary products	DW0006	DW0006	DW0008	DW0008	DW0008
Rigid dipper	DW0006	DW0006	DW0008	DW0008	DW0008
Spare cork cover	DW0001	DW0012	DW0013	DW0014	DW0015
Manual discharge device	-	-	DW0018	DW0019	DW0020
Roller base	DW0021	-	DW0023	DW0024	DW0025



1. Specifications are subject to change without prior notice due to continuous product development/upgrades.
2. Measurements exclude canisters, side rings and bottom pads.
3. Static evaporation loss rate and static holding time are nominal. Actual rate may be affected by manufacturing tolerances, atmospheric condition and container history. These evaporation rates are specified at standard atmospheric conditions and without canisters. Net evaporation rate is calculated based on weighing method.

Cryosurgery products	Product	Part number
	Brymill Cry-Ac® Cryogun 300ml	B-800
	Brymill Cry-Ac® Cryogun 500ml	B-700
	Aperture A (0.04 in)	B102-A
	Aperture B (0.031 in)	B102-B
	Aperture C (0.022 in)	B102-D
	Aperture D (0.016 in)	B102-D
	Bent spray	B103
	Straight spray 1.0 in x 20 g	B107-20
	Brymill cryoplate	B303



For more information contact BOC

Australia
1800 050 999 | healthcare@boc.com | www.boc.com.au/healthcare

New Zealand
0800 699 2273 | healthcareNZ@boc.com | www.boc.co.nz/healthcare

BOC: Living healthcare

Always read and follow instructions for use.

© BOC Limited 2020. BOC is a trading name of BOC Limited. Reproduction without permission is strictly prohibited. Details given in this document are believed to be correct at the time of printing. Whilst proper care has been taken in the preparation, no liability for injury or damage resulting from its improper use can be accepted.

MP20-0328 CO 1120