

Anti-Fog Solution

Anti-fog solution for fog reduction during laparoscopic surgery. Simply apply to the end of the scope using the foam pad provided.

Anti-Fog is safe, non-toxic and non-irritating, and conveniently packaged in a 6ml bottle with adhesive foam pad for single patient use.

The foam pad is X-ray detectable.

Presentation:

- Sterile, single use
- 6ml bottle with foam pad application
- Box 25

Sterilisation: Radiation

ARTG No: 315548



PRODUCT CODE

MAAF-IBB10	Anti Fog Solution, 6ml bottle with foam pad, sterile, single patient use, Box 25
------------	--



For more information, call
Customer Service on 1300 616 366
or your local Product Specialist