


Bed Wedges Reference Guide


Small

Item	Description
117	Standard
117V	With Hook Tape
117HS	Heat Sealed
117HSV	Heat Sealed with Hook Tape




Large

Item	Description
122	Standard
122V	With Hook Tape
122HS	Heat Sealed
122HSV	Heat Sealed with Hook Tape




Small - 1/2 Length

Item	Description
118	Standard
118V	With Hook Tape
118HS	Heat Sealed
118HSV	Heat Sealed with Hook Tape




Large - 1/2 Length

Item	Description
123	Standard
123V	With Hook Tape
123HS	Heat Sealed
123HSV	Heat Sealed with Hook Tape




Small - 3/4 Length

Item	Description
119	Standard
119V	With Hook Tape
119HS	Heat Sealed
119HSV	Heat Sealed with Hook Tape




Large - 3/4 Length

Item	Description
124	Standard
124V	With Hook Tape
124HS	Heat Sealed
124HSV	Heat Sealed with Hook Tape




Small - Extra Long

Item	Description
117X	Standard
117XV	With Hook Tape
117XHS	Heat Sealed
117XHSV	Heat Sealed with Hook Tape




Large - Extra Long

Item	Description
122X	Standard
122XV	With Hook Tape
122XHS	Heat Sealed
122XHSV	Heat Sealed with Hook Tape




Small - Angled

Item	Description
117A	Standard
117AV	With Hook Tape
117AHS	Heat Sealed
117AHSV	Heat Sealed with Hook Tape




Large - Angled

Item	Description
122A	Standard
122AV	With Hook Tape
122AHS	Heat Sealed
122AHSV	Heat Sealed with Hook Tape




Extra Large

Item	Description
1245	Standard
1245V	With Hook Tape
1245HS	Heat Sealed
1245HSV	Heat Sealed with Hook Tape



Sidelying

Item	Description
116	Standard
116V	With Hook Tape
116HS	Heat Sealed
116HSV	Heat Sealed with Hook Tape



Bariatric - Heat Sealed

Item 1360

Length 48cm (19") x Width 30cm (12") x Height 15cm (6")

