



## Microlite® Nitrile Disposable Gloves

### Intended Use

To protect the user from incidental chemical and substance contamination within food and medical industry and trade based work environments.

Trade Name	Materials	Type	Colour	UOM
Microlite®	100% synthetic nitrile. No powder lubricant added.	Nitrile Examination Glove, powder free, online single chlorinated, non-sterile, ambidextrous, single-use only.	Blue	100 gloves

Surface Finish	Shelf Life	Storage Condition
Textured fingers, beaded cuff.	A shelf life of five years from the date of manufacture with the correct storage conditions.	Gloves shall maintain their properties when stored in a dry condition. Avoid direct sunlight.

Dimensions	Size	Palm width	Length (mm)
Millimetres (Inches)	Small	85 ± 5mm (3.34" ± 0.2")	240 (minimum) (9")
	Medium	95 ± 5mm (3.74" ± 0.2")	240 (minimum) (9")
	Large	106 ± 5mm (4.17" ± 0.2")	240 (minimum) (9")
	Extra Large	113 ± 5mm (4.45" ± 0.2")	240 (minimum) (9")

Thickness	Palm	Finger
Millimetres (Inches)	>0.05mm (>1.9mils)	>0.06mm (>2.3mils)

### Packaging Configuration & Labelling

Packaging	100 gloves per inner box. (by weight) 10 Boxes per Outer Carton
Box dimension	210mm(L) x 120mm(W) x 55mm(D)
Labelling	The inner boxes and outer cartons have the following information: Product description, brand, size, and net quantity of gloves. Country of origin, name and address of distributor, and single use only symbol.
Barcodes	GS1 format

### Glove Standard



Microlite® meets ATSM D6319 for Medical Exam use

### 100 Box Part No.

Small	230031
Medium	230032
Large	230033
Extra Large	230034

### Australian Government Department of Health

TGA approved Class 1 medical device

