

BLACK AIR® Nitrile Disposable Gloves

Intended Use

To conduct medical examination, diagnostic and therapeutic procedures to protect patient and user from cross contamination or infection.

Trade Name	Materials	Type	Colour	UOM	Surface Finish
BLACK AIR®	100% Nitrile (Carboxylated Butadiene Acrylonitrile)	Single Use, Non-sterile, Powder-Free Nitrile Examination Gloves.	Black	100 gloves	Ambidextrous, finger textured, beaded cuff.

Shelf Life

The gloves shall have shelf life of five years from the date of manufacture with the above storage condition.

Storage Condition

The gloves shall maintain their properties when stored in a dry condition. Avoid direct sunlight.

Glove Standard



ISO 374-5:2016 VIRUS
Marking of gloves protecting against virus, bacteria and fungi. Protection against bacteria & fungi- Pass Protection against viruses - Pass

Part No.

Small	120001
Medium	120002
Large	120003
Extra Large	120004

Dimensions	Size	Palm width	Length (mm)
Millimetres (Inches)	Small	85 ± 5mm (3.34" ± 0.2")	240 (minimum) (9")
	Medium	95 ± 5mm (3.74" ± 0.2")	240 (minimum) (9")
	Large	105 ± 5mm (4.13" ± 0.2")	240 (minimum) (9")
	Extra Large	115 ± 5mm (4.53" ± 0.2")	240 (minimum) (9")

Thickness	Cuff	Palm	Finger
Millimetres (Inches)	>0.07mm (>2.7mils)	>0.08mm (>3.1mils)	>0.09mm (>3.5mils)

Quality Specification

Water Tight Test	EN 455 Standard
Residual Powder	G1 - AQL 1.5 n = 5 (Average)

Packaging Configuration & Labelling

Packaging	100 gloves per inner box (by weight) 10 Boxes per Outer Carton
Box dimension	210mm(L) x 120mm(W) x 55mm(D)
Barcodes	GS1 format

